MicroMega K-Files / MicroMega K-Reamers / MicroMega H-Files

Instructions for use

EN

1. Indication

The medical devices (K-Files, H-Files, K-Reamers) are used for catheterisation (n°006 to 015) and root canal shaping (all other diameters) during a non-surgical endodontic treatment. K-Files 008 and 010 are also used in the exploration, initial penetration and root canal permeabilisation.

2. Characteristics

Ø	Maximum recommended number of uses (if the file is not visually damaged)
< 20	2
≥ 20	5

3. Reprocessing instructions

General recommendations	<ul style="list-style-type: none"> For all metal devices, it is recommended to use anticorrosion disinfecting and cleaning agents. For your own safety, please wear personal protective equipment (gloves, glasses and mask). Do not use cleaning or disinfecting agents containing phenol, aldehyde and alkaline composition. Always follow the instructions for use provided by the manufacturer of the products. Symbol: <table border="1" style="margin-left: 20px;"> <tr> <td style="text-align: center;">STERILE EO</td> <td style="text-align: center;">Sterilised using ethylene oxide</td> </tr> </table>	STERILE EO	Sterilised using ethylene oxide
STERILE EO	Sterilised using ethylene oxide		
Limitations on reprocessing	<ul style="list-style-type: none"> Due to the product design and the materials used, no definite limit to the maximum number of performable processing cycles can be specified. The service life of the medical devices is determined by their function and careful handling. Multiple use disinfection and re-sterilisation cycles may lead to increased risk of file separation. The user must ensure that the processing method used, including resources, materials and personnel, is appropriate and meets the applicable requirements. State of the art and national laws require that validated procedures be followed. 		
Materials needed	<ul style="list-style-type: none"> Gloves, masks, gown as recommended by the manufacturer of the cleaner and detergent Running or deionised water Disinfectant (neodisher® Septo Active) Detergent (neodisher® MediZym) Small soft brushes Container Ultrasonic bath or washer-disinfector Class B sterilisation apparatus <p><i>Remark: All materials used should be cleaned and replaced regularly. Identify materials used for each step of the process (initial treatment, cleaning or rinsing).</i></p>		

1	<p>Initial Treatment</p> <p>Place used products in a container or a wipe with tap water at 20-40°C and 1.0% neodisher® Septo Active for between 5 and 15 min.</p> <p>Rinse the products with tap water at 20-40°C for 1 min.</p> <p>If there is a pause before the next step, make sure the device stays moist by packing it in a wet wipe. Do not exceed 1 hour waiting time.</p> <p><i>Remarks:</i></p> <ul style="list-style-type: none"> Do not use fixing agents or hot water (> 40°C), as this causes fixation of residues and can impair successful cleaning. Follow instructions and observe concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or other defects of the devices).
2	<p>Preparation Before Cleaning</p> <p>If the devices have visible impurities, it is recommended to manually pre-clean by brushing under tap water at 20-40°C at least 1 min with a soft brush until the impurities disappear.</p> <p><i>Remark: Follow instructions and observe concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or other defects of the devices).</i></p>
3	<p>Visual Inspection</p> <p>Inspect the used products and discard any damaged products (broken, untwisted or abnormally bent)</p>
Manual	<p style="text-align: center;">Cleaning</p> <p>Insert the products in the ultrasonic bath.</p> <p>Run ultrasonic device for 10-30 min with tap water and 0.5-2.0% neodisher® MediZym</p> <p><i>Remarks: Follow instructions and observe water quality, concentrations and cleaning time stated by the manufacturer of the cleaning solution.</i></p>
5	<p style="text-align: center;">Rinsing</p> <p>Rinse the products with tap water at 20-40°C for 1 min.</p> <p><i>Remark: It is recommended to use deionised water.</i></p>
6	<p style="text-align: center;">Drying</p> <p>Dry the products with compressed air until product is visibly dry.</p>

Automatic	<p style="text-align: center;">Cleaning/Rinsing/Drying</p> <p>Place the instruments in the tray of the slide-in trolley of the washer-disinfector. Perform cleaning cycle with 0.2-1.0% neodisher® MediZym.</p> <p>Perform drying.</p> <p><i>Remarks:</i></p> <ul style="list-style-type: none"> Disinfection (thermal or chemico-thermal) is not required, as the products are sterilised after cleaning. Follow the instructions and concentrations given by the manufacturer of the detergent solution. Follow the instructions of the washer-disinfector and check that the success criteria after each cycle have been met as stated by the manufacturer. The final rinsing step should be performed with deionised water. For other steps follow the water quality defined by the manufacturer. Use only approved washer-disinfectors according to EN ISO 15883 that are maintained and validated regularly.
7	<p style="text-align: center;">Visual Inspection</p> <p>Inspect the used products.</p> <p>Repeat steps 4-5-6 if the product is not visibly clean or discard the damaged products.</p>
8	<p style="text-align: center;">Packaging</p> <p>Place the instruments in a paper-plastic pouch for steam sterilisation compliant with ISO 11607 and EN 868.</p> <p><i>Remarks:</i></p> <ul style="list-style-type: none"> For sharp devices that are not contained within a box, silicone tubes should be placed around the devices to prevent piercing of the packaging. Seal the pouches according to the recommendation of the pouch manufacturer. If a thermo-sealer is used, the process must be validated and the thermo-sealer must be calibrated and qualified.
9	<p style="text-align: center;">Sterilisation</p> <p>Sterilise the products with steam vapour</p> <ul style="list-style-type: none"> Apparatus: class B Minimum temperature: 132°C Minimum time: 3 min Absolute pressure: 2.2 bar Minimum drying: 20 min <p>Check the physico-chemical indicators and cycle parameters.</p> <p>The settings (temperature at 134 °C and duration of 18 min) are mandatory for prion inactivation according to French regulations.</p> <p><i>Remarks:</i></p> <ul style="list-style-type: none"> When sterilising multiple instruments in one autoclave cycle, ensure that the steriliser's maximum load is not exceeded. Place the pouches in the steam steriliser according to the recommendation given by the steriliser manufacturer. Use only pre-vacuum air removal steam sterilisers that meet the requirements of EN 13060 (class B, small steriliser) and EN 285 (full size steriliser), with saturated steam.
10	<p style="text-align: center;">Storage</p> <p>Store the products in a dry, clean and dust-free environment at temperatures specified for paper-plastic pouches by the steam sterilisation manufacturer.</p> <p><i>Remark: Check the packaging and the medical devices before use (packaging integrity, no humidity and use-by date). In case of damage, a complete rework should be performed.</i></p>

4. Disposal

After use, instruments must be placed in a secure container used to collect cutting or piercing instruments (like needles or disposable blades) as per good dentistry practices.

5. Storage and transport condition

6. Symbols

	Stainless steel material		Root canal preparation
	Quantity		Sterilizable in a steam sterilizer (autoclave) at the temperature specified
	K Files		Assortment
	K Reamers		Do not use if packaged is damaged
	H Files		Medical device

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