Indication	Increment Depth	All halogens or LED lights (with output of 550-1000 mW/cm²)	LED lights (with output 1000-2000 mW/cm²)
Core Build-Up and Class II Direct Restorations	5 mm	20 sec occlusal, 20 sec buccal, 20 sec lingual	10 sec occlusal, 10 sec buccal, 10 sec lingual
All indications listed (except Core Build-Up) and deeper Class II Direct Restorations)	4 mm	40 sec	20 sec
Anterior or shallow Class I restorations	≤ 3 mm	20 sec	10 sec

**Note:** For Class II restorations, remove the matrix band prior to the buccal and lingual curing steps.

- Contouring:** Contour restoration surfaces with fine finishing diamonds, bur or stones. Contour proximal surfaces with 3M™ SoLx™ Finishing Strips.

- Adjust Occlusion:** Check occlusion with a thin articulating paper. Examine centric and lateral excursion contacts. Carefully adjust occlusion by removing material with a fine polishing diamond or stone.

- Finish and Polishing:** Polish with the 3M™ SoLx™ Finishing and Polishing System.

#### Indirect Procedure for Inlays, Onlays or Veneers

- Dental Operatory Procedure**
  - 1.1 Shade selection: Choose the appropriate shade(s) of Filtek One Bulk Fill Restorative prior to isolation.
  - 1.2 Preparation: Prepare the tooth.
  - 1.3 Impressioning: After preparation is complete, make an impression of the prepared tooth by following the manufacturer’s instructions of the impressioning material chosen. An impressioning material, such as manufactured by 3M, may be used.

- Laboratory Procedure**
  - 2.1 Flour the impression of the preparation with die stone. Place pins at the preparation site at this time if a “triple tray” type of impression was used.
  - 2.2 Separate the cast from the impression after 45 to 60 minutes. Place pins in die and base, then cast as for a typical crown and bridge procedure. Mount or articulate the cast to its counter model on an adequate articulator.
  - 2.3 If a second impression was not sent, pour a second cast using the same impression registration. This is to be used as a working cast.
  - 2.4 Section out the preparation with a laboratory saw and trim away excess or, expose the margins so they can be easily worked. Mark the margins with a red pencil if needed. Add a spacer at this time if one is required.
  - 2.5 Soak the die in water, then with a brush, apply a very thin coat of separating medium to the preparation, let it dry somewhat, and then add another thin layer.
  - 2.6 Add the first increment of composite to the floor of the preparation, stay short of the margins, and follow the cure recommendations described in the Direct Restoration section (Step 10).
  - 2.7 Place and cure additional increments of composite. Allow for the last increment (incisal) to include the contact areas.
  - 2.8 Place the die back into the articulated arch. Add the last increment of composite to the occlusal surface. Overfill very slightly mesially, distally, and occlusally. This will allow for the mesiodistal contacts and the proper occlusal contact when the opposing arch is brought into occlusion with the uncured increment. Light cure for only ten seconds, then remove the die to prevent adhering to adjacent surfaces. Finish the curing process following the cure times in the Direct Restoration section (Step 10).
  - 2.9 With the occlusal contacts already established, begin removing the excess composite from around the points of contact. Develop the incisal and ridges as per remaining occlusal anatomy.
  - 2.10 Care must be taken when removing the prosthesis from the die. Break off small amounts of the die from around the restoration, the die stone should break away cleanly from the cured restoration, until all of the restoration is recovered.

- 2.11 Using the master die, check the restoration for flash, undercuts, and fit. Adjust as necessary, and then polish as noted above in Direct Restorations steps 11-13.

#### Dental Operatory Procedure

- 3.1 Roughen the interior surfaces of the indirect restoration.

- 3.2 Clean the prosthesis in a soap solution in an ultrasonic bath and rinse thoroughly, using water with soap and then re-cure.

- 3.3 Cementation: Cement the prosthesis using a 3M™ resin cement system.

#### Storage Conditions

1. This product is designed to be used at room temperature. If desired, the product may be warmed in a commercial warmer prior to use (no higher than 70°C/158°F, no longer than 1 hour); for capsules only.
2. The product is best stored at room temperature. If stored in cooler allow product to reach room temperature prior to use. Shelf life at room temperature is 36 months. Ambient temperatures routinely higher than 27°C/80°F may reduce shelf life. See outer package for expiration date.
3. Do not expose restorative materials to intense light.
4. Do not store materials in proximity to eugenol containing products.

#### Handling after Use

The multiple-use syringe is not intended for direct patient contact. Use clean gloves when handling the syringe. Directions for cleaning and low-level disinfection of the syringe are provided below.

#### Step 1 (Cleaning):

Use a Cav/Wipes™, or equivalent cleaning wipe, and wipe the entire surface of the device thoroughly for at least 30 seconds and until no visible soil remains on the device.

#### Step 2 (Disinfecting):

Use a new Cav/Wipes, or equivalent alcohol-quaternary ammonium disinfectant wipe, to disinfect the entire surface of the device by keeping wet for the contact time listed on the disinfectant label.

Refer to the separate restorative capsule dispenser instructions for reprocessing instructions.

#### Ingredients

Silane Treated Ceramic, Aromatic Urethane Dimethacrylate, Diurethane Dimethacrylate, Ytterbium Fluoride, Silane Treated Silica, 1,12-Dodecane Dimethycrylate, Silane Treated Zirconia, Water

#### SDS and Disposal

\*See the Safety Data Sheet (available at www.3M.com) for disposal and ingredient information.

\*See SDS for disposal information or to be disposed in accordance with the local or regional regulations.

- Please pay special attention to the disposal of contaminated waste to avoid health risks due to improper handling.

#### Customer Information

**Report a Safety Incident**
Please report a serious incident occurring in relation to the device to 3M or local regulatory authority.

\*No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

#### Warnings

3M MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user’s application. If this product is defective within the warranty period, your exclusive remedy and 3M Company’s sole obligation shall be repair or replacement of the 3M product.

#### Limitation of Liability

Except where prohibited by law, 3M will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Symbol Glossary		
Symbol Title	Symbol	Description and Reference
<b>Manufacturer</b>		Indicates the medical device manufacturer. Source: ISO 15223, 5.1.1
<b>Date of manufacture</b>		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
<b>Use-by date</b>		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
<b>Batch code</b>		Indicates the manufacturer’s batch code so that the batch or lot can be identified. Source: ISO 15223, 5.1.5
<b>Catalogue number</b>		Indicates the manufacturer’s catalogue number so that the medical device can be identified. Source: ISO 15223, 5.1.6
<b>Temperature limit</b>		Indicates the temperature limits to which the medical device can be safely exposed. Source: ISO 15223, 5.3.7
<b>Do not re-use</b>		Indicates a medical device that is intended for one single use only. Source: ISO 15223, 5.4.2

#### Caution

Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences. Source: ISO 15223, 5.4.4

#### Contains nano materials

Indicates a medical device that contains nano materials. Source: ISO 15223-1, 5.4.11.

#### Medical Device

Indicates the item is a medical device. Source: ISO 15223, 5.7.7

#### Unique device identifier

Indicates a carrier that contains Unique Device Identifier information. Source: ISO 15223, 5.7.10

#### Flammable

Product is classified as flammable. Source: Globally Harmonized Symbols for Classifying Hazardous Chemicals

#### Corrosive

Product is classified as corrosive. Source: Globally Harmonized Symbols for Classifying Hazardous Chemicals

#### Green Dot

Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62/EEC concerning the national law. Packaging Recovery Organization Europe.

#### Non-recorrged fiber

Indicates cardboard material is recyclable. Official Journal of the EC, Commission Decision (97/129/EC)

For more information see, HCBGregulatory.3M.com
---

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