General Information

3M™ ESPE® Filtek™ Supreme XTE Flowable Restorative, is a low viscosity, visible-light-activated, radiopaque, flowable nanocomposite. The restorative is packaged in capsules and syringes. It is available in a variety of tooth-colored shades. The shades offered with Filtek™ Supreme XTE flowable restorative were designed to coordinate with shades offered with Filtek™ Supreme XTE Universal Composite. Filtek™ Supreme XTE flowable restorative is a bisGMA, TEGDMA and Phenylone resin. The fillers are a combination of vitrified fluoride filler with a range of particles sizes from 0.1 to 5.0 microns, a non-agglomerated/non-aggregated surface modified 20nm silica filler, a non-agglomerated/non-aggregated surface modified 75nm silica filler, anti-a surface modified aggregated zirconia/silica cluster filler (composed of 20 nm silica and 4.1 nm zirconia particles). The aggregate has an average cluster particle size of 0.6 to 10 microns. The inorganic filler loading is approximately 65% by weight (46% by volume).

Indications

- Class III and V restorations
- Restoration of minimally invasive cavity preparations (including small, non stress-bearing occlusal restorations)
- Baseline under direct restorations
- Repair of small defects in esthetic indirect restorations
- Pit and fissure sealant
- Undercut blockout
- Repair of resin and acrylic temporary materials

Precautionary Information for Patients

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If allergic reaction occurs, seek medical attention as needed. Remove the product if necessary and discontinue future use of the product.

Precautionary Information for Dental Personnel

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If product contacts glove, remove and discard glove, wash hands immediately with soap and water and then re-glove. If allergic reaction occurs, seek medical attention as needed.

3M ESPE MSDS information can be obtained from www.3MESPE.com or contact your local subsidiary.

Instructions for Use**Preparation**

1. **Prophy:** Teeth should be cleaned with pumice and water to remove surface stains.

2. **Shade Selection:** Prior to isolation of tooth, select the appropriate shade(s) of Filtek Supreme XTE flowable restorative using a standard VITAPAN® classical shade guide. Shade selection accuracy may be enhanced via a mock-up; place the chosen shade of the restorative material on the unetched tooth. Manipulate the material to approximate the thickness and site of the restoration. Cure. Evaluate the shade match under different lighting sources. Remove the restorative material from the unetched tooth with an explorer. Repeat process until an acceptable shade match is achieved. For sealants, a contrasting shade may be desirable to enhance detection.

3. **Isolation:** A rubber dam is the preferred method of isolation. Cotton rolls and an evacuator can also be used.

Directions

1. **Pulp protection:** If a pulp exposure has occurred and the situation warrants a direct pulp capping procedure, use a minimum amount of calcium hydroxide on the exposure followed by an application of 3M™ ESPE™ Vitebond™ or Vitebond™ Plus Light Cure Glass Ionomer. Vitrebond or Vitrebond Plus liner/bases may also be used to line areas of deep cavity excavation.

2. **Adhesive System:** Follow manufacturer's instructions regarding placement, Note: Silane treatment is recommended for repair of ceramic restorations followed by the adhesive application.

3. Dispensing:

3.1. **Syringe delivery:** Filtek Supreme XTE flowable restorative can be delivered directly from the dispensing tip.
3.1.1. Protective eyewear for patients and staff is recommended when using the dispensing tip.

3.1.2. Prepare the dispensing tip: Remove cap and save. Twist the present disposable dispensing tip securely onto the syringe. Holding the tip away from the patient and any dental staff, express a small amount of Filtek Supreme XTE flowable restorative to assure that the delivery system is not plugged.

3.1.3. If plugged, remove the dispensing tip and express a small amount of Filtek Supreme XTE flowable restorative directly from the syringe. Remove any visible gel, if present, from the syringe opening. Replace dispensing tip and again express restorative. Filtek Supreme XTE flowable restorative may be extruded onto a dispensing pad and applied with a brush or other appropriate instrument.

3.2. **Single-Dose Capsule:** Insert capsule into 3M™ ESPE™ Restorative Dispenser. Refer to separate restorative dispenser instructions for full instructions and precautions. Extrude restorative directly into cavity.

4. Placement:

4.1. Place and light cure restorative in increments as indicated in Section 5.
4.2. Avoid intense light in the working field.

4.3. To prevent oozing of paste from the syringe tip after placement, pull back on syringe plunger.

5. **Curing:** This product is intended to be cured by exposure to a halogen or LED light with a minimum intensity of 400 mW/cm² in the 400-500nm range. Cure each increment by exposing its entire surface to a high intensity visible light source, such as a 3M™ ESPE curing light. Hold the light guide tip as close to the restorative as possible during light exposure.

Current Time			
Shades	Increment depth	All halogen lights LED lights (with output 400-1000 mW/cm ²)	Elipar™ S10 and Elipar™ Freelite 2 (LED lights with output 1000-2000 mW/cm ²)
Opaque shades	1.5 mm	40 sec.	20 sec.
All other shades	2.0 mm	20 sec.	10 sec.

6. Complete the restoration:

6.1. **Base/liner application:** Place a composite restorative material, such as Filtek Supreme XTE Universal restorative directly over the cured Filtek Supreme XTE flowable restorative. Follow the manufacturer's instructions regarding placement, curing, finishing, occlusal adjustment, and polishing.

6.2. **Direct Restorative Application:**

6.2.1. Contour restoration surfaces with fine finishing diamonds, burs or stones.

6.2.2. Check occlusion with a thin articulating paper. Examine centric and lateral excursion contacts. Carefully adjust occlusion by removing material with a fine polishing diamond or stone.

6.2.3. Polish with 3M™ ESPE™ Sol-Lex™ Finishing and Polishing System or with white stones, rubber points or polishing paste where discs are not suitable.

6.3. **Pit and Fissure Sealants:** Gently remove the inhibited layer remaining after light curing with slurry of pumice or polishing paste.

Storage and Use

1. This product is designed to be used at room temperature. If stored in cooler allow product to reach room temperature prior to use. Shell life at room temperature is 24 months. Ambient temperatures routinely higher than 27°C/80°F may reduce shelf life. See outer package for expiration date.

2. Do not expose restorative materials to elevated temperatures, or to intense light.

3. Do not store materials in proximity to eugenol containing products.

4. Discard used syringe tip and replace with syringe storage cap.

Disinfect this product using an intermediate level disinfection process (liquid contact) as recommended by the Centers for Disease Control and endorsed by the American Dental Association. Guidelines for Infection Control in Dental Health-Care Settings - 2003 (Vol. 52; No. RR-17), Centres for Disease Control and Prevention.

Disposal – See the Material Safety Data Sheet (available at www.3MESPE.com or through your local subsidiary) for disposal information.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Caution: U.S. Federal Law restricts this device to sale or use on the order of a dental professional.

Warranty: 3M ESPE warrants this product will be free from defects in material and manufacture.

3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be to repair or replace the 3M ESPE product.

Limitation of Liability: Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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