3M | Impregum™ Penta[®] ESPE

- (en)** Polyether Impression Material
- (es)** Material de impresión a base de poliéter
- (pt)** Material de poliéter para impressões
- (ko)** 폴리에테르 인상재
- (zh)** 聚乙醚印模材料
- (th)** วัสดุพิมพ์ปากชนิดโพลีเอทาร์
- (tr)** Matériau pour empreintes, à base de polyéther

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3M Deutschland GmbH
Health Care Business
Carl-Schurz-Str. 1
41453 Ness - Germany

3M, ESPE, Garant, Impregum, Penta, Pentamix, Pentamic and Permadyne are trademarks of 3M or 3M Deutschland GmbH. © 2020, 3M. All rights reserved.

(en) ENGLISH
Product Description
Impregum™ Penta [®] is a medium-bodied consistency polyether impression material for the Pentamix™ mixing device. The mixing ratio is 5 volumes base paste : 1 volume catalyst.
⚠️ These Instructions for Use should be kept for the duration of product use. The product may only be used when the product labeling is clearly readable. Refer to all additionally mentioned products please refer to the corresponding Instructions for Use.
Intended Purpose
Intended purpose: dental impression material. Intended users: educated dental professionals who have theoretical and practical knowledge on users of dental products. Patient target group: all patients requiring dental treatment, unless the patient’s condition limits the use. Clinical benefit: create a precise reproduction of the clinical situation within dental workflows.
Indications
<ul style="list-style-type: none">Impression of inlay, onlay, crown, bridge, and veneer preparations Functional impressions Fixation impressions Implant impressions
Contraindications
Polyether impression materials are contraindicated for use with patients who have a known allergy to polyether impression materials or contained allergens, as listed in section Precautionary Measures.
Precautionary Measures
For Patients
The base paste and catalyst paste of this product contain substances that may cause an allergic reaction by skin contact in certain individuals. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If an adverse reaction is observed in a patient discontinue use of the product in this patient and seek medical attention as needed. Known allergens are 1-Dodecylimidazole, Grapefruit flavors (different flavors) and Sulfonium salt.
For Dental Personnel
The base paste and catalyst paste of this product contain substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a no-touch technique is recommended.
Potential Undesirable Side Effects and Complications
Potential effects are local irritations, or local or systemic allergic reactions. In rare cases impression material can be swallowed or aspirated. Swallowing of impression material can cause harm to the gastrointestinal tract. As the impression material cures in the mouth of the patient, it is possible that loose adhesive restorations are removed with the impression material. Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority. 3M Safety Data Sheets/Safety Information Sheets can be obtained from www.3M.com, or contact your local subsidiary.
Preparation
<ul style="list-style-type: none">Block-out undercuts or areas where gingival recession exists to prevent the impression material from “backing” into the tooth structure. Failure to block-out may make tray removal difficult, or cause extraction of natural teeth or prosthesis. The block-out is especially important when using impression materials with high Shore hardness.

Impression trays: Any impression trays generally used for precision impressions are suitable. For sufficient adhesion, apply a thin layer of Polyether Adhesive or Polyether Contact Tray Adhesive to the tray and allow to dry completely according to the corresponding Instructions for Use.

Pentamix/Penta cartridge/poly bag:

- Place Impregum Penta in the designated cartridge only!
- Poly bags equipped with a **PentaMatic** sealing cap exclusively must be used with Pentamix™ Mixing Ties Press.
- Discard the first 3 cm of paste extruded from newly filled cartridges prior to the first use using impression tubing. The color of the dispensed paste must be uniform.
- If a new mixing tip is installed when the cartridge is inserted, prior to mixing check whether the drive shaft is engaged.

Retraction
Solutions based on aluminum hydroxide chloride or aluminum sulphate are suitable retraction agents. Retraction with threads or rings containing epinephrine (adrenaline), 8-hydroxyquinoline sulphate or iron (III) sulphate may impair the setting of polyether impression materials.

- Areas from which impressions are to be taken should be kept dry.
- In subgingival preparations, hemostatic threads or rings may be used.
- Prior to taking the impression, completely remove all residual retraction agent by rinsing and drying.

Dosing and Mixing

- Dosing and mixing are performed automatically in the Pentamix.

Times

	Working Time from start of mixing* min:sec	Setting Time from start of mixing* min:sec	Intraoral Setting Time min:sec
Impregum Penta	02:45	06:00	03:15

The working times shown are in conformity with ISO 4823. Higher temperatures will shorten the working time, while lower temperatures will prolong it. We recommend a working time of 01:45 min:sec and compliance with the indicated setting time from the start of mixing to assure an ideal impression result of all times under clinical conditions.
* Start of mixing = entry of paste into the mixing cannula.

Impression Taking

- If necessary for removal of the set impression from the mouth appropriately block-out deep undercuts in interdental areas.

Nono-phase technique

- For application according to the preparation attach the Penta Elastomer syringe to the Penta Mixing Tie Rod of the Pentamix device, and fill the syringe.
- Then load the tray previously prepared with adhesive. The mixing tip should be kept immersed in the material at all times.
- Apply the material into the sulcus or cavity from bottom up, keeping the tip of the application nozzle immersed in the material and in contact with the tooth surface at all times.
- Depending on tray size, the application of the material should be initiated only once the tray is half-loaded in order to ensure that the tray can be completely filled. In situations difficults handling, use the syringe material may set more rapidly than the tray material causing distortion of the impression.

- Immediately after application around the preparation, position the loading tray in the mouth, and hold it there without applying pressure until setting is complete.
- To remedy initial adhesion (‘setting the valve’), especially with upper jaw impressions, remove the tray from the gingiva on one side in a posterior position. If this proves difficult it may be necessary to carefully blow some air or water between the impression and the gingiva.

Functional impression

- After preparing the custom tray with adhesive, load the tray with paste mix, position it in the mouth, and ask the patient to perform functional exercises.
- If required, use Permadyne Penta L or Permadyne Garant 2:1 for precise impression taking (correction).

After Impression Taking

- Thoroughly examine and explore the sulcus of the prepared teeth and surrounding dentition. Remove any residual cured impression material from the mouth.

Hygiene

- Place the impression in a standard disinfectant solution commonly used for impressions, for the period of time recommended by the manufacturer. Excessive disinfection may damage the impression.
- After disinfection, rinse the impression under running water for approx. 15 sec.

Model Preparation

- Prepare a cast from the impression with a specialized stone plaster no earlier than 30 min and no later than 14 days after impression taking.
- To avoid introducing bubbles into the model, briefly pre-rinse with water and dry with air before pouring. Do not use surfactants as these impair the quality of polyether impressions and are not recommended.
- Polyether impressions can be silver-coated, whereas copper-coating is not feasible.

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