<b>3M   Impregum™ Soft</b>	<b>ESPE</b>
<b>☞ Polyether Impression Material</b>	
<b>☞ Material de Impresión de Poliéter</b>	
<b>☞ Material de moldagem a base de poliéter</b>	
<b>☞ 폴리에스테르 치형재</b>	
<b>☞ 聚乙 醜印模材料</b>	
<b>☞ วัสดุพิมพ์ยางชนิดโพลีเอเทอร์</b>	

Instructions for Use
Información de empleo
Instruções de Uso

3M Deutschland GmbH
Health Care Business
Carl-Schurz-Str. 1
41453 Neuss - Germany

3M, ESPE and Impregum are trademarks of 3M or 3M Deutschland GmbH.

© 2020, 3M. All rights reserved.

<b>☞ ENGLISH</b>	
------------------	--

**Product Description**
Impregum Soft is a polyether impression material that is mixed by hand and is available in heavy-bodied, medium-bodied and light-bodied consistency. The base paste is cavity-filling, medium to 7:1 by volume.
Please refer to the instructions for use, use Polyether Adhesive/Polyether Contact Tray Adhesive. This adhesive can be used on metal or methacrylate trays, but not with vinyl polysiloxane impression materials.
☞ These Instructions for Use should be kept for the duration of product use. The product may only be used when the product labeling is clearly readable. For details on all additionally mentioned products please refer to the corresponding instructions for use.

**Intended Purpose**
Intended purpose: dental impression material.
Intended users: educated dental professionals who have theoretical and practical knowledge on usage of dental products.
Patient target group: all patients requiring dental treatment, unless the patient's condition limits the use.
Clinical benefit: create a precise reproduction of the clinical situation within dental workflows.

**Indications**
• All types of precision impressions

**Contraindications**
Polyether impression materials are contraindicated for use with patients who have a known allergy to polyether impression materials or contained allergens, as listed in section Precautionary Measures.

**Precautionary Measures**
**For Patients**
The base paste and catalyst paste of this product contain substances that may cause an allergic reaction by skin contact in certain individuals. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If an adverse reaction is observed in a patient discontinue use of the product in this patient and seek medical attention as needed.
Known allergens are 1-Dodecylimidazole, mint flavors and Sulfonium salt.

**For Dental Personnel**
The base paste and catalyst paste of this product contain substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a no-touch technique is recommended.

**Potential Undesirable Side Effects and Complications**
Potential effects are local irritations, or local or systemic allergic reactions. In rare cases impression material can be swallowed or aspirated. Swallowing of impression material can cause harm to the gastrointestinal tract. As the impression material cures in the mouth of the patient, it is possible that loose teeth or restorations are removed with the impression material.

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.
3M Safety Data Sheets/ Safety Information Sheets can be obtained from www.3M.com or contact your local subsidiary.

**Preparation**
Block-out undercuts or areas where gingival recession exists to prevent the impression material from "locking" onto the tooth structure. Failure to block-out may make tray removal difficult, or cause extraction of natural teeth or prostheses.

**Impression Tray:**
Any impression trays generally used for precision impressions are suitable.
☞ For sufficient adhesion, apply a thin layer of Polyether Adhesive/Polyether Contact Tray Adhesive to the tray and allow to dry completely according to the Instructions for Use.

**Retraction**
Aluminum hydroxychloride- or aluminum sulfate-based solutions are suitable retraction agents. Retraction using coils or rings containing epinephrine (adrenalin), 8-hydroxy quinoline sulfate and iron (III) sulfate can interfere with the setting of the polyether impression materials.
► Keep the areas from which the impression is to be taken dry.
► For subgingival preparations use cords or rings as necessary.
► Before taking the impression, thoroughly rinse and dry the area to remove the rest of the retraction agent.

**Dosing and Mixing**
► Place strands of base paste and catalyst of equal length on a mixing pad. Using too much or not enough catalyst will not affect the working time, but will adversely affect the quality of the impression.
► Mix the paste with a spatula for about 45 seconds. Wipe any unmixed paste remaining on the spatula on the edge of the pad, pick up the paste again and spread it until it is a uniform color. Never mix by stringing, as bubbles may form.

**Times**
The working times apply at 23°C/74°F. Higher temperatures decrease the total working time, lower temperatures increase it.

	Heavy-Bodied Consistency min:sec	Medium-Bodied Consistency min:sec	Light-Bodied Consistency min:sec
Mixing	00:45	00:45	00:30
Working from beg. of mixing	02:30	02:30	03:00
Setting from beg. of mixing	05:30	06:00	06:00
Time in the mouth*	03:00	03:30	03:00

\* When products are combined, the longer time applies.

- Taking the Impression**
**Monophase technique**
► Using the filling device, fill the Elastomer Syringe with the mixed paste of Impregum Soft medium-bodied consistency.
► Then fill the tray prepared with Polyether Adhesive/Polyether Contact Tray Adhesive.
► Springe the sulcus or cavity from the deep part out. Always keep the tip of the intraoral tip immersed in the material and touching the tooth surface during application.
► Immediately after syringing, position the filled tray in the mouth and hold it there without pressure until the material has set.
- Double-Mix Technique**
► Using the filling device, fill the Elastomer Syringe with the mixed paste of Impregum Soft light-bodied consistency.
► Then fill the tray prepared with Polyether Adhesive/Polyether Contact Tray Adhesive using Impregum Soft heavy-bodied consistency.
► Springe the sulcus or cavity from the deep part out. Always keep the tip of the intraoral tip immersed in the material and touching the tooth surface during application.
► Immediately after syringing, position the filled tray in the mouth and hold it there without pressure until the material has set.

**Functional Impression**
► Fill the individual tray prepared using adhesive with Impregum Soft heavy-bodied consistency or Impregum Soft medium-bodied consistency. Position the tray in the mouth and have the patient make functional movements.
► If necessary, use Impregum Soft light-bodied consistency subsequently for detailed impressions.

**After Impression Taking**
Thoroughly examine and explore the sulcus of the prepared teeth and surrounding dentition. Remove any residual cure impression material from the mouth.

- Higiene**
► Place the impression in a standard disinfectant solution commonly used for impressions for the period of time recommended by the manufacturer. Excessive disinfection may damage the impression.
► After disinfection, rinse the impression under running water for approx. 15 sec.

- Pouring a Model**
► Prepare a cast from the impression with a commercial specialized stone plaster no earlier than 30 min and no later than 14 days after taking the impression.
► To avoid introducing bubbles into the model, briefly pre-rinse with water and dry with air before pouring. Do not use surfactants as these impair the quality of polyether impressions and are not necessary!
► Polyether impressions can be silver-coated, whereas copper-coating is not feasible.

**Cleaning and Disinfection**
**Multiple Use and Cleaning of Medical Device**
If the tube is not visibly clean after use, i.e., if, contrary to expectations, there is contamination with blood or saliva during the treatment, for example, discard the tube immediately.
**Specifications for Product Use**
Before use of the tube, put on new, unused disposable gloves during all clinical applications and the tube to reduce the risk of contamination. The material should be applied in advance to prevent contamination during the treatment.

**Cleaning**
The tube must not be placed into a disinfection bath or a cleaning-disinfection device (washer-disinfector). A proper cleaning of the tube can be carried out as part of hygiene management, i.e., in the event of soiling from the environment, e.g., dust (soiling does not mean contamination with saliva, blood, secretion, etc.), CavWipes™ have been tested for use as a cleaning agent. Always observe all applicable legal and hygiene regulations for dental offices and/or homes. Cleaning is done using ready-to-use disinfectant wipes (e.g., CavWipes™) based on the active substances alcoholols and quaternary compounds; alternatively, use equivalent materials. Follow the cleaning agent manufacturer's directions for use. Clean the medical device by wiping it for one minute with a new cloth.

**Inspection, Maintenance, Testing**
Examine the medical device for damage, discoloration, and contamination before every use. Do not use damaged and contaminated devices under any circumstances.

**Cleaning the Elastomer Syringe**
Please refer to the pertinent Instructions for Use for information about cleaning and disinfecting the Elastomer Syringe.

**Cleaning**
► Paste that has not set can be removed using ethanol or by washing with soap and water.
► Acetone can be used to clean the adhesive from metal trays.

- Notes**
• Elastomer Syringe Tips are single-use devices and would be damaged when removing soft impression material.
• Please keep the single-use accessory in the original bag to guarantee the traceability of the batch.
• Due to increased viscosity, the paste cannot be taken from the tube at temperatures below 18°C/64°F. When returned to room temperature, they become workable again with no loss of quality.
• Direct sunlight and humid storage conditions limit the impression.
• Polyether impressions should never come in contact with liquids containing solvents. This could result in swelling and make the model inaccurate.
• 3M polyether materials may only be combined with other 3M polyethers. The working and setting times of the individual products must necessarily be strictly observed.

**Disposal**
Dispose of the contents or container in accordance with the applicable regulations. Please pay special attention to the disposal of contaminated waste to avoid health risks because of improper handling.

**Storage and Stability**
Keep the product at 18-25°C/64-77°F. **Do not refrigerate!**
Do not use after the expiration date.
Store impressions dry and below 30°C/86°F in the dark.

**Customer Information**
No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

**Warranty**
3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

**Limitation of Liability**
Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

<b>Reference Number and Symbol</b>	<b>Symbol</b>	<b>Description of Symbol</b>
ISO 15223-1 5.1.1 Manufacturer		Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/CE.
ISO 15223-1 5.1.3 Date of Manufacture		Indicates the date when the medical device was manufactured.
ISO 15223-1 5.1.4 Use-by date		Indicates the date after which the medical device is not to be used.
ISO 15223-1 5.1.5 Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified.
ISO 15223-1 5.1.6 Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified.
ISO 15223-1 5.3.7 Temperature limit		Indicates the temperature limits to which the medical device can be safely exposed.
ISO 15223-1 5.4.4 Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
CE Mark		Indicates conformity to European Union Medical Device Regulation or Directive with notified body involvement.
Medical Device		Indicates the item is a medical device.
Recycle		Indicates packaging is recyclable with mixed paper.)
Aluminum		Indicates metal component made of aluminum is recyclable.
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national Law. Packaging Recovery Organization Europe.

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