## Imprint™ 3

(en) **VPS Impression Material**

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(es) **Material de Impresión VPS**

Instructions for Use
Mode d’emploi
Instrucciones de uso

Imprint™ 3 Penta™ Quick Step Heavy BodyImprint™ 3 Quick Step Light BodyImprint™ 3 Quick Step Regular BodyImprint™ 3 Quick Step Ultra-Regular BodyImprint™ 3 Penta™ PuttyImprint™ 3 Penta™ Heavy BodyImprint™ 3 Light BodyImprint™ 3 Regular BodyImprint™ 3 Ultra-Regular Body

**3M Deutschland GmbH**
Dental Products
Carl-Schurz-Str. 1
41453 Neuss – Germany

Imprint™ 3 Quick Step Heavy BodyImprint™ 3 Heavy Body

**3M ESPE**
Dental Products
2510 Conway Avenue
St. Paul, MN 55144-1000 U.S.A.

3M ESPE Customer Care/MSDS Information:
U.S.A. 1-800-634-2249 and Canada 1-888-363-3685.

3M, ESPE, Garant, Imprint, Penta and Pentamix are trademarks of 3M Company or 3M Deutschland GmbH.
Used under license in Canada. © 2015, 3M. All rights reserved.

(en) **ENGLISH**

**Product Description**

These Instructions for Use apply to the subsequently listed VPS Impression materials, which are arranged for One-Step Impression Technique.

Product Name	ISO 4823	ISO 4823 Consistency	Paste Color	Mixing Ratio (Vol.)	Mixing Device
Imprint™ 3 Penta™ Quick Step Heavy Body	Type 1	heavy-bodied	blue green	5 <span> </span> : 1	all Pentamix versions
Imprint™ 3 Quick Step Heavy Body	Type 1	heavy-bodied	blue green	1 <span> </span> : 1	Garant™ Dispenser
Imprint™ 3 Quick Step Light Body	Type 3	light-bodied	pink	1 <span> </span> : 1	Garant Dispenser
Imprint™ 3 Quick Step Regular Body	Type 3	light-bodied	pink	1 <span> </span> : 1	Garant Dispenser
Imprint™ 3 Quick Step Ultra-Regular Body	Type 2	medium-bodied	caramel	1 <span> </span> : 1	Garant Dispenser

Product Name	ISO 4823	ISO 4823 Consistency	Paste Color	Mixing Ratio (Vol.)	Mixing Device
Imprint™ 3 Penta™ Putty	Type 0	putty	frog green	5 <span> </span> : 1	Pentamix 2 and younger versions*
Imprint™ 3 Penta™ Heavy Body	Type 1	heavy-bodied	burgundy	5 <span> </span> : 1	all Pentamix versions
Imprint™ 3 Heavy Body	Type 1	heavy-bodied	burgundy	1 <span> </span> : 1	Garant Dispenser
Imprint™ 3 Light Body	Type 3	light-bodied	yellow	1 <span> </span> : 1	Garant Dispenser
Imprint™ 3 Regular Body	Type 3	light-bodied	yellow	1 <span> </span> : 1	Garant Dispenser
Imprint™ 3 Ultra-Regular Body	Type 2	medium-bodied	petrol	1 <span> </span> : 1	Garant Dispenser

\* Imprint 3 Penta Putty can only be mixed in Pentamix versions from Pentamix 2 (September 1999) on, earlier Pentamix units (before September 1999) extrude the material too slow.

Imprint 3 Penta Heavy Body and Imprint 3 Penta Quick Step Heavy Body can be mixed in all versions of Pentamix mixing units.

All other products are statically mixed and directly applied using the Garant™ Dispenser.

Not all products are available in all countries.

(fr) These Instructions for Use should be kept for the duration of product use. For all additional mentioned products please refer to the corresponding Instructions for Use.

**Areas of Application**

- All dual phase impressions, e.g. crown, bridge, inlay and onlay preparations, implant abutments, orthodontic impressions and related impressions such as a matrix for provisional restorations, tooth reduction guides, study models, bleaching trays, etc.

The quick setting products are especially suited for single-unit and two-unit impressions.

**Precautionary Measures**
3M MSDSs can be obtained from www.mmm.com or contact your local subsidiary.

**Preparation**

- Block out undercuts or areas where gingival recession exists to prevent the impression material from “locking” onto the tooth structure. Failure to block-out may make tray removal difficult, or cause extraction of natural teeth or prosthesis. This is especially necessary when using Imprint 3 Penta Putty or any other rigid-setting impression material.

**Impression Tray:**
All impression trays generally used for precision impressions are suitable, e.g. full arch, quadrant and dual-arch impression trays. Rigid impression trays are recommended when using Imprint 3 Penta Putty.

- For sufficient adhesion, apply a thin layer of VPS Tray Adhesive, on the tray and allow it to dry completely (5–15 minutes drying time are ideal).

**Using Penta Material:**

- Place the Penta impression material foil bags into the designated metal-reinforced Penta cartridges.
- Insert the loaded cartridge into the Pentamix mixing unit.
- Attach a red Penta™ Mixing Tips.
- Turn the plungers of the mixing unit downward until resistance is noticed using the hand-wheel, do not force.
- When using new foil bags: Initiate mixing and extrude, bleed the first few amounts of the mix until it is uniform in color, then discard.
  - Please see “Product Description” for paste color.

- Material that has set inside the Penta Mixing Tip should not be extruded by force as this may damage the Penta Cartridge and Penta Mixing Tip and result in the formation of leaks.

- Please follow the section “Impression Taking” for tray filling procedure.

**Using the Garant Dispenser:**

- Place the cartridge of the tray and wash material in the Garant Dispenser each. When using Penta material for tray filling, only place the cartridge of the wash material in the Garant Dispenser.
- Check before application that both Garant cartridge openings are not clogged.
- Using new Garant cartridges: Extrude a small amount of paste until both, the base and catalyst paste flow out evenly.
- Attach a green Garant™ Mixing Tip, to the cartridge of the heavy-bodied impression material and a yellow Garant Mixing Tip to the cartridge of the wash impression material.
- For intraoral application of the wash impression material attach a yellow Garant Intraoral Tip onto the yellow Garant Mixing Tip.
- If needed: Enlarge the exit opening for the intraoral tip by cutting it with a scalpel.
- Make sure that base paste and catalyst are mixed completely and is being extruded in a uniform color.
  - Please see “Product Description” for paste color.
- Material that has set inside the Garant Mixing Tip should not be extruded by force as this may damage the Garant cartridge and Garant Mixing Tip and result in the formation of leaks.

**Retraction**
Suitable retraction agents include aluminum chloride, aluminium sulfate or ferric sulfate solutions.

- Keep the preparation area clear of any pooled or excess fluids while the impression is being taken.
- Hemostatic cords may be used for sub-gingival preparations.
- Thoroughly remove the residue of the retraction/hemostatic agent by rinsing before taking the impression.

	Working Time at 23° C/74° F from start of mixing* min:sec	Intraoral Syringe Time min:sec	Intraoral Setting Time min:sec
Symbols on product			
Quick Setting			
Imprint 3 Penta Quick Step Heavy Body	1:15	–	2:30
Imprint 3 Quick Step Heavy Body	1:15	–	2:30
Imprint 3 Quick Step Light Body	1:30	0:40	2:30
Imprint 3 Quick Step Regular Body	1:30	0:40	2:30
Imprint 3 Quick Step Ultra-Regular Body	1:30	0:40	2:30
Intraoral setting time for all combinations of quick setting materials:			2:30
Regular Setting			
Imprint 3 Penta Putty	1:30	–	3:00
Imprint 3 Penta Heavy Body	2:00	–	3:30
Imprint 3 Heavy Body	2:00	–	4:00
Imprint 3 Light Body	2:00	1:00	3:30
Imprint 3 Regular Body	2:00	1:00	3:30
Imprint 3 Ultra-Regular Body	2:00	1:00	3:30
Intraoral setting times for combinations of regular setting materials:			
Imprint 3 Penta Putty	Imprint 3 Light Body, Imprint 3 Regular Body or Imprint 3 Ultra-Regular Body		3:00 **
Imprint 3 Penta Heavy Body	Imprint 3 Light Body, Imprint 3 Regular Body or Imprint 3 Ultra-Regular Body		3:30
Imprint 3 Heavy Body	Imprint 3 Light Body, Imprint 3 Regular Body or Imprint 3 Ultra-Regular Body		4:00

\* Start of mixing = Entry of paste into the mixing tip.

\*\* Due to the higher temperature of the Penta Putty paste extruded from the Pentamix mixing unit the intraoral setting of the wash material is accelerated.

These are the processing times for products that have been kept at a temperature of 23° C/74° F and 20–80% relative humidity. These times are shortened due to higher temperatures of the product or prolonged if the temperature of the product is lower.

The working time and intraoral setting time are indicated on the cartridge label of the product accompanied by the respective symbol:

Symbol for <b>maximum working time</b> for tray impression materials	
Symbol for <b>maximum working time</b> for syringing wash impression materials	
Symbol for <b>minimum intraoral setting time</b> (residence time in the mouth)	

**Impression Taking**

**One-Step Technique using Penta Material:**

Recommended product combinations	
Imprint 3 Penta Putty	Imprint 3 Ultra-Regular Body or Imprint 3 Regular Body
Imprint 3 Penta Quick Step Heavy Body	Imprint 3 Quick Step Light Body or Imprint 3 Quick Step Regular Body
Imprint 3 Penta Heavy Body	Imprint 3 Light Body or Imprint 3 Regular Body

Symbol for <b>maximum working time</b> for tray impression materials	
Symbol for <b>maximum working time</b> for syringing wash impression materials	
Symbol for <b>minimum intraoral setting time</b> (residence time in the mouth)	

Recommended product combinations	
Imprint 3 Penta Putty	Imprint 3 Ultra-Regular Body or Imprint 3 Regular Body
Imprint 3 Penta Quick Step Heavy Body	Imprint 3 Quick Step Light Body or Imprint 3 Quick Step Regular Body
Imprint 3 Penta Heavy Body	Imprint 3 Light Body or Imprint 3 Regular Body

- At the start of mixing the Penta impression material, be sure to observe the entry of the paste into the mixing tip, and ensure that base paste and catalyst paste flow into the mixing tip evenly. The mixed paste must be extruded in a uniform color free of streaks.
  - Please see “Product Description” for paste color.
- Load the impression tray prepared with adhesive, keeping the mixing tip constantly immersed in the paste to avoid air bubble formation.
- While the tray is being loaded, apply the wash material from the bottom up around the cleaned and dried preparation. Depending on the number of units to be syringed, initiate application, so the tray loading and application of material around the preparation are completed at the same time.
  - Be careful not to exceed the working time and intraoral syringe time of the respective material to avoid incomplete or inaccurate impressions.
- Slowly place the loaded tray in the mouth parallel to the long axes of the prepared teeth, and hold it in position passively.
- Allow the impression material to fully set, then remove the impression from the mouth.

Recommended product combinations	
Imprint 3 Quick Step Heavy Body	Imprint 3 Quick Step Light Body or Imprint 3 Quick Step Regular Body
Imprint 3 Heavy Body	Imprint 3 Light Body or Imprint 3 Regular Body

- Load the impression tray prepared with adhesive with the tray material, keeping the Garant Mixing Tip constantly immersed in the paste to avoid air bubble formation.
- While the tray is being loaded, apply the wash material from the bottom up around the cleaned and dried preparation. Depending on the number of units to be syringed, initiate application, so the tray loading and application of material around the preparation are completed at the same time.
  - Be careful not to exceed the working time and intraoral syringe time of the respective material to avoid incomplete or inaccurate impressions.
- Slowly place the loaded tray in the mouth parallel to the long axes of the prepared teeth, and hold it in position passively.
- Allow the impression material to fully set, then remove the impression from the mouth.

**After Impression Taking**

- Thoroughly examine and explore the sulcus of the prepared teeth and surrounding dentition. Remove any residual cured impression material from the mouth.

**Hygiene**

- Place the impression in a standard disinfectant solution for the period of time recommended by the manufacturer. Do not leave the impression submerged in the disinfection solution for more than the specified amount of time as otherwise the impression may be damaged.
- After disinfection, clean the impression of any remaining disinfection solution by rinsing under running water for approx. 15 seconds.

**Model Preparation**

- If Imprint 3 Penta Putty or Imprint 3 Quick Step Heavy Body have been used to take an impression, the impression can be poured at the earliest two hours after removal from the mouth. Impressions of all other product combinatios can be poured 30 minutes after removal from the mouth at the earliest. There are no other time limitations for pouring.
- Pre-rinse the impression with water and remove excess standing water in order to obtain a bubble-free model. A silicone surfactant can be used before pouring.
- Prepare the cast for the impression using a conventional die or model stone. The impressions can be silver or copper plated.

**Cleaning**

- Garant Dispenser: Remove paste that has not set using an alcohol-soaked cloth. The dispenser handle and plunger can be autoclaved up to a temperature of 135° C/275° F. Disassemble the dispenser before autoclaving. Glutaraldehyde-based solutions can be used for disinfection.
- Impression Tray: The adhesive can be removed from re-usable impression trays using acetone.

**Notes**

- The oxygen inhibition layer of composite materials, e.g. fillings or core buildups, may impair or retard the setting of silicone impression materials, and should be removed completely.
- Keep the filled Penta or Garant Mixing Tip on top of the cartridge to serve as a closure until the next use.** Removing and replacing the spent Penta or Garant Mixing Tip may lead to carry-over of paste and ensuing formation of clogs.
- Vinyl polysiloxane impression materials must not be combined with condensation polymerization silicones or polyether impression materials, trace amounts will impair the setting process.
- Disposable latex gloves impair the setting of silicone impression materials. Vinyl gloves are more suitable.
- Impressions should not come into contact with solvent-containing liquids, as this may result in swelling and imprecise modelling.

**Storage and Stability**
Store the product at 15–25° C/59–77° F. Do not use after the expiration date. Store set impressions at temperatures below 30° C/86° F.

**Customer Information**
No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

**Warranty**
3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user’s application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH’s sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

**Limitation of Liability**
Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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