3M Ketac™ Cem Easymix

ESPE

-  **Glass Ionomer Luting Cement – Radiopaque**
-  **Glassionomer-Befestigungszement – Röntgensichtbar**
-  **Ciment verre ionomère pour scellement – Radiopaque**
-  **Cemento vetro ionomero per fissaggio – Radiopaque**
-  **Cemento de vidrio ionómero para fijación – Radiopaque**
-  **Cimento de vidro ionómero para fixação – Radiopaque**

Instructions for Use / Informationen für den Gebrauch / Instructions for Use / Informações de uso

3M Deutschland GmbH
Health Care Business
Carl-Schurz-Str. 1
41453 Neuss – Germany
3M.com

3100019065/01

3M, ESPE and Ketac are trademarks of 3M or 3M Deutschland GmbH.
Unauthorized use prohibited. All rights reserved.

ENGLISH

Product Description
Ketac™ Cem Easymix is a luting cement for mixing by hand.
These instructions for Use should be kept for the duration of product use.
The product may only be used when the product labeling is clearly readable.

Intended Purpose
Intended purpose: dental luting agent for cementation of indirect restorations and orthodontic bands.
Intended users: educated dental professionals who have theoretical and practical knowledge on use of dental products.
Clinical benefit: holds indirect restorations and orthodontic bands in place in the context of prosthetic or orthodontic therapy.
Patient target group: all patients requiring dental treatment, unless the patient's condition limits the use.

- Indications**
- Cementation of inlays, onlays, crowns, and bridges made of metal, with metal frameworks, made of high-strength ceramic or composite, provided they are suitable for conventional cementation
 - Cementation of crowns and bridges made of metal, with metal frameworks, made of high-strength ceramic or composite, provided they are suitable for conventional cementation, on implant abutments
 - Cementation of posts and screws made of metal or high-strength ceramic, provided they are suitable for conventional cementation
 - Cementation of orthodontic bands.

Contraindications
None

Composition
The ingredients of each component of the product are indicated in decreasing order of their concentration.
Liquid: water, tartaric acid, powder: glass powder, polymeric acid.
For more information please consult the Safety Data Sheet/Safety Information Sheet (www.3m.com) or contact your 3M subsidiary.

Precautionary Measures
Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

Potential Undesirable Side Effects and Complications
Potential effects are local irritations. Irritation loss can result in loss of breakage of the restoration, which can then lead to swallowing or aspiration of the restoration or parts of it in rare cases.
The summary of safety and clinical performance for the device can be retrieved in the European database for medical devices (EuDamEd) under the following link (as soon as the corresponding database module is available):
<https://ec.europa.eu/tools/eudamed>
The instructions for Use and the summary of safety and clinical performance are also available at: <https://hcbregulatory.3m.com>
Basic UDI-DI: 06082232761020000000044W

- Preparation**
- Provide for optimal adhesion, carefully clean the surfaces of the tooth and restoration to be cemented.
 - Blow the preparation dry in only 2–3 short intervals with water and oil-free air or dab dry with cotton pellets. Do not dry over-dry!
 - **The preparation should be as dry as dry that the surface has a matt shiny appearance.** Excessive drying with compressed air can negatively affect adhesion or result in post-operative sensitivity after seating the restoration.

► The surfaces of the restoration should be dried well with compressed air.
Avoid further contamination by saliva.
Pulp Protection
The glass ionomer cement must not be applied directly to dentin – e.g. on inlays – situated close to the pulp or to the exposed pulp. Prior to taking the impression, cover areas close to the pulp using a hard-setting calcium hydroxide preparation.

- Dosage**
- Shake bottle to loosen powder.
 - Screw off cap and replace by dispensing dropper.
 - For cementing work, the standard powder-to-liquid mixing ratio is 3.8 : 1 w/w, this corresponding to 1 level spoonful of powder to two drops of liquid.
 - Skim off spoon on plastic seaghter. Do not compress the powder.
 - Dose powder and liquid next to each other on a mixing pad or glass slab.
 - Keep bottle upright while dispensing.
 - Crystallized liquid must not be allowed to form on the dispensing dropper.
 - Carefully seal bottles after use.
 - If the mix is too thin (due to insufficient dosing of powder) this may produce symptoms of pulpitis.

- Mixing**
Working of Ketac Cem Easymix should proceed at a room temperature of 20–25°C.
- Mix using a metal or plastic cement spatula.
 - Add the powder to the liquid in one portion.
 - On principle a sufficient quantity of powder should be added to the liquid to produce viscous consistency.
 - The quantity of cement required for a crown should just about drip from the spatula.
 - Continue to smooth out the paste until a homogeneous mixture is obtained.

Application

- During application, water and saliva must be prevented from entering the working area.
- Apply a thin coat of cement to the surface of the restoration to be cemented and to the preparation.
- Subsequently seat the restoration.
- Avoid overfilling the restoration. Applying hydrostatic pressure buildup in tightly fitting restorations may cause pulp discomfort.

Times
The following time scale applies at a room temperature of 23°C and 50% relative atmospheric humidity:

Mixing	min:sec	0:30
Application incl. mixing	3:10	
Setting from begin of mix	7:00	
Higher temperatures shorten the working time, whereas lower temperatures will extend this period (e.g. mixing on a cooled glass slab). A higher amount of powder also brings about a shorter working time. Exceeding the working time causes loss of adhesion to enamel and dentine.		

Removal of Excess

- Remove excess using a Heidemann spatula and/or probe 6–8 min. after beginning of mix.

After Cementing

- Carefully check the sulcus of the treated teeth and the surrounding areas; remove any cement or the restoration if still remaining. This is especially important when cementing to implant abutments.

Cleaning and Disinfection

Caution
Make sure that the cleaning and disinfecting agents you have chosen do not contain any of the following materials:

- Oxidizing agents (e.g., hydrogen peroxides)
- Chlorine
- Glutaraldehydes

Use only validated methods for cleaning and disinfection.
Use safety cleaning and disinfecting agents tested for efficacy and compatibility with the medical device (e.g., FDA certification). CavliWipes™ have been tested for use as a cleaning and disinfecting agent. Always observe all applicable legal and hygiene regulations for dental offices and/or hospitals.

General
The following procedure should be applied after use or before any subsequent use. The bottle must first be cleaned manually and, after that, be manually disinfected in accordance with the following instructions.

Manual Cleaning by Wiping

- Clean the bottle after or before every application with ready-to-use disinfectant wipes (e.g., CavliWipes™ based on the active substances alcohol and quaternary compounds) or equivalent cleaning agents. Follow the cleaning agent manufacturer's directions for use.
- Thoroughly clean the surface of the bottle with a new ready-to-use CavliWipes™ for at least 30 seconds or until no more dirt is visible.
- Special care must be exercised in cleaning of the grooves on the lid.

Manual Disinfection by Wiping (Intermediate Level Disinfection)

- Disinfect the bottle with ready-to-use disinfectant wipes (e.g., CavliWipes™ based on the active substances alcohol and quaternary compounds) or equivalent disinfecting agents. Follow the disinfecting agent manufacturer's directions for use, observing in particular the contact times.
- Disinfect the surface of the bottle for three minutes using at least two new disinfectant wipes. Keep all surfaces moist during the three-minute period. Special care must be exercised in disinfection of the grooves on the lid.

Inspection, Maintenance, Testing

- Examine the bottle for damage, discoloration, and contamination before every use.
- Do not use damaged devices under any circumstances. If the bottle is not visibly clean, repeat the cleaning and disinfection procedure.

Notes

Ketac Cem Easymix is a multiple patient multiple use device.





Disposal
Dispose of the contents or container in accordance with the applicable regulations. Please pay special attention to the disposal of contaminated waste to avoid health risks because of improper handling.

Storage and Shelf Life
Keep the powder in a dry place away from damp.
Carefully reseal the powder bottle after each use.
Store the product at 15–25°C/59–77°F.
Do not use after the expiry date.

Customer Information
No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty
3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

Limitation of Liability
Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Reference Number and Symbol Title	Symbol	Description of Symbol
ISO 15223-1 5.1.1 Manufacturer		Indicates the medical device manufacturer.
Swiss Authorised Representative		Indicates the authorised representative in Switzerland.
ISO 15223-1 5.1.3 Date of Manufacture		Indicates the date when the medical device was manufactured.
ISO 15223-1 5.1.4 Use-by date		Indicates the date after which the medical device is not to be used.
ISO 15223-1 5.1.5 Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified.
ISO 15223-1 5.1.6 Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified.
ISO 15223-1 5.3.7 Temperature limit		Indicates the temperature limits to which the medical device can be safely exposed.
ISO 15223-1 5.4.3 Consult instructions for use or consult electronic instructions for use		Indicates the need for the user to consult the instructions for use.
CE Mark		Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement.
Medical Device		Indicates the item is a medical device.
Non-corrugated fiberboard		Indicates product packaging is made of non-corrugated fiberboard. Official Journal of the EC, Commission Decision (97/129/EC)
Colorless Glass		Indicates product packaging is made of colorless glass. Official Journal of the EC, Commission Decision (97/129/EC)
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.
ISO 15223-1 5.1.8 Importer		Indicates the entity importing the medical device into the locale.

For more information, see HCBRegulatory.3M.com

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Καυτόόηοηετ
Brugsanvisning
Bruksveiledning

			
3M Deutschland GmbH <div>Health Care Business Carl-Schurz-Str. 1 41453 Nüssa - Germany 3M.com</div>			
3M, ESPE and Ketac are trademarks of 3M or 3M Deutschland GmbH. © 2021, 3M. Unauthorised use prohibited. All rights reserved.			

			
3M Deutschland GmbH <div>Health Care Business Carl-Schurz-Str. 1 41453 Nüssa - Germany 3M.com</div>			
3M, ESPE and Ketac are trademarks of 3M or 3M Deutschland GmbH. © 2021, 3M. Unauthorised use prohibited. All rights reserved.			

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