3M Ketac™ Cem Apilcap™

ESPE

- Glass Ionomer Luting Cement – Radiopaque**
- Glasionomer-Befestigungszement – Röntgensichtbar**
- Ciment verre ionomère pour scellement – Radiopaque**
- Cemento vetroso ionomero per fissaggio – Radiopaque**
- Cemento de vidrio ionómero para fijación – Radiopaque**
- Cimento de vidro ionómero para fixação – Radiopaque**

Instructions for Use
Gebrauchsinformation
Notice d'utilisation

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3M, ESPE, Apilcap, CapMix, Ketac, Maxicap and RotoMix are trademarks of 3M or 3M Deutschland GmbH. Used and/or reprinted in Canada.
© 2021, 3M. Unauthorized use prohibited. All rights reserved.

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ENGLISH

Product Description
Ketac™ Cem Apilcap™ (yellow, similar to dentin) and Ketac™ Cem Maxicap™ (light pink) are luting cements in capsule form. The minimum amount being dispensed from an Apilcap capsule is 0.07 ml; the minimum amount from a Maxicap capsule is 0.36 ml.

These Instructions for Use should be kept for the duration of product use. The product may only be used when the product labeling is clearly readable. For details on all additionally mentioned products please refer to the corresponding instructions for Use.

Intended Purpose
Intended purpose: dental luting agent for cementation of indirect restorations and orthodontic bands.
Intended users: educated dental professionals who have theoretical and practical knowledge on usage of dental products.
Clinical benefit: holds indirect restorations and orthodontic bands in place in the context of prosthetic or orthodontic therapy.
Patient target group: all patients requiring dental treatment, unless the patient's condition limits the use.

- Indications**
- Cementation of inlays, onlays, crowns, and bridges made of metal, with metal frameworks, made of high-strength ceramic or composite, provided they are suitable for conventional cementation.
 - Cementation of crowns and bridges made of metal, with metal frameworks, made of high-strength ceramic or composite, provided they are suitable for conventional cementation, on implant abutments
 - Cementation of posts and screws made of metal or high-strength ceramic, provided they are suitable for conventional cementation
 - Cementation of orthodontic bands

Contraindications
None.

Composition
The ingredients of each component of the product are indicated in decreasing order of their concentration.
Liquid: water, tartaric acid, Powder: glass powder.
For more information please consult the Safety Data Sheet/Safety Information Sheet (www.3m.com) or contact your 3M subsidiary.

Precautionary Measures
Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

Potential Undesirable Side Effects and Complications
Potential effects are local irritations. Retention loss can result in loss or breakage of the restoration, which can then lead to swallowing or aspiration of the restoration or parts of it in rare cases.
The summary of safety and clinical performance for the device can be retrieved in the European database for medical devices (Edumed) under the following link (as soon as the corresponding database module is available):
<https://ec.europa.eu/tools/edumed>

The instructions for Use and the summary of safety and clinical performance are also available at: <https://hcbregulatory.3m.com>
Basic UDI-DI: 0608223276102000000044DW

Preparation
► To provide for optimal adhesion, carefully clean the surfaces of the tooth and restoration to be cemented.
► Blow the preparation dry in only 2–3 short intervals with water and oil-free air or dab dry with cotton pellets. Do not over-dry!
– The preparation should be just dry so that the surface has a matt shiny appearance. Excessive drying with compressed air can negatively affect adhesion or result in post-operative sensitivity after seating the restoration.

► The surfaces of the restoration should be dried well with compressed air.
► Avoid further contamination by saliva.

Pulp Protection
The glass ionomer cement must not be applied directly to dentin – e.g. on inlays – situated close to the pulp or to the exposed pulp.
► Prior to taking the impression, cover areas close to the pulp using a hard-setting calcium hydroxide material.

Activation
► Place the Apilcap or Maxicap Activator on a sturdy surface and insert the Apilcap or Maxicap capsule in the activator.
► Blow the ball of your hands, depress the activator lever firmly and completely to the stop and hold it down for 2 to 4 seconds.
– Firmly and completely depressing the lever to the stop and holding it in this position is the only way to ensure that all of the liquid is pressed or mixed into the powder. Inadequate pressure is applied or if the lever is not depressed completely, it is possible that too little liquid will be mixed with the powder, increasing the viscosity and changing the product characteristics.

Mixing
► Mix the capsules at approx. 4,300 rpm, in a high frequency mixing device such as CapMix™ for 10 sec or the RotoMix™ rotating mixer (for 8 sec), and apply with the Apilcap or Maxicap Applier. Both the Apilcap Activator and the Maxicap Applier are color-coded in orange; both the Maxicap Activator and the Maxicap Applier are color-coded in blue.

Application
► After mixing, insert the capsule into the applier and secure it in place by pressing the lever once. Before applying the paste intracavity, open the nozzle to the stop and press the applier several times until the mixed paste becomes visible in the intra-pulp.
► Ensure that there is no contamination from water, saliva, etc., entering the working area for the entire duration of application; the working area must be kept dry.
► Apply a thin coat of cement to the inside of the restoration and to the preparation. Subsequently seat the restoration. Avoid overfilling the restoration.

Attention
Capsules must be mixed immediately after activation. Pauses between mixing and application must be avoided; otherwise, as soon as the material starts to set, removal from the capsule will be rendered difficult or even impossible. Hydrostatic pressure buildup in tightly fitting restorations may cause pulp discomfort.

Times
The following time scale applies at a room temperature of 23°C/73°F:

	min:sec
Activation	0:08
Mixing in RotoMix	0:10
Mixing in high-speed mixer, e.g. CapMix	0:10
Working (from beginning of mix)	3:00
Setting (from beginning of mix)	7:00

The setting time is accelerated at room temperatures above 23°C/73°F. Exceeding the working time will cause loss of adhesion to enamel and dentin. Shorter mixing times will prolong the working and setting times, whereas longer mixing times will shorten them. On high-speed mixers the mixing time can be set from 8–15 sec and on the RotoMix from 7–12 sec.

Removal of Excess
► Remove excess using a Heidemann spatula and/or probe 7 min after beginning of mix.

After Cementing
► Carefully check the sulcus of the treated teeth and the surrounding areas; remove any cement residues still remaining. This is especially important when cementing to implant abutments.

Cleaning the Apilcap and Maxicap Activator and Applier
Please refer to the pertinent Instructions for Use for information about cleaning and disinfecting the Apilcap/Maxicap Activator and Applier.

Notes
Apilcap and Maxicap capsules are single-use devices. Please keep the single-use capsule in the original blister to guarantee the traceability of the batch (LOT).

Disposal
Dispose of the contents or container in accordance with the applicable regulations. Please pay special attention to the disposal of contaminated waste to avoid health risks because of improper handling.

Storage and Shelf Life
Store the product at 15–25°C/59–77°F. Do not use after the expiry date.

Customer Information
No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty
3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy is that 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Symbol Glossary

Reference Number and Symbol	Symbol	Description of Symbol
ISO 15223-1 5.1.1 Manufacturer		Indicates the medical device manufacturer.
ISO 15223-1 5.1.3 Swiss Authorized Representative		Indicates the authorised representative in Switzerland.
ISO 15223-1 5.1.3 Date of Manufacture		Indicates the date when the medical device was manufactured.
ISO 15223-1 5.1.4 Use-by date		Indicates the date after which the medical device is not to be used.
ISO 15223-1 5.1.5 Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified.
ISO 15223-1 5.1.6 Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified.
ISO 15223-1 5.3.7 Temperature limit		Indicates the temperature limits to which the medical device can be safely exposed.
ISO 15223-1 5.4.2 Do not re-use		Indicates a medical device that is intended for one single use only.
ISO 15223-1 5.4.3 Consult instructions for use or consult electronic instructions for use		Indicates the need for the user to consult the instructions for use.
CE Mark		Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement.
Medical Device		Indicates the item is a medical device.
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a dental professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1).
Non-corrugated fiberboard		Indicates product packaging is made of non-corrugated fiberboard. Official Journal of the EC, Commission Decision (97/129/EC)
Paper		Indicates product packaging is made of paper (blister). Official Journal of the EC, Commission Decision (97/129/EC)
Low-density polyethylene		Indicates plastic component made of low-density polyethylene is recyclable (invisibly marked with recycling symbol).
Other resins		Indicates product packaging is made of other plastic component (blister). ASTM International (American Society for Testing and Materials) 7611/D 7611M
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.
ISO 15223-1 5.1.8 Importer		Indicates the entity importing the medical device into the locale.

For more information, see [HCBRegulatory.3m.com](https://hcbregulatory.3m.com)

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