3M RelyX™ Unicem 2 Automix

Self-adhesive Resin Cement

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Used under license in Canada. © 2022, 3M. Unauthorized use prohibited. All rights reserved.

Do not use substances such as desensitizers, disinfectants, stringents, dentin sealants, rinsing solutions containing EDTA, scanning powder, etc., after the final curing of water. Their residues may have a detrimental effect on the bonding strength and setting reaction of the cement.

Preparation of Root Canals

- Remove the root canal filling, leaving at least 4 mm of root canal filling (guttaepatch) apically.
- Prepare the root canal to the post.
- Clean the root canal with sodium hypochlorite solution (NaOCl solution).
- Immediately rinse with water, and dry with paper tips.
- The use of a rubber dam during cementation of posts is recommended.

Preparation for Maryland and Inlay/Onlay Bridges

Maryland and inlay/onlay bridges are not recommended for patients with bruxism or periodontitis. Adjacent teeth must have an adequate enamel surface for bonding. They should be healthy or only slightly restored and the periodontal conditions should be good. It is the sole responsibility of the dentist to ensure proper selection of indication and technique. The guidelines of the relevant national professional associations must be observed for such indications.

- Prepare retentive elements such as angular rests and/or approximal grooves.
- Use a rubber dam and keep the cavity free of any contamination during cementation.
- Etch the enamel surface of the cavity with 37% phosphoric acid for 15–20 sec. Then rinse thoroughly with water and dry with water-free and oil-free air. In the case of exposed dentin, make sure to selectively etch the enamel to avoid post-operative sensitivity.

Preparatory Measures

- Remove the temporary restoration and thoroughly clean the tooth surface of any residual cement or other contaminants.
- Make a trial insertion of the final restoration and check its fit and contact points.
 - If a low viscosity silicone is used for the trial fit, the silicone residue must be cleaned by the dentist with atmospheric oxygen during the polymerization process.
 - For glass ceramic restorations, do not check the occlusion until the restoration has been cemented. Breakage could occur if done prior to being cemented in place.
- Avoid any contaminant whatsoever of the surfaces to be treated during pre-treatment and until the final cementation.

Pre-treatment of Restorations

Pre-treatment of Metal Restorations Please follow the instructions for use of the restoration material. In the absence of deviating instructions, we recommend the following procedure:

- Sandblast surfaces to be cemented. Use 30 or 50 µm aluminum oxide at a pressure of 2 bar (30 psi) to create a roughened matte surface appearance.
- Clean the blasted surface with alcohol and dry it with water-free and oil-free air.
- Surfaces must be cleaned after try-in (sodium hypochlorite (NaOCl) solution is recommended), rinsed with water, and dried.

For Maryland and inlay/onlay bridges, please refer to "Pre-treatment of Maryland and Inlay/Onlay Bridges".

Pre-treatment of Etchable Glass Ceramic Restorations Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:

- Use hydrofluoric acid to etch the glass ceramic restoration surface to be cemented.
- Rinse thoroughly with water for 15 seconds and dry with water-free and oil-free air.

For Maryland and inlay/onlay bridges, please refer to "Pre-treatment of Maryland and Inlay/Onlay Bridges".

Pre-treatment of Zirconia and Aluminum Oxide Ceramic Restorations Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:

Alternative 1:

- Sandblast surfaces to be cemented. Use 30 or 50 µm aluminum oxide at a pressure of 2 bar (30 psi) to create a roughened matte surface appearance.
- Clean the blasted surface with alcohol and dry it with water-free and oil-free air.

Alternative 2:

- Coat (silicate) the restoration surface to be luted with the micro-blasting device CoJet™ Prap and the blast-coating agent CoJet™ Sand from a distance of 30 mm and vertically to the surface for 15 sec, see instructions for use for CoJet Prap and CoJet Sand, respectively.
- Blow away any residues of the blasting agent with water-free and oil-free air.
- Apply an appropriate silane in accordance with the instructions for use.

For Maryland and inlay/onlay bridges, please refer to "Pre-treatment of Maryland and Inlay/Onlay Bridges".

Pre-treatment of Composite Restorations Please follow the instructions for use of the restoration material. If the manufacturer has not provided deviating instructions, we recommend the following procedure:

- Sandblast surfaces to be cemented. Use 30 or 50 µm aluminum oxide at a pressure of 2 bar (30 psi) to create a roughened matte surface appearance.
- Clean the blasted surface with alcohol and dry it with water-free and oil-free air.
- Surfaces must be cleaned after try-in (sodium hypochlorite (NaOCl) solution is recommended), rinsed with water, and dried.

Pre-treatment of Glass Fiber-reinforced Posts

Please follow the instructions for use of the post to be used. If the manufacturer has not provided deviating instructions, we recommend the following procedure:

- Clean the post with alcohol and dry it with water-free and oil-free air.
- When using RelyX™ Fiber Post application of a silane is not required. Other fiber-reinforced posts should be treated in accordance with the relevant instructions for use.

Pre-treatment of Implant Abutments

Please follow the recommendations for the relevant restorative material (e.g., metal, zirconia).

Pre-treatment of Maryland and Inlay/Onlay Bridges

Metal, Zirconia and Aluminum Oxide Surfaces

- Coat (silicate) the restoration surface to be luted with the micro-blasting device CoJet™ Prap and the blast-coating agent CoJet™ Sand from a distance of 30 mm and vertically to the surface for 15 sec, see instructions for use for CoJet Prap and CoJet Sand, respectively.
- Blow away any residues of the blasting agent with water-free and oil-free air.
- Apply an appropriate silane in accordance with the instructions for use.

For etchable glass ceramics, please refer to "Pre-treatment of Etchable Glass Ceramic Restorations".

The processing and setting times depend on the ambient and oral temperature. The times shown are based on conditions relevant for practice. As is the case with every composite cement, the setting of RelyX Unicem 2 Automix cement slows down significantly at room temperature. RelyX Unicem 2 Automix is not so light sensitive and therefore also sensitive to natural or artificial light. The working time is significantly reduced during application under operating lights!

	minutes
Working time from start of mixing	02:30
	seconds
Light-curing in a lightwave range of 400–500 nm	
– Single surface, from occlusal	00:20
– Any other surface additional	00:40
– RelyX Fiber Post posts, from occlusal	00:40
Self-curing	
– Setting time after start of mixing	06:00

Application

- Select a RelyX Unicem 2 Automix syringe containing the desired color and place it ready for use.

When using a new RelyX Unicem 2 Automix syringe:

- Remove the automix syringe from the foil bag and discard the bag. Note down the date and the lot number on the syringe label.
- Remove and discard the sealing cap from the automix syringe. Check the syringe openings for blockage; remove any paste plugs.

The **sealing cap must not be used to close the syringe** because doing so might introduce air bubbles into the material. Squeeze out a small quantity of paste onto a block to equal the base paste and calibrate the syringe in the automix syringe. Discard the paste which has been squeezed out.

During every application:

Protect the working area from contamination with water, blood, saliva and sputum by using a rubber dam or other suitable barrier.

- Remove and discard the mixing tip remaining on the automix syringe from the previous application.

- Attach one of the two types of mixing tips ("Standard" or "Wide") and secure it by turning it to the right.
- The mixing tip "Standard" is used throughout an Intraoral Tip for application in easily accessible areas.

- The mixing tip "Wide" is used for the application with the Intraoral Tip or the Endo Tip.

When using the mixing tip "Wide", attach an Intraoral Tip or an Endo Tip.

- Endo Tips may not be shortened due to the risk of deformation or blockage from shavings.

- Squeeze out and discard a peppercorn-size quantity of RelyX Unicem 2 Automix until an evenly mixed paste in a homogeneous color flows out of the tip.
 - The paste requires a certain amount of time to flow through the mixing tip. The flowthrough speed cannot be accelerated by increasing the pressure on the plunger.
 - As a result of the pressure on the plunger decreases, the material flow stops and the paste begins to set up. Do not use force to press out paste that has set as this could cause damage to the mixing tip and the automix syringe.
- Apply RelyX Unicem 2 Automix evenly to the entire cavity and as appropriate to the bottom side of the inlay/onlay, or fill the crown with cement.
 - Keep the opening of the mixing tip "Standard" or of the Intraoral or Endo Tip immersed in the material during the entire application to prevent the inclusion of any air bubbles.
- Seat the restoration firmly and stabilize long enough for the cement to set fully.
- Leave the used mixing tip on the automix syringe as a cap until the next application.

Application in the Root Canal See instructions under "Application" for handling of the automix syringe.

- Do not use Lento-Loops to insert the cement in the root canal as this can excessively accelerate setting.
- Attach an Endo Tip to the mixing tip "Wide" for application in the root canal.

Insert the Endo Tip as deeply as possible in the root canal and apply RelyX Unicem 2 Automix, beginning apically. Keep the tip of the Endo Tip immersed in the cement and slowly move the Endo Tip upwards as the level of the paste rises.

- Do not remove the Endo Tip from the cement until the root canal has been completely filled.
 - This so-called minimum filling should not be carried out in less than 5 sec.; this will minimize the entrapment of air bubbles.
- Place the post in the root canal filled with cement; apply moderate pressure to hold it in position. We recommend rotating the post slightly during insertion to avoid the inclusion of air bubbles.

Removal of Excess

Sores on oxygen inhibition: As is the case with all composite materials, RelyX Unicem 2 Automix is subject to oxygen inhibition, i.e., the upper layer (about 50 µm) which is in contact with atmospheric oxygen during the polymerization does not harden. Leaving sufficient excess permits the removal of the uncured layer during shaping/polishing without leaving behind a deficit.

From Restorations

Using light-curing:

- After brief curing (about 2 sec using a standard polymerization device), use a suitable instrument (e.g., scaler) to remove the excess. Larger volume excess can be removed more easily.

Immediately light-cure after removal [see Polymerization and Shaping].

Using self-curing:

- Use a suitable instrument (e.g., scaler) to remove the excess during the self-curing phase (about 3 min after beginning the mixing in the "gel phase").
- Immediately cover the edges of the restoration with a glycerin gel to prevent oxygen inhibition.

Alternatively:

- Remove the excess before polymerization, e.g., using a sponge pellet, while using a suitable instrument to hold the restoration in position.

- Cover the edges of the restoration with a glycerin gel to prevent oxygen inhibition.
- Immediately light-cure or wait for the self-curing to finish [see Polymerization and Shaping].

From Posts

- Remove the cement with a suitable instrument or a cotton pellet.

Polymerization and Shaping

We recommend light-curing the cement through the restoration when doing ceramic and composite work. Select the exposure times appropriate for the number of surfaces (please refer to "Times"). Polymerize the cement through the post when using translucent posts. The exposure time depends on the translucency of the post being used; for RelyX Fiber, it is 40 sec.

- Shape any remaining uneven areas on the edges of restorations, and polish each area with a suitable device, aluminum oxide coated disc (e.g., Sof-Lex™), and diamond polishing paste.
- Then check the occlusion.

After Shaping

- Carefully check the gutlets of the treated teeth and the surrounding areas; remove any cement residues still remaining. This is especially important when cementing on implant abutments.

Hygiene and Disinfection

Handling and Cleaning Before use of the Automix syringe, put on new, unused disposable gloves during all clinical applications of the Automix syringe to reduce the risk of cross-contamination.

For hygienic reasons, use a new mixing tip and protective sleeve for every patient to avoid risk of infection.

Place the protective sleeve over the Automix syringe with the previously attached mixing tip and use the tip of the mixing tip to carefully push through the end of the protective sleeve.

Avoid contact between the reusable parts (e.g., the body of the syringe) and the patient's mouth.

A mixing tip (disposable product) is used to apply the product to the patient and/or restoration; it is discarded after being used once on the patient.

After using the Automix syringe in the protective sleeve, remove the sleeve carefully so that there is no contamination from the outer surface of the protective sleeve. Discard the used protective sleeve.

After removal of the protective sleeve, wipe the Automix syringe thoroughly for 30 seconds with a ready-to-use CavWipe™ or an equivalent agent.

Leave the used mixing tip on the Automix syringe as a cap until the next application.

Always observe all applicable legal and hygiene regulations for dental offices and/or hospitals during use and reprocessing.

Note: Dispose the Automix syringe if contaminated, e.g., with blood, saliva, spray of body fluids, or contaminated hands. Do not place the Automix syringe into a disinfection bath or a cleaning-disinfection device (washer-disinfector).

Inspection, Maintenance, Testing

Examine the material for damage, discoloration, and contamination before every use. Do not use damaged and contaminated syringes under any circumstances.

Notes

- Automix Mixing Tips Regular and Wide, Intraoral Tips and Endo Tips are single-use devices and would be damaged when removing the set paste. When reattaching Automix Tips the proper fit to the tip on the Automix cartridge cannot be ensured.
- Examine the Automix Mixing Tips Regular and Wide, Intraoral Tips and Endo Tips for damage, discoloration, and contamination before every use. Do not use damaged, discolored, or contaminated tips under any circumstances. Attach the Mixing Tips Regular or Wide, Intraoral Tips or Endo Tips immediately before use and do not detach and reattach to avoid accidental detachment.
- Please keep the single-use accessory in the original bag to guarantee the traceability of the batch.
- RelyX Unicem 2 Automix is a multiple patient multiple use device.
- RelyX Unicem 2 Automix can be prematurely polymerized when exposed to natural or artificial lighting. Therefore avoid intensive light exposure during application.

Disposal

Dispose of the contents or container in accordance with the applicable regulations to avoid health risks because of improper handling.

Storage and Stability

Store RelyX Unicem 2 Automix in the pouch. After removal from the pouch, use the product within 6 months and before the lapse of the expiration date. Constantly high humidity accelerates the setting and must be avoided. Always store at 15°–25 °C/59°–77 °F. Do not use after expiration date.

Customer Information No product information provides any information which deviates from the information provided in this instruction sheet.

Warnings

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH's product.

Limitation of Liability Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty contract, negligence, or strict liability.

Symbol Glossary

Reference and Symbol Title	Symbol	Description of Symbol
ISO 15223-1 5.1.1. Manufacturer		Indicates the medical device manufacturer.
Swiss Authorized Representative		Indicates the authorised representative in Switzerland.
ISO 15223-1 5.1.3. Date of Manufacture		Indicates the date when the medical device was manufactured.
ISO 15223-1 5.1.4. Use-by date		Indicates the date after which the medical device is not to be used.
ISO 15223-1 5.1.5. Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified.
ISO 15223-1 5.1.6. Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified.
ISO 15223-1 5.3.7. Temperature limit		Indicates the temperature limits to which the medical device can be safely exposed.
ISO 15223-1 5.4.2. Do not re-use		Indicates a medical device that is intended for one single use only.
ISO 15223-1 5.4.3. Consult instructions for use or consult electronic instructions for use		Indicates the need for the user to consult the instructions for use.
CE Mark		Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement.
ISO 15223-1 5.7.7.		

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3M, ESPE, CoJet, Lava, RelyX and Sof-Lex are trademarks of 3M or 3M Deutschland GmbH.
3M, ESPE, CoJet, Lava, RelyX und Sof-Lex are trademarks of 3M or 3M Deutschland GmbH.
3M, ESP, CoJet, Lava, RelyX and Sof-Lex are trademarks of 3M or 3M Deutschland GmbH.
Este uso no autorizado, suministrado en el nombre de 3M. No autorizado use prohibited. All rights reserved.