

RelyX™ Veneer Cement RelyX™ Try-In Paste Adper™ Single Bond Plus Adhesive

- Verify that the dental lab has etched the veneers.
- · Remove provisional restoration if present.

Preparation for Bonding:

- Clean the teeth using a plain flour of pumice slurry, rinse and lightly air dry.
- Place matrix strips between teeth to prevent unwanted bonding to adjacent teeth.
- Isolate the bonding area to prevent contamination.

Try-In/Shade Selection:

- Try-in veneer with water or Translucent shade of the watersoluble RelyX Try-In Paste to verify fit and esthetics.
- If esthetics and fit are acceptable, remove veneer, rinse tooth and veneer thoroughly, dry and proceed with bonding procedure.
- If color adjustment is needed, select appropriate shade of RelyX Try-In paste, apply to veneer, seat, and examine for color and fit.
- When the correct shade has been chosen, remove veneer, rinse tooth and veneer thoroughly, dry and proceed with bonding procedure.

Silane Prime Veneer:

- Clean the bonding surface of the veneer by applying Scotchbond[™] Etchant (35% phosphoric acid) for 15 seconds, rinse and dry.
- Apply a single coat of the RelyX[™] Ceramic Primer to the bonding surface of the veneer and dry.

Etchant:

- Apply Scotchbond[™] Etchant to enamel and dentin for 15 seconds.
- Rinse for 10 seconds.
- · Blot excess water leaving tooth moist.





RelyX[™] Veneer Cement RelyX[™] Try-In Paste Adper[™] Single Bond Plus Adhesive



Adhesive:

- Using a fully saturated brush tip for each coat, apply 2 consecutive coats of Adper Single Bond Plus Adhesive to enamel and dentin.
- Dry gently for 2-5 seconds.
- Do Not Light-cure!



Treat the Veneer:

- Apply 1 coat of adhesive to the silane treated bonding surface of the veneer.
- Dry gently for 2-5 seconds.
- Do Not Light-cure!
- Apply selected shade of RelyX Veneer Cement to the veneer



Seat Veneers:

- Seat the veneer with gentle pressure.
- Spot cure the veneer to secure in place by light-curing on the facial surface with a small diameter light guide for 20 seconds. Avoid curing excess cement.
- Clean excess cement from margins.
- Alternatively, the excess cement can be "tack-cured" for easier clean-up.
- Light-cure each area and margin of the veneer for 30 seconds (40 seconds for the A5 opaque/Dark shade).
- Remove matrix strips and finish restoration.
- See product instructions for detailed curing and finishing information.



