ENGLISH

General information:

3M™ ESPE™ Adper™ Multi-Purpose Plus System, is a versatile system for bonding all classes of direct composite restorations as well as for indirect procedures involving metal, porcelain, or composite crowns, inlays, and onlays. The Adper Scotchbond multi-purpose plus system also bonds amalgam and self cure composite and can be used to bond orthodontic brackets to crowns. Bonding light cure composite requires only the Adper Scotchbond multi-purpose etchant, primer and adhesive components. This system is recommended for all classes of restorations. Direct composite placement involves etching, priming, and application of a light cure adhesive. When bonding light cure composite to porcelain or alloys 3M™ ESPE™ RelyX™ ceramic primer is used instead of the Adper Scotchbond multi-purpse primer.

The addition of Adper Scotchbond multi-purpose plus activator and Adper Scotchbond multi-purpose plus catalyst extends the applications to include indirect bonding procedures and bonding amalgam and self cure composite. For indirect bonding procedures, oral surfaces are prepared and etched. Adper Scotchbond multi-purpose plus activator is applied and dried. Next, Adper Scotchbond multi-purpose primer is applied and dried. Lastly, Adper Scotchbond multi-purpose plus catalyst is applied. Using only the catalyst rather than a mixture of adhesive and catalyst ensures that the adhesive will not be setting while the indirect restoration is being prepared and placed. The indirect restorative is surface treated as appropriate for the material chosen, followed by application of the Adper Scotchbond multi-purpose plus catalyst and the dual cure or self cure resin-based luting material. Set time is dictated by the choice of resin-based luting material.

Note: When bonding indirect restorations, the adhesive components that have been placed on the tooth will not set until they come in contact with the mixed resin-based luting material. When the luting material on the indirect restoration contacts the treated tooth the time available for adjustment of the restoration will be reduced by both the elevated oral temperature and by interaction with the adhesive system components. For example, a 5 minute room temperature adjustment time can be reduced to approximately 2 minutes at body temperature, and further reduced to approximately 30 seconds by interaction with the adhesive system. This effect can be tested in the office using two microscope slides.

Treat one slide with Adper Scotchbond multi-purpose plus activator, primer, and catalyst to represent the preparation. Treat the other slide with Adper Scotchbond multi-purpose plus catalyst and the mixed luting cement to represent the indirect restoration. Place the two treated surfaces in contact and move the slides at 10 second intervals to measure the adjustment time available after placement. Elevating the temperature using an incubator or by holding the slides in gloved hands will add to the realism of the test. If a dual cure luting material is used, it is recommended that the margins be light cured to yield the highest physical properties to the luting material. Bonding amalgam or self cure composite begins with an etch of enamel and dentin. Adper Scotchbond multi-purpose plus activator is then applied, followed with an application of Adper Scotchbond multi-purpose primer. Adper Scotchbond multi-purpose adhesive and catalyst are mixed and applied, and the restorative material is placed. In this case adhesive and catalyst are mixed because no luting cement is present.

Indications

Bonding all classes of direct composite restorations as well as for indirect procedures involving metal, porcelain, or composite crowns, inlays, and onlays. The Adper Scotchbond multi-purpose plus system also bonds amalgam and self-cure composite and can be used to bond orthodontic brackets to crowns.

Recommendations:

Use 3M™ ESPE™ Vitrebond™ light cure glass ionomer liner/base in areas of deep cavity excavation such as Class I and II restorations. If a pulp exposure has occurred, use a minimum amount of calcium hydroxide followed by an application of Vitrebond liner/base.

Adper Scotchbond multi-purpose adhesive will bond to Vitrebond light cure glass ionomer liner/base whether or not the ionomer was treated with etchant and primer. If a temporary material was used on the Vitrebond glass ionomer liner/ base, an alcohol wipe is required after removal of the temporary cement.

The Adper Scotchbond multi-purpose system includes an etch of enamel and dentin as part of the procedure. It is recommended that the surfaces be left moist after the rinse. Excess surface moisture can be removed with a brief exposure to an air syringe or by blotting.

The Adper Scotchbond multi-purpose primer can be dried immediately without agitation or waiting time. Do not rinse off the primer. The primed surface should appear shiny. With gentle drying this will occur with only one application of the primer.

Light curing time assumes the use of a 3M™ ESPE™ light curing unit or other dental visible light curing unit of comparable intensity. Air used for drying the various components should be free of oil and water contaminants.

Components:

3M™ ESPE™ Scotchbond™ Universal Etchant contains 32% phosphoric acid by weight etches the enamel and removes the dentinal smear layer. Use of an etchant is critical on both enamel and dentinal surfaces. If etched surfaces are contaminated it is recommended that the surfaces be re-etched and rinsed.

Adper Scotchbond multi-purpose plus activator is needed only for self cure or dual cure applications such as bonding amalgam or indirect restorations. The activator reacts with the primer to facilitate curing.

Adper Scotchbond multi-purpse primer facilitates the wetting of the adhesive onto the prepared tooth structure. Use on dentinal surfaces is mandatory for optimum adhesion.

Adper Scotchbond multi-purpose adhesive is the light cure component of the system. It bonds to etched enamel and to dentin when conditioned using the etchant and primer. It will not self cure without the addition of Adper Scotchbond multi-purpose plus catalyst.

Adper Scotchbond multi-purpose plus catalyst, when combined with Adper Scotchbond multi-purpose adhesive, results in a dual cure resin suitable for bonding amalgam and self cure composite. The catalyst is also used alone in indirect bonding applications.

Precautionary Information for Patients

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of this product in patients with known acrylate allergies. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If allergic reaction occurs, seek medical attention as needed, remove the product if necessary and discontinue future use of the product.

Precautionary Information for Dental Personnel

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If product contacts glove, remove and discard glove, wash hands immediately with soap and water and then re-glove. If allergic reaction occurs, seek medical attention as needed.

Adper Scotchbond multi-purpose plus Activator and RelyX Ceramic Primer are highly flammable. Avoid static discharge. Keep away from heat, sparks, open flame, pilot light, and other sources of ignition.

Scotchbond Universal etchant: Contact with eyes or skin may cause burns. Wear suitable protective clothing, gloves and eye/face protection.In case of contact, flush area with plenty of water and seek medical advice.

3M ESPE MSDSs can be obtained from www.3MESPE.com or contact your local subsidiary.

Sensitivity

Some patients may experience transitory postoperative sensitivity. The risk of sensitivity can be minimized by the following measures: Remove minimal tooth structure. Use proper isolation. Use of a rubber dam is highly recommended. Use adequate pulp protection. Use a glass ionomer cavity liner/base in areas of deep

excavation. Place restorative material in increments, curing each increment separately. Adequately cure restorative according to instructions for shade and thickness of restorative and light exposure time. Adjust occlusion carefully.

Check for hyperocclusion, particularly in lateral excursion contacts.

Instructions for Use

Etchant syringe assembly

- Protective eyewear for patients and dental staff is recommended when using the delivery system.
- Prepare delivery system: Remove cap from etchant syringe and SAVE. Twist a blue disposable tip securely onto the syringe. Holding the tip away from the patient and any dental staff, express a small amount of etchant onto a dispensing pad or a 2 x 2 gauge to assure that the delivery system is not clogged. If clogged, remove the dispensing tip and express a small amount of etchant directly from the syringe. Remove any visible plug, if present, from the syringe opening.

Replace dispensing tip and again express etchant. If clog remains, discard dispensing tip and replace with a new one. Bend the dispensing tip to a desired angle. Place bend midway along tip. Do not bend dispensing tip at its hub. The small chamfered hole in the three-well tray was designed for bending syringe tips.

- Delivery system storage: Remove used dispensing tip and discard. Twist on storage cap. Storage of the delivery syringe with a used dispensing tip or without the storage cap will allow drying of the etchant and consequent clogging of the system. Replace storage cap with a new dispensing tip at next use.

- If desired, the etchant may be extruded onto a dispensing pad and applied with a brush or other appropriate instrument.

5. If a liquid etchant is desired, the etchant may be dispensed into a dappen dish and stirred to increase its fluidity.

- Disinfect: Discard used dispensing tip. Replace syringe storage cap.

Disinfect this product using an intermediate level disinfection process (liquid contact) as recommended by the Centers for Disease Control and endorsed by the American Dental Association. Guidelines for Infection Control in Dental Health-Care Settings – *MMWR*, December 19, 2003:52(RR-17), Centers for Disease Control and Prevention.

Direct light cure restorations in enamel and dentin

- Isolation: Rubber dam is the preferred method of isolation.

- Cavity preparation: Prepare cavity with minimal tooth reduction. Bevel cavosurface enamel margins.

- Etching: Apply etchant (phosphoric or maleic) to enamel and dentin. Wait 15 seconds. Rinse for 15 seconds. Dry for 5 seconds.

- Priming: Apply Adper Scotchbond multi-purpose primer to etched enamel and dentin. Dry gently for 5 seconds.

- Adhesive application: Apply Adper Scotchbond multi-purpose adhesive to the primed enamel and dentin.

- Adhesive curing: Light cure for 10 seconds.

7. Restorative placement, cure and finishing: Refer to manufacturer’s instructions for placement, cure and finishing of restorative material.

Direct light cure restoration of noncarious cervical lesions

- To restore noncarious cervical lesions, clean the lesion and adjacent enamel with a plain pumice/water slurry. Do not use prophy pastes containing oils. Remove excess surface moisture with a brief exposure to an air syringe or by blotting.

- Leave dentin moist. Follow the procedural steps 1 through 7 above limiting step #2 to beveling enamel adjacent to the erosion.

Direct light cure restorations in enamel only

1. In all-enamel preparations, use of the primer is not required. Application of the primer will neither adversely affect nor significantly enhance the bond of the adhesive to etched enamel surfaces. For restorations in enamel only, follow the procedural steps above omitting step #4.

Bonding amalgam or self cure composite to tooth structure and existing amalgam or composite

- Isolation: Rubber dam is the preferred method of isolation.

- Cavity preparation: Prepare a standard amalgam preparation. If previous amalgam or composite has not been totally removed, roughen the residual restorative with a sandblast technique or a bur.

- Matrix application: Lightly rub the internal surface of the matrix band with wax before placement.

- Etching: Apply etchant (phosphoric or maleic) to enamel, dentin, and the residual restorative. Wait 15 seconds. Rinse for 15 seconds. Remove excess water with an air syringe or by blotting. Appearance should be shiny. If surface has matte appearance it should be remoistened.

- Activator/Primer: Mix one drop each of Activator and Primer. Apply to etched enamel and dentin – wait 15 seconds. Dry gently for 5 seconds.

- Adhesive application: Mix one drop each of Adper Scotchbond multi-purpose adhesive and catalyst. Triturate amalgam. Apply the mixed adhesive to the primed surfaces. Thin the mixture using a brush tip, if desired. Avoid excessive air thinning.

- Condense and burnish the amalgam or place the self cure composite in the usual way.

Bonding porcelain veneers

- Silane Treatment: Porcelain bonding surfaces should have been etched with hydrofluoric acid by the dental laboratory. Apply RelyX ceramic primer to the bonding surface of the restoration. Dry for 5 seconds.

- Clean the prepared teeth in preparation for seating and bonding using a plain flour of pumice slurry. Rinse and dry thoroughly and isolate from moisture.

- Try in veneer with 3M™ ESPE™ RelyX™ Try-In Paste. After try in, isolate from adjacent teeth with clear matrix strip.

Etching: Apply etchant (phosphoric or maleic) to both enamel and dentin. Wait 15 seconds. Rinse for 15 seconds. Dry for 2 seconds. Leave moist. The etched enamel should display a “frosted white” appearance. An additional 15 second etch time may be appropriate for teeth that were not prepared using a diamond or bur. Residual organic matter can also require additional etch time.

- Priming: Apply Adper Scotchbond multi-purpose primer to etched enamel and dentin. Dry gently for 5 seconds.

- Adhesive application to the tooth: Apply Adper Scotchbond multi-purpose adhesive to the primed enamel and dentin. Do not light cure.

- Adhesive application to veneer: Apply Adper Scotchbond multi-purpose adhesive to the silane treated bonding surface of the veneer.

- Luting material application to veneer: Apply 3M™ ESPE™ RelyX™ Veneer Cement restorative to the bonding surface of the veneer.

- Seating and curing: Carefully seat the veneer. Clean excess luting materials from the veneer margins. Cure each area of the veneer for times recommended by the luting cement manufacturer. We recommend curing the gingival margin first, followed by the body and the incisal margin. Avoid direct contact with the light guide during initial curing.

Bonding inlays, onlays, crowns and Maryland (adhesive) bridges

- Prepare the bonding surface of the indirect restoration and the core buildup, if applicable. Porcelain bonding surfaces should have been etched with hydrofluoric acid by the dental laboratory. Metal and amalgam bonding surfaces should be roughened, preferably using a sandblast technique. Composite surfaces also benefit from roughening with a diamond. Light cure glass ionomer buildups should be pumiced.

- Silane Treatment: Apply RelyX ceramic primer to the bonding surfaces of porcelain or metal indirect restorations. Dry for 5 seconds.

- Clean the prepared teeth in preparation for seating and bonding using a plain flour of pumice slurry. Rinse and dry thoroughly and isolate from moisture.

- Etching: Apply etchant (phosphoric or maleic) to enamel, dentin, and the residual restorative. Wait 15 seconds. Rinse for 15 seconds. Remove excess water with an air syringe or by blotting. Appearance should be shiny. If surface has matte appearance it should be remoistened.

- Activator/Primer: Mix one drop each of Activator and Primer. Apply to etched enamel and dentin – wait 15 seconds. Dry gently for 5 seconds.

- Catalyst application to the tooth: Apply Adper Scotchbond multi-purpose plus catalyst to the primed enamel, dentin and core material.

- Catalyst application to the indirect restoration: Apply Adper Scotchbond multi-purpose plus catalyst to the treated bonding surface of the indirect restoration.

- Apply dual cure or self cure luting material to the bonding surface of the restoration.

- Seat the restoration. Remove excess paste from the margins. If using a dual cure luting material, light curing the margins is recommended to ensure maximum physical properties. Set time is dictated by the choice of luting material.

Bonding cast and prefabricated endodontic posts

- Prepare the endodontically treated tooth to receive the post. Trial fit, and adjust the post as needed. Bond to cast posts will be enhanced if the metal is surface treated as recommended under the section “Instructions for bonding inlays, onlays, crowns and Maryland (adhesive) bridges,” numbers 1 and 2.

- Etching: Apply etchant (phosphoric or maleic) to the prepared tooth. Wait 15 seconds. Rinse for 15 seconds. Dry for 2 seconds. Use a paper point to remove any excess water in the canal.

- Activation: Apply Adper Scotchbond multi-purpose plus activator to the canal using a paper point. Dry for 5 seconds.

- Priming: Apply Adper Scotchbond multi-purpose primer to the canal using a paper point. Dry for 5 seconds.

- Catalyst: Apply Adper Scotchbond multi-purpose plus catalyst using a paper point.

- Post Preparation: Apply a coating of Adper Scotchbond multi-purpose plus catalyst to the post. Mix and apply the luting agent to the post.

- Seat the post. If a dual cure cement was used, cure from the occlusal surface to allow immediate preparation of the post. Set time is dictated by the choice of luting material. Core material can now be bonded to post and tooth structure.

Bonding porcelain and metal orthodontic brackets to porcelain, metal and composite crowns

Porcelain or metal orthodontic brackets may be bonded to crowns using the Adper Scotchbond multi-purpose system. Bonding to silane treated porcelain can result in very high bonds, and debonding may remove areas of porcelain from a porcelain fused to metal crown. In this event the crown surface can be repaired using the porcelain repair procedure.

- Clean the crown in preparation for seating and bonding using a plain flour of pumice slurry. Rinse and dry thoroughly and isolate from moisture.

- Etching: Apply etchant (phosphoric or maleic) to the prepared surface of the crown. Wait 15 seconds. Rinse for 15 seconds. Dry for 2 seconds. Leave moist.

- Silane priming: Apply RelyX ceramic primer to the prepared surface and dry.

- Activation: Apply Adper Scotchbond multi-purpose plus activator to the etched surface of the crown. Dry the activator with a gentle air stream for 5 seconds.

- Priming: Apply Adper Scotchbond multi-purpose primer to etched surface of the crown. Dry gently for 5 seconds.

- Catalyst application to the crown: Apply Adper Scotchbond multi-purpose plus catalyst to the primed surface of the crown.

- Catalyst application to the orthodontic bracket: Apply Adper Scotchbond multi-purpose plus catalyst to the bonding surface of the bracket.

- Apply dual cure or self cure luting material to the bonding surface of the orthodontic bracket.

- Seat the orthodontic bracket. Remove excess paste from the margins. If using a dual cure luting material, light cure the margins to ensure maximum physical properties. Set time is dictated by the choice of luting material.

Porcelain repair using light cure composite

- Clean the surface to be repaired with a slurry of plain flour of pumice. Rinse and dry thoroughly and isolate from moisture.

Composite repair using light cure composite

- Clean the composite surface and tooth structure using a plain flour of pumice slurry. Rinse and dry thoroughly and isolate from moisture. Note: If bonding only to composite, pumice the composite, apply primer and adhesive as indicated below (etching is not necessary in this case).

- Roughen composite surfaces with a bur or diamond.

- Etching: Apply etchant (phosphoric or maleic) to both enamel, dentin and composite. Wait 15 seconds. Rinse for 15 seconds. Dry for 2 seconds.

- Priming: Apply Adper Scotchbond multi-purpose primer to etched enamel, dentin and composite. Dry gently for 5 seconds.

- Adhesive application: Apply Adper Scotchbond multi-purpose adhesive to the primed tooth structure and composite.

- Adhesive curing: Light cure for 10 seconds.

- Restorative placement, cure and finishing: Refer to manufacturer’s instructions for placement, cure and finishing of restorative material.

Composite repair using self cure composite

Follow steps 1-3. Apply Adper Scotchbond multi-purpose plus activator and dry for 5 seconds. Apply Adper Scotchbond multi-purpose primer and dry.

Apply a 1:1 mixed solution of adhesive and catalyst. Place self cure composite.

Bonding porcelain or metal to existing porcelain, metal or composite

- Prepare the new bonding surfaces. Porcelain bonding surfaces should have been etched with hydrofluoric acid by the dental laboratory. Metal and amalgam bonding surfaces should be roughened, preferably using a sandblast technique. Composite surfaces also benefit from roughening. Apply RelyX ceramic primer to the new bonding surfaces and dry. The surface will appear dull.

- Clean the existing oral bonding surfaces in preparation for seating and bonding using a plain flour of pumice slurry. Rinse and dry thoroughly and isolate from moisture.

- Etching: Apply etchant (phosphoric or maleic) to existing oral bonding surfaces. Wait 15 seconds. Rinse for 15 seconds. Leave moist.

- Apply RelyX ceramic primer to the oral bonding surfaces. Dry for 5 seconds.

- Activation: Apply Adper Scotchbond multi-purpose plus activator to existing oral bonding surfaces. Dry the activator with a gentle air stream for 5 seconds.

- Priming: Apply Adper Scotchbond multi-purpse primer to activated surfaces. Dry gently for 5 seconds. The surface will appear shiny.

- Catalyst: Apply Adper Scotchbond multi-purpose plus catalyst to the primed existing oral bonding surfaces.

- Catalyst application to the indirect restoration: Apply Adper Scotchbond multi-purpose plus catalyst to the bonding surface of the indirect restoration.

- Apply dual cure or self cure luting material to the bonding surface of the restoration.

- Seat the restoration. Remove excess paste from the margins. If using a dual cure luting material, light cure the margins. Set time is dictated by the choice of luting material.

Bonding light cure composite to set amalgam

- Prepare tooth as desired.

- Roughen the set amalgam surface using a sandblast technique. The remaining steps duplicate the use of Adper Scotchbond multi-purpose system for direct composite placement.

- Etching: Apply etchant (phosphoric or maleic) to both enamel, dentin and set amalgam. Wait 15 seconds. Rinse for 15 seconds. Dry for 2 seconds. Leave moist.

- Priming: Apply Adper Scotchbond multi-purpose primer to etched enamel, dentin and set amalgam. Dry gently for 5 seconds.

- Adhesive application: Apply Adper Scotchbond multi-purpose adhesive to the primed tooth structure and amalgam.

- Adhesive curing: Light cure for 10 seconds.

- Masking: Mask amalgam with an appropriate masking agent. Cure.

- Restorative placement, cure and finishing: Refer to manufacturer’s instructions for placement cure and finishing of restorative material.

Storage and use

This product is designed to be stored and used at room temperature. Shelf life at room temperature is 36 months. Ambient temperatures routinely higher than 27°C/80°F, or lower than 10°C/50°F may reduce shelf life. See outer package for expiration date.

Adper Scotchbond multi-purpose plus activator and RelyX ceramic primer are ethanol solutions and should be recapped immediately after dispensing to reduce evaporation. Do not expose materials to elevated temperature or intense light. Do not store materials in proximity to eugenol-containing products.

For cleanup, Adper Scotchbond multi-purpose etchant and primer can be removed with water, while the activator, uncured adhesive and catalyst can be removed with alcohol.

Disinfect this product using an intermediate level disinfection process (liquid contact) as recommended by the Centers for Disease Control and endorsed by the American Dental Association. Guidelines for Infection Control in Dental Health-Care Settings – *MMWR*, December 19, 2003:52(RR-17), Centers for Disease Control and Prevention.

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user’s application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE’s sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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- Porcelain preparation: Use a diamond to remove all weakened porcelain and roughen the surfaces to be bonded. Bevel the margin and remove 1 mm of surface glaze beyond the margin.

- Metal preparation: Roughen areas of exposed metal with diamond or sandblast technique.

- Surface treatment: Cleanse the metal and porcelain surfaces by applying etchant for 15 seconds, followed by rinsing and drying.

- Silane priming: Apply RelyX ceramic primer to the prepared surface and dry.

- Adhesive application: Apply Adper Scotchbond multi-purpose adhesive. Light cure for 10 seconds.

- Masking: Mask any metal surfaces with an appropriate masking agent. Cure.

- Restorative placement, cure and finishing: Refer to manufacturer’s instructions for placement, cure and finishing of restorative material.

Porcelain repair using self cure composite:

Follow steps 1-5. Apply Adper Scotchbond multi-purpose plus activator and dry for 5 seconds. Apply Adper Scotchbond multi-purpose primer and dry for 5 seconds.

Apply a 1:1 mixed solution of adhesive and catalyst. Place self cure composite.

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