## ENGLISH

### General Information

The **3M™ ESPE™ Sof-Lex™ Diamond Polishing System** is a 2-step, multi-use polishing system. The Spirals use the same mandrel as 3M™ ESPE™ Sof-Lex™ Finishing and Polishing Discs. The Spirals are made with either aluminum oxide or diamond particles impregnated in a thermoplastic elastomer.

### Indications

The **3M™ ESPE™ Sof-Lex™ Pre-Polishing Spiral** (beige) is indicated to smooth the surface of restorations. The **3M™ ESPE™ Sof-Lex Diamond Polishing Spiral** (pink) is indicated for final high gloss polishing.

Spirals can be used to polish surfaces of:

- Direct:
- Composite restorations
- Resin modified glass ionomers
- Bisacrylic temporary materials
- Indirect:
- Composite
- 3M™ ESPE™ Lava™ Ultimate CAD/CAM Restorative
- Precious and semi-precious metal

### Precautinary Instruction for Dental Personnel


Wear eye protection when using this product.

**DO NOT** use with high speed hand pieces.

Operating hand piece at speeds greater than 25,000 rpm may cause the Spiral and mandrel to separate or fragment resulting in injury.

Excessive, heavy pressure (over 100 grams force) may cause a Spiral to separate or fragment resulting in injury.

Contacting the hub (center of the wheel near the mandrel side) of the Spiral with teeth (especially occlusal) for incisal edge) may result in instruments separating from mandrel causing injury.
The useful life of the Sof-Lex Pre-Polishing and Diamond Polishing Spirals depend on many factors, including the method and duration of each use, and the handling between uses. Careful inspection before use is the best method to determine the end of life. The Spirals are safe for reuse according to the instructions below until bristles are notched, or missing. Inspect each Spiral carefully before use and discard if any bristles are notched (A) or missing (B). Refer to the instructions for reprocessing the Spirals.

Spirals are designed to be used **only** with 3M™ ESPE™ Sof-Lex™ mandrels. Some resistance should be felt when placing a Spiral on a mandrel. If fit appears loose, replace the instrument.

Heads of Sof-Lex mandrels will wear over time. Replace mandrel if Spirals persistently feel loose.

3M ESPE MSDS information can be obtained from [www.3MESPE.com](http://www.3MESPE.com) or contact your local subsidiary.

### Instructions for Use

- Gross reduction. Use 3M™ ESPE™ Sof-Lex™ Coarse Discs or trimming and finishing diamonds or carbide burs to remove excess (or air inhibited layer) restorative material.
- Final contour: Complete and refine final restoration contours and margins with 3M™ ESPE™ Sof-Lex™ Medium Discs, fine or very fine diamond or finishing burs. The restoration should have desired form and a smooth, clean surface.
- Pre-Polish: Smooth surface with Pre-Polishing Spiral (beige).
  - Insert a Sof-Lex mandrel into a slow speed hand piece that is designed for RA mandrels.
  - Attach the Pre-Polishing Spiral (beige) by pushing it onto the Sof-Lex mandrel.
  - A forward handpiece direction is recommended for operation of the Spirals. Turning the Spirals in reverse direction can alter the shape of the instruments significantly and reduce polishing effectiveness.
  - Optimum results are obtained when the Spirals are used on wet surfaces, with moderate pressure, and a low-speed handpiece rpm range of 15,000-20,000.
    - For best results, use medium (30 grams force bends the bristle tips) to moderate pressure (60 grams force causes a significant bend in the bristles) with hand piece speeds of 15,000-20,000 rpm (typically the maximum speed on a slow speed handpiece) and keep the instrument in motion over the restoration surface. If the hub of the wheel is being deformed (bent), too much pressure is being applied.
    - Intermittent to continuous water spray while polishing will keep the polishing surface wet, optimize results and prolong tool life.
- All surfaces of the Spiral (top, bottom, edge) can be utilized for maximum access to the restorative surface.
- When satisfied with results of the Pre-Polishing Spiral (beige), remove wheel from the mandrel by grasping the edge and pulling upwards.
- Rinse debris from restoration and tooth surfaces.
- Polish: High gloss polish with Diamond Polishing Spiral (pink).
  - For highest gloss, place the Diamond Polishing Spiral (pink) onto the mandrel and repeat steps 3.3 to 3.6.
  - After bringing surface to high gloss polish, rinse with water and wipe off restoration surface as necessary.

5. No other additional steps are necessary.

### Reuse Instructions:

Reprocessing Instructions
<b>General</b> This medical device is non-sterile when delivered. It may be used non-sterile for the initial use, but must be cleaned, sterilized and inspected before every re-use.
<b>Warnings/Precautions</b> Make sure that the cleaning agents you have chosen do not contain any of the following materials: <ul style="list-style-type: none"><li>Organic solvents</li> <li>Oxidation agents</li> <li>Halogens</li> <li>Oils</li> <li>Glutaraldehydes</li></ul> Never use metal brushes or steel wool to clean the device. The medical device must not be exposed to temperatures higher than 138°C (280°F).
<b>Limitations on reprocessing</b> The useful life of the Spiral depends on many factors, including the method and duration of each use, and the handling between uses. Careful inspection of each Spiral before use is the best method to determine the end of life. Inspect each spiral carefully before use. The Polishing Spirals are safe for reuse according to the instructions below until a bristle is notched or missing, or the Spiral appears loose on the mandrel.
<b>Point of Use Pre-Treatment</b> The device should be placed into a lidded container for transportation to the reprocessing location.
<b>Preparation for decontamination</b> As soon as practical after use, remove gross contamination from the device under running water or with a wet cloth. <ul style="list-style-type: none"><li>Rinse the device under running water (temperature &lt;35°C/95°F) for 1 minute.</li> <li>Use a clean, soft, lint-free cloth to manually remove all visible contamination.</li> <li>Rinse again for 1 minute under running water.</li></ul>

**Cleaning: Automated (Washer-Disinfector)**  
Cleaning may be performed using a washer-disinfector using a mild (pH 7.5-11) enzymatic detergent, such as Cidezyme® Xtra or Neodisher® MedClean Forte.

- The cleaners being used must be compatible with the device (see "Warnings").
- Follow the detergent manufacturer's instructions for concentrations, temperatures, exposure times, and rinsing instructions exactly.
- Place the device in the washer basket, making sure that there is no contact between the devices.
- Recommended Cycle:

Phase	Recirculation Time	Water Temperature
Pre-Wash	4 minutes	45°C
Wash	3 minutes	60°C
Rinse	1 minute	20-25°C

**Cleaning: Automatic (Ultrasonic Washer)**  
Cleaning may be performed ultrasonically using an enzymatic detergent (pH 7.5-11), such as Enzo®®, Cidezyme® Xtra or Neodisher® MedClean Forte.

- Follow the detergent manufacturer's instructions for concentrations, temperatures, exposure times, and rinsing instructions exactly. Pour the detergent into the ultrasonic washer.
- Place the Spiral in a wire container and rinse under running tap water for 10 seconds.
- Immerse the container with the Spiral in the detergent and sonicate for 5 minutes.
- Remove the container with the Spiral from the detergent and rinse in tap water for 1 minute.
- Remove the Spiral from the container and immerse in deionized water for 30 seconds.

### Cleaning: Manual

Cleaning may be performed manually using an enzymatic detergent (pH 7.5-11) enzymatic detergent, such as Enzo®®, Cidezyme® Xtra or Neodisher® MedClean Forte.

- Follow the detergent manufacturer's instructions for concentrations, temperatures, exposure times, and rinsing instructions exactly.
- Rinse the Spiral under running tap water for 10 seconds.
- Immerse the Spiral in the detergent solution and soak for 1 minute.
- Remove the Spiral from the detergent solution and brush all surfaces for a total of 15 seconds using a non-abrasive nylon toothbrush-style brush.
- Rinse the Spiral under running tap water for 10 seconds. Immerse the Spiral in deionized water for 30 seconds.

### Disinfection

A final rinse at 93°C for 3 minutes, or 90°C for 5 minutes (A<sub>1</sub> 3000) can be used to achieve a thermal disinfection.

### Drying

Dry with a clean, soft, lint free cloth as needed. When drying is performed as part of the washer-disinfector cycle, do not exceed 95°C.

### Maintenance, Inspection and Testing

- Examine the Spiral for notched or missing bristles before every use and dispose of damaged Spirals.
  - Any Spirals which is still contaminated must be cleaned a second time.
  - Store the Spirals under dry and dust-free conditions.
  - Never use instrument oils or greases.
- Packaging**  
Package the Spiral in disposable sterilization packaging (single or double packaging). The packaging must meet the following requirements:
- DIN EN ISO/ANSI AAMI ISO 11607
  - For USA: FDA clearance
  - Suitable for use in steam sterilization (temperature stability of at least 142°C (288°F), adequate steam permeability).

### Sterilization

Steam sterilization is the only sterilization method which may be used. Optimum results are obtained when the Spirals are used on wet surfaces, with moderate pressure, and a low-speed handpiece rpm range of 15,000-20,000.

- Do not exceed the maximum sterilization temperature of 134°C (273°F; tolerances pursuant to DIN EN ISO 17665).
- For best results, use medium (30 grams force bends the bristle tips) to moderate pressure (60 grams force causes a significant bend in the bristles) with hand piece speeds of 15,000-20,000 rpm (typically the maximum speed on a slow speed handpiece) and keep the instrument in motion over the restoration surface. If the hub of the wheel is being deformed (bent), too much pressure is being applied.
- Intermittent to continuous water spray while polishing will keep the polishing surface wet, optimize results and prolong tool life.

- All surfaces of the Spiral (top, bottom, edge) can be utilized for maximum access to the restorative surface.
- When satisfied with results of the Pre-Polishing Spiral (beige), remove wheel from the mandrel by grasping the edge and pulling upwards.
- Rinse debris from restoration and tooth surfaces.
- Polish: High gloss polish with Diamond Polishing Spiral (pink).
  - For highest gloss, place the Diamond Polishing Spiral (pink) onto the mandrel and repeat steps 3.3 to 3.6.
  - After bringing surface to high gloss polish, rinse with water and wipe off restoration surface as necessary.

The instructions provided above have been validated by 3M ESPE as being CAPABLE of preparing this medical device for re-use. It remains the responsibility of the processor to ensure that the reprocessing is actually performed using equipment, materials and personnel in the reprocessing facility to achieve the desired result. This normally requires validation and routine monitoring of the process.

### Disposal:

Properly dispose of used Spirals in accordance with your healthcare facility's policy.

### Sof-Lex Mandrel re-use:

Sof-Lex mandrels must be cleaned and high level disinfected or sterilized between patients. See Sof-Lex mandrel Instructions for Use for further detail.

### Customer Information

No person is authorized to provide any information which deviates from the intation provided in this instruction sheet.

**Caution:** U.S. Federal Law restricts this device to sale or use on the order of a dental professional.

### Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the products for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be the repair or replacement of the 3M ESPE product.

### Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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