BIO-C[®] SEALER

CHOOSE YOUR LANGUAGE

ENGLISH.....	4	HRVATSKI	55
ESPAÑOL.....	8	SLOVENŠČINA	59
PORTUGUÊS	12	SLOVENSKY	63
FRANÇAIS	16	POLSKI	67
ITALIANO	21	EESTI KEEL.....	71
ROMÂNĂ	25	LIETUVIŠKAI	75
NEDERLANDS.....	29	ČESKY	79
DEUTSCH	33	MAGYAR	83
SVENSKA	38	TÜRKÇE	88
NORSK.....	42	БЪЛГАРСКИ	92
DANSK	46	СРПСКИ.....	97
SUOMI	51	ΕΛΛΗΝΙΚΑ	101

ENGLISH

INTENDED PURPOSE

BIO-C® SEALER is a bioceramic root canal sealer, intended for exclusive use by dentists during endodontic obturation procedures. It is designed to fill and seal previously cleaned and shaped root canal systems, aiming to prevent reinfection and ensure long-term obturation. The device is supplied ready-to-use.

INDICATIONS

- Obturation of root canals of permanent teeth.
- Use during internal resorption treatment.

INTENDED POPULATION

Patients with permanent teeth with closed apices, diagnosed with pulpal or periapical disease.

INTENDED USERS

BIO-C® SEALER should be used by dentists trained to perform endodontic procedures.

INTENDED ENVIRONMENT

Dental offices.

CONTRAINDICATIONS

- Do not use the product after the expiration date.
- Do not use the product if the packaging is damaged.
- Do not use in patients with reported sensitivity to any of the formula components.
- Do not reuse the applicator tips.

COMPOSITION

Calcium silicate, calcium aluminate, calcium oxide, zirconium oxide, iron oxide, silicon dioxide, and dispersing agent.

TECHNICAL INFORMATION

- Biocompatible product.
- Ready-to-use presentation, no mixing required.
- Setting time: 30–240 minutes after insertion into the canal.

Important: Dry with paper points but avoid over-drying the canal. Moisture from dentinal tubules will

initiate the material's setting reaction.

- High alkalinity.
- High radiopacity.
- No setting shrinkage.
- Can be used with conventional techniques, single cone, and thermoplastic obturation.

WARNINGS

- Avoid contact with eyes or skin. In case of contact, rinse immediately with water.
 - Avoid contact with oral mucosa. In case of contact, rinse with water and avoid ingestion of the product.
- If any sensitivity persists, seek medical attention immediately.
- Safety and efficacy have not been established in vulnerable populations.
 - Keep out of reach of children.

Note: Any serious incident related to the device must be reported to the manufacturer and the competent health authority in the country where the user and/or patient is located.

PRECAUTIONS

- Wear protective glasses, mask, and gloves when handling the product. In case of contact with eyes or skin, wash with water.
- Use absolute isolation in all indications.
- The product is moisture-sensitive, so close the packaging properly with adequate pressure to prevent drying of the product.
- Use the product only after infection control, as the material will not set in an acidic environment.
- Avoid excessive extrusion of the product into the periapical region.
- Disinfect the applicator tip prior to use. After use, the tip should be discarded.
- Clean and disinfect the syringe for subsequent uses.
- Return the syringe properly closed, along with the silica sachet, inside the packaging, which must be properly sealed and stored according to the "storage" section.
- For subsequent use, ensure the silica sachet inside the envelope remains blue.

TECHNIQUE OF USE

Use of Applicator Tips

1. A disposable plastic cover on the syringe is recommended to prevent cross-contamination.

2. Confirm the material flow from the syringe before direct application into the canal.
3. Position the tip and syringe assembly properly in the canal and press the plunger gently to avoid excessive extrusion of the product.
4. Before closing the syringe, pull back the plunger to stop the flow.
5. Clean excess material from the syringe tip and close with adequate pressure after use to prevent moisture entry that causes drying.
6. Clean and disinfect the syringe for subsequent uses to avoid cross-contamination; if contaminated with saliva or blood, discard it.

A. OBTURATION OF ROOT CANALS OF PERMANENT TEETH

Traditional Technique – Lateral Compaction

1. Administer anesthesia, place absolute isolation, then perform biomechanical preparation of the canal.
2. Remove the sealing tip and attach the applicator tip to the syringe.
3. Dry the canal only with paper points without over-drying.
4. Apply BIO-C® SEALER while gradually retracting the syringe to achieve complete canal filling
Important: press the plunger gently to avoid excessive extrusion of the product.
5. Insert the main gutta-percha cone according to the last file used, followed by accessory cones, all coated with BIO-C® SEALER.
6. Take a radiograph to verify proper canal filling.
7. Cut the cone at the desired height with heated instruments, followed by vertical compaction.
8. Remove excess material from canal walls with water, perform coronal sealing, and restoration.

Important:

- BIO-C® SEALER does not undergo shrinkage setting.
- BIO-C® SEALER can be used with conventional, single cone, and thermoplastic obturation techniques.
- Use conventional gutta-percha removal techniques.

Single Cone Technique

1. Administer anesthesia, place absolute isolation, then perform biomechanical preparation of the canal.
2. Remove the sealing tip and attach the applicator tip to the syringe.
3. Dry the canal only with paper points without over-drying.

4. Apply BIO-C® SEALER while gradually retracting the syringe to achieve complete canal filling. Important: press the plunger gently to avoid excessive extrusion of the product.
5. Insert the gutta-percha cone.
6. Take a radiograph to verify proper canal filling.
7. Cut the cone at the desired height with heated instruments, followed by vertical compaction.
8. Remove excess material from canal walls with appropriate irrigating solutions, perform coronal sealing, and restoration.

B. USE DURING INTERNAL RESORPTION TREATMENT

1. Administer anesthesia and place absolute isolation.
2. Remove granulation tissue from the resorption area with sharp curettes.
3. Neutralize the environment with calcium hydroxide paste or BIO-C® TEMP.
4. Remove calcium hydroxide or BIO-C TEMP in the following session.
5. Dry the canal only with paper points without over-drying.
6. Insert BIO-C® SEALER with the applicator tip throughout the canal, prioritizing the resorption site.
7. Insert accessory gutta-percha cones according to the lateral compaction technique.
8. Cut cones at the top of the resorption area with heated instruments. Perform vertical compaction to improve cement flow at the resorption site.
9. Take a radiograph to verify correct filling of the resorption site and canal.
10. Obturate the remainder of the canal with gutta-percha.
11. Perform coronal sealing with glass ionomer or another preferred material and restore the tooth.
12. Take a radiograph and maintain follow-up for at least two years.

Obturation removal

1. Use conventional gutta-percha removal techniques.

STORAGE

- Keep the product in a dry, ventilated place between 15 and 30°C with relative humidity below 60%.
- Do not store in a refrigerator.
- Do not store the product near ammonia, ammonium nitrate, or chlorine-containing products. Avoid using disinfectant solutions containing any of these ingredients.

DISPOSAL METHOD

Dispose of the product in accordance with national laws and local authority recommendations for proper disposal. Improper disposal may present biological or physical risks; handle and dispose of as hazardous biological material if contaminated.

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