R-Endo / R-Endo assorted

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EN

1. Indication

R-Endo: Desobturation/removal and shaping of root canal during non-surgical endodontics retreatment which intervene after an endodontic treatment.

For use by a dental professional only.

2. Contraindications

There is no contraindication to the use of R-Endo for endodontically treating a tooth by an orthograde route.

3. Complications

In cases of complex canal anatomy, per-operative risks (instrumental breakage, ledge, stripping, zipping, false path, perforation, etc.) could occur and lead to a risk of infectious processes.

4. Warning and precaution

- The decision to use the endodontic instrument must be based on the clinical case expertise, particularly where the canal anatomy is considered to be too complex.
- Patients identified as having a risk of infectious endocarditis.
- Contains nickel and titanium: should not be used in persons with a known allergy to these metals.
- Respect good dental practices, in particular by using a dental dam and gloves.
- Use in continuous rotation at the recommended speed.
- In cases of complex anatomy, the maximum torque may be different from the manufacturer's recommended value.
- Use according to the procedure recommended (§7).
- Instruments supplied non-sterile: Follow the reprocessing instructions (§8) before use.
- Keep UDI information present on the labelling only until the last use.
- Year of manufacture: see labelling.
- Do not use the instrument if there is any doubt about compliance with the storage conditions.
- If there is any doubt about the identification of the product, do not use it.
- Check the instrument integrity before using between each canal. If the instrument is damaged or shows signs of wear, do not use it.
- Inform the manufacturer and the regulatory authority in the country of any serious incident relating to the instrument.

5. Clinical claims

R-Endo in normal conditions.

Clinical performance:

- Shaping ability (debris removal and cutting efficiency)
- Respect of canal anatomy
- Geometry adapted to irrigation
- Maintains the apical foramen in its original position

Safety:

- Limited risk of breakage
- Limited peroperative risks
- Limited risk of screwing effects

6. Characteristics

Type 2 continuous rotating instruments under EN ISO 3630-1: 2019 standard:

- Maximum recommended torque: see table 'Recommended torque'.

	Ø	%	L	Stop color	Ring color	RPM	Recommended torque (N.cm)
R-Endo R1	25	8	15	●	●	300-400	1.5 - 2
R-Endo R2	25	6	19	●	●	300-400	1 - 1.5
R-Endo R3	25	4	23	●	●	300-400	0.6 - 1


Material of operative part: Nitinol.

Use in connection with an endodontic contra-angle according to standard EN ISO 1797: 2017 (Type 1).

Number of uses: maximum of 10 canals recommended if the instrument is not visually impaired.

Medical device class according to Directive 93/42/EEC and MDR 2017/745 : IIa.

7. Protocol

R1 : Shaping with removal of stresses in the coronal third. Picking motion to progress.

R2 : Shaping with removal of stresses in the median third.

R3 : Shaping up to the working length with removal of stresses in the apical third.

Irrigate abundantly between each step and each passage of the instrument.

8. Reprocessing instructions

General recommendations	- For all metal devices, the use of anticorrosion disinfecting and cleaning agents is recommended. - For your own safety, please wear personal protective equipment (gloves, glasses and mask). - Do not use cleaning or disinfecting agents containing phenol, aldehyde and alkaline composition. - Always respect the instructions for use provided by the manufacturer of the products.
Limitations on reprocessing	- Due to the product design and the materials used, no definite limit to the maximum number of performable processing cycles can be specified. The service life of the medical devices is determined by their function and careful handling. Multiple-use disinfection and re-sterilisation cycles may lead to increased risk of yile separation. - The user must ensure that the processing method used, including resources, materials and personnel, is appropriate and meets the applicable requirements. - The state of the art and national laws require that validated procedures be followed.
Material needed	- Gloves, masks, gown as recommended by the manufacturer of cleaner and detergent - Running or deionized water - Disinfectant (neodisher® Septo Active) - Detergent (neodisher® MediZym) - Small soft brushes - Container - Ultrasonic tub or washer-disinfector - Class B sterilisation apparatus <i>Remark: All material used should be cleaned and replaced regularly. Identify material used for each step of the process (initial treatment, cleaning or rinsing).</i>

1	Initial treatment Place used products in a container or a wipe with tap water at 20–40 °C and 1.0 % neodisher® Septo Active for 5 to 15 min. Rinse the products with tap water at 20–40 °C for 1 min. If there is waiting time before the next step, make sure the device stays moist by placing it in a wet wipe. Do not exceed 1 hour waiting time. <i>Remarks:</i> - Do not use fixating agents or hot water (> 40 °C), as this causes fixation of residues and can impair successful cleaning. - Follow instructions and respect concentrations and immersion times specified by the manufacturer (an excessive concentration may cause corrosion or other defects on devices).
2	Preparation before cleaning If the devices have visible impurities, it is recommended to manually pre-clean by brushing under tap water at 20–40°C at least 1 min with a soft brush until all impurities have been removed. <i>Remark: Follow instructions and respect concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or others defects on devices).</i>
3	Visual inspection Inspect the used products and discard the damaged products (broken, untwisted or abnormally bent).
4	Cleaning Insert the products into an ultrasound apparatus beaker. Run ultrasound apparatus for 10–30 min with tap water and 0.5–2.0% neodisher® MediZym. <i>Remarks: Follow the instructions, observe the water quality, concentrations and cleaning time stated by the manufacturer of the cleaning solution.</i>
5	Rinsing Rinse the products with tap water at 20–40°C for 1 min. <i>Remark: the use of deionised water is recommended.</i>
6	Drying Dry the products with compressed air until products are visibly dry.
7	Cleaning/Rinsing/Drying Place the instruments in the tray on the slide-in trolley of the washer/disinfector. Perform cleaning cycle with 0.2–1.0% neodisher® MediZym Perform drying. <i>Remarks:</i> - Disinfection (thermal or chemical-thermal) is not required since the products are sterilised after cleaning. - Follow instructions and concentrations given by the manufacturer of the detergent solution. - Follow the instructions of the washer-disinfector and verify the success criteria have been met after each cycle as specified by the manufacturer. - The final rinse step should be with deionized water. For other steps follow the water quality defined by the manufacturer. - Only use approved washers-disinfectors that comply with EN ISO 15883, and are maintained and validated regularly.
7	Visual inspection Inspect the used products. Re-do steps 4-5-6 if the product is visibly not clean or discard any damaged products.

8	Packaging Place the instruments in a paper-plastic pouch for steam sterilisation compliant with ISO 11607 and EN 868. <i>Remarks:</i> - For sharp devices that are not contained within a box, silicon tubes should be placed around the devices to prevent packaging being pierced. - Seal the pouches according to the recommendations of the pouch manufacturer. If a thermo-sealer is used, the process must be validated and the thermo-sealer must be calibrated and qualified.
9	Sterilisation Sterilise the products using steam vapor: - Apparatus: class B - Minimum temperature: 132°C - Minimum time: 3 min - Absolute pressure: 2.2 bar - Minimum drying time: 20 min Control the physico-chemical indicators and cycle parameters. The settings for temperature at 134°C and duration of 18 min are mandatory for prion inactivation according to French regulations. <i>Remarks:</i> - When sterilising multiple instruments in one autoclave cycle, ensure that the steriliser's maximum load is not exceeded. - Place the pouches in the steam steriliser according to the recommendations specified by the steriliser manufacturer. - Use only pre-vacuum air removal steam sterilisers that meet the requirements of EN 13060 (class B, small steriliser) and EN 285 (full size steriliser), with saturated steam.
10	Storage Store the products in a dry, clean and dust-free environment at the temperature specified on the paper-plastic pouch by the steam sterilisation manufacturer. <i>Remark: Check the packaging and the medical devices before using them (packaging integrity, no humidity and use-by date). In case of damage, a complete rerun should be performed.</i>

9. Storage and transport conditions

There are no special storage and transport conditions for non-sterile devices.

10. Disposal

After use, instruments must be placed in a secure container, used to collect cutting or sticking instruments (like needles or disposable bistouries) as per good dentistry practices.

11. Symbols

	Nitinol raw material		Root canal preparation
	Quantity		Assortment
	"CLASSICS" handle		Medical device
	Sterilizable in a steam sterilizer (autoclave) at the temperature specified		

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