Your practice is *our* inspiration.®

Instructions For Use

Demi^{Plus}

LED Light Curing System
Powered by Demetron™ Technology


Table of Contents

| | |
|-------------------------|-----|
| ENGLISH | 5 |
| FRANÇAIS | 9 |
| ESPAÑOL..... | 13 |
| DEUTSCH | 17 |
| ITALIANO | 21 |
| DANSK | 25 |
| PORTUGUÊS (BRASIL)..... | 29 |
| PORTUGUÊS | 33 |
| NEDERLANDS | 37 |
| SUOMI..... | 41 |
| ΕΛΛΗΝΙΚΑ | 45 |
| NORSK | 50 |
| SVENSKA | 54 |
| ČESKÝ | 58 |
| EESTI | 62 |
| MAGYAR | 66 |
| LIETUVIU..... | 70 |
| LATVIEŠU..... | 74 |
| POLSKI | 78 |
| ROMÂNEȘTE..... | 82 |
| РУССКИЙ | 86 |
| БЪЛГАРСКИ | 91 |
| SLOVENSKÝ | 96 |
| SLOVENSKO..... | 100 |
| SRPSKI | 104 |
| TÜRKCÈ | 108 |
| TIẾNG VIỆT | 112 |
| 中文体..... | 116 |
| 中文版..... | 120 |
| 日本語 | 124 |
| 한국어 | 129 |
| ภาษาไทย..... | 134 |
| عربي | 138 |

Illustrations

Fig. 1

Fig. 2


Fig. 3

Fig. 4

Fig. 5

Fig. 6


Fig. 7

Fig. 8

Fig. 9

Demi^{Plus} L.E.D Curing Light

Operator Manual

MADE IN USA

DESCRIPTION

The Demi^{Plus} L.E.D. (Light Emitting Diode) curing light is a visible curing unit intended for polymerization of light-cured materials by dental professionals.

IMPORTANT: Batteries must be fully charged prior to use. Refer to the CHARGE THE BATTERY(IES) section. There are no user serviceable parts inside. Do not open or tamper with the batteries.

SPECIFICATIONS AND EQUIPMENT CLASS

| | |
|---|--|
| AC Supply connection: | 100-240V AC/0.8-0.4A 47-63 Hz |
| Power input: | 12VA, nominal |
| Equipment class: | Class II  |
| Safety: | IEC 60601-1 |
| EMC (Electro-Magnetic Compliance): | IEC 60601-1-2 |
| Protection from electric shock: | Type BF  |
| Protection from ingress of liquids: | IPX0 (ordinary equipment) |
|  Operation: | Continuous operation patient application, duty cycle 20 seconds ON/1 minute OFF. |

OPERATING ENVIRONMENT

Note: Take caution when using this product in the presence of a flammable anesthetic mixture or an oxidizer such as oxygen or nitrous oxide.

- Ambient temperature: 60° F to 104° F (16° C to 40° C)
- Relative humidity: 10% to 80%, non-condensing
- Atmospheric pressure: 0.5-atm to 1.0-atm (500 hPa to 1060 hPa)

TRANSPORT AND STORAGE ENVIRONMENT

- Ambient temperature: -4° F to 104° F (-20° C to 40° C)
- Relative humidity: 10% to 85%
- Atmospheric pressure: 0.5 atm to 1.0 atm (500 hPa to 1060 hPa)
- Battery performance: Handpiece will generate 500, five-second cures with a fully charged battery
- Battery over-current protection: Electronic and Resettable polyfuse
- Light Emitting Diodes: An array of four LEDs die-mounted on a substrate.
- Output peak wavelength range: 450 - 470 nanometers
- Baseline output: 1,100 mW/cm²
- Standard light guide(s): Part Number 902495, 13mm to 8mm Extended Light Guide

Unit handpiece dimensions: Length: 9.25 in (23.5 cm); Width: 1.2 in (2.9 cm)

Charger dimensions: Length: 6.0 in (16.0 cm); Width: 4.63 in (11.1 cm)

Unit weight: Handpiece and battery: 5.5 oz (155 g); Battery charger: 12.4 oz (352 g)

SYMBOLS


Attention, consult accompanying documents: 

DC power:  

Waste Electrical and Electronic Equipment (WEEE): 

SAFETY PRECAUTIONS

- The fiber-optic light guide is provided in a non-sterile condition and must be sterilized prior to patient contact. For sterilization instructions, refer to the MAINTENANCE: CLEANING/DISINFECTING/STERILIZATING section.
- This curing unit produces high-output curing energy! A significant increase in curing energy is possible compared with equipment previously used. It is important to observe the following precautions and procedures:
 - Do not place light directly on or towards unprotected gingiva or skin.
 - Adjust curing techniques in accordance with the increase in curing energy.
 - Do not look directly at the light emitted from this curing unit. Do not use this device without suitable eye protection for the operator, assistant and patient. Refer to the ACCESSORIES AND REPLACEMENT PARTS section for suitable eye protection.
 - Persons having a history of photosensitive reactions or who are using photosensitizing drugs should not be exposed to light from this unit.
- **NOTE:** Take caution when using this product in the presence of a flammable anesthetic mixture or an oxidizer such as oxygen or nitrous oxide.
- All Demetron models have been designed to meet worldwide electrical safety standards, including U.S., Canada, and Europe, to be safe and effective for all dental applications.
- The Demi^{Plus} requires special precautions regarding EMC and must be installed and put into service according to the Guidance and Manufacturer's Declaration for EMC included in this manual.
- Portable and mobile RF communications equipment can affect the Demi^{Plus}. Reference the Guidance and Manufacturer's Declaration for EMC included in this manual.
- We are certain that you will receive many years of exceptional service from the Demi^{Plus} product. However, it is inevitable that at some time in the future you will need to dispose of it. When that time comes, we encourage you to return it to Kerr in its original packaging for environmentally sound recycling. Kerr will be glad to bear the cost of the return shipping. Please contact your sales rep for details. The Waste Electrical and Electronic Equipment (WEEE) Directive forbids the disposal of waste electrical and electronic equipment as unsorted municipal waste and requires that they be collected and recycled or disposed of separately. Hazardous substances are present in electrical and electronic equipment and present potential risks to human health and the environment if disposed of in municipal landfills which are not designed to prevent migration of substances into the soil and groundwater. By properly disposing of this waste, you will be contributing to the advancement of reuse, recycling, and other forms of recovery, as well as the prevention of pollution at a very local level.

NOTE: The dimensions in the figure represent the minimum patient environment (1.5 meters). The Demi^{Plus} charger base should not be stored within this patient environment.

ASSEMBLY/MOUNTING

The box contains:

| Qty. | Description |
|------|--|
| 1 | Demi ^{Plus} curing light handpiece |
| 1 | Protective light shield |
| 1 | Plug-in power supply |
| 1 | Battery charger |
| 1 | Battery pack |
| 1 | Operator manual |
| 1 | Demi/Demi ^{Plus} disposable barrier |

OPERATION

BATTERY CHARGER/BATTERIES

The battery pack contains Lithium ion cells. A new battery pack should be charged for 16 hours prior to its first use.

ATTACH THE ADAPTER TO THE POWER SUPPLY

The power supply is a universal-input unit that accepts 100V - 240V AC nominal. **Note:** Only Kerr-supplied power supplies (part number 921656) are to be used to power the Demi^{Plus} battery charger. There are two types of power supplies available as shown in Figure 1 and 3. The power supplies differ slightly in appearance, but operate the same both electrically and mechanically.

- I. Select the proper adapter (as shown in Figure 2 and 3). Their applications are as follows:
 - 120V US — Type A
 - Euro — Type C
 - UK — Type G
 - Australian — Type I
- II. Each adapter also has a slot into which the retaining latch on the power supply fits.
- III. With the slot oriented toward the power supply, place the adapter onto the power supply so that the tabs on the adapter are aligned with their respective slots on the power supply.
- IV. Pressing down on both ends of the adapter (Figure 4), slide the adapter into the power supply until a click is heard and the adapter is locked in place. **Note:** Pressing down on the adapter during installation assures that all the tabs are held by the retainers on the power supply. If the adapter feels loose, remove the adapter and re-install.
- V. Figure 5 shows the adapter locked in place on the power supply.
- VI. To remove the adapter from the power supply, press on the latch retainer button (7.), slide the adapter away from the power supply, and remove the adapter.
- VII. Retain unused adapters for future use.

INSTALL THE BATTERY ONTO THE WAND

- I. The battery will only fit onto the wand in one orientation. If you have difficulty fitting the battery to the wand, do not force it.
- II. Orient the battery (Figure 6) so that the broad flat area (1.) is toward the side of the wand that contains the switches.
- III. Slide the battery and wand together until a click is heard and the battery is locked in place.
- IV. To remove the battery (Figure 6), press the textured indents (2.) on both sides of the battery and slide the battery free.

CHARGE THE BATTERY(IES)

Plug the power supply into the charger base. The socket is located at the back of the charger base (Figure 7).

- I. Plug the power supply into an electrical outlet.
- II. Insert the battery and/or curing unit with battery into the opening(s) of the battery charger.
 - a. The charger can charge any combination of two batteries or wands (in differing states of charge) at the same time.
 - b. The lighted rings (Figure 7, 2.) that surround the openings in the charger base indicate the state of charge of each battery as follows:
 - Green light—battery is charged.
 - Yellow light—battery is charging.
 - Blinking Yellow Light (with battery or wand in the charger) — indicates a fault condition. Please call Customer Service.

Note: If the battery is below room temperature, allow the battery to come up to room temperature and try again; otherwise, replace the battery.

Battery pack life: 300 full charge/recharge cycles. Output: 3.7-Volts nominal @ 2.1A-H capacity.

BATTERY DISPOSAL: The battery pack contains a Lithium ion (Li-ion) battery. Recycle or dispose of batteries according to national, state and local regulations.

RECOMMENDED CURING TIME

The Demi^{Plus} is a high-output curing unit. Kerr testing indicates the following recommended curing times:

| Material | Recommended Curing Time |
|--|-------------------------|
| Universal composite shades A3 and lighter (2 mm depth) | 5 seconds |
| Universal composite shades A3.5 and darker | 10 seconds |
| Bonding agents, universal resin luting cements and Self-etch, self-adhesive cements. | 5 seconds |

Note: It is highly recommended that curing capability be tested using a Hardness Disk to ensure a complete cure.

HANDPIECE: The Demi^{Plus} is intended for all-day use without changing the battery. A fully charged battery will yield 500, five second cures.

The Demi^{Plus} is a wand-style curing light used for the polymerization of light-cured dental material and is powered by a rechargeable Lithium ion battery pack.

Note: When new, the battery pack must be charged prior to using the Demi^{Plus} Curing Light. The unit employs proprietary Periodic Level Shifting (PLS), patent pending. PLS is a constant periodic cycling of the output from a predetermined base level to a higher level for a short duration of time. The PLS technology provides faster, deeper cures without overheating.

DEMI^{Plus} PARTS AND CONTROLS

The following components are used by the operator during the curing process (Figure 8):

- I. **Curing Mode LEDs (1.)** Three LEDs represent the three curing modes; from left to right these are: 5 second, 10 second, and 20 second cycles.
- II. **Mode Selection Switch (2.)** Pressing this switch sequentially toggles the unit through the three curing modes. If the unit is in Sleep mode, pressing this switch awakens the unit in the curing mode last used.
- III. **Trigger Switch (3.)** Pressing (and releasing) this switch initiates the chosen curing cycle. If the unit is in Sleep mode, pressing this switch awakens the unit in the curing mode last used.
- IV. **Battery Power LED (4.)** This LED indicates the state of charge of the battery in the following way:
 - a) If the light is **unlit**, the battery is either **CHARGED** or the unit is in Sleep mode.
 - b) If the light is **yellow**, the battery has a **LOW** charge.
Note: The curing light will operate normally in this mode. Recharge the battery at the next available opportunity.
 - c) If the light is **red**, the battery is **DISCHARGED** (dead-battery) and will need to be recharged.

So that no curing cycles are interrupted, when the battery voltage drops below the dead-battery threshold during a curing cycle, the curing cycle will complete and the battery status LED will not change to red until the curing cycle is completed.

- V. **Lithium-Ion Battery Pack (5.)** Use to power the curing wand. Can be removed by pressing the indents and pulling the battery free.
- VI. **Light Guide Assembly (6.)** The output of the Curing Light LED is emitted from the end of this piece.*
- VII. **Light Shield (7.)** This shield protects the eyes of the operator from the light emitted from the Light Guide when in use.

*Other light guides that may yield different outputs are available for the Demi^{Plus}.

DEMI^{Plus} CURING LIGHT OPERATING STATES

There are three distinct operational states when operating the wand:

- I. **Idle:** The unit is Idle when it is not Curing or in Sleep mode and the battery is CHARGED or LOW. The user may cycle between 5 second, 10 second, and 20 second curing modes by pressing the mode switch. Curing cycles are also initiated when the unit is in the Idle state.
- II. **Curing:** A Curing cycle is initiated by pressing, and releasing, the trigger switch while the unit is Idle. A 5 second, 10 second, or 20 second cycle will be started, depending on the status indicated by the Curing Mode LEDs. Once a Curing cycle is started, the Curing Light LED will be turned on and a tone sounds. Depending upon which mode was selected, a tone will sound every five seconds until the cycle is complete, at which time a final, longer duration, tone sounds.
Note: Pressing the Trigger Switch for a second time during a Curing cycle immediately interrupts the Curing cycle and causes the final tone to sound, regardless of how much time has elapsed since beginning the cycle (pressing the Mode Selection Switch during a Curing cycle has no effect).
- III. **Sleep:** The Demi^{Plus} goes to Sleep mode to conserve battery life after approximately 8 minutes of no activity. It can be awakened by pressing either the mode or trigger switch once, at which time the unit will be in the Idle state. The unit will reflect the last cure mode chosen. The mode or trigger switches can then be pressed to switch between modes or to initiate a cure. While in Sleep mode, all LEDs are turned off and the curing light goes into a low power state. Battery life (if left uncharged and unused after receiving a full charge) is approximately 6 months whether the battery is in the handpiece or not.

OPERATE THE DEMI^{Plus} CURING WAND

- I. Select the desired curing mode by pressing the Mode Selection Switch (Figure 9, 1.)
- II. Place the Light Guide (2.) over the tooth, as close to the tooth as possible without touching it. Press and release the Trigger Switch (3.). The unit will sound a tone at 5 second intervals as previously described.

MAINTENANCE:

CLEANING / DISINFECTING / STERILIZING

- I. Fiber-Optic Light Guide:
 - a) Immerse/thoroughly drench light guide in a surface decontaminate cleaner (soap and water cleaning solution).
 - b) After 10 minutes, scrub with an instrument cleaning brush until visibly free of debris.
 - c) Rinse all parts with warm water.
- II. Fiber-Optic Light Guide - Sterilizing:
 - a) Sterilize in a steam autoclave using distilled water at a temperature of 270°F (132°C) for 20 minutes.
 - b) Do not add chemical solutions to the autoclave.
 - c) Do not use a chemical autoclave.
- III. Optics Maintenance Kit (PN 21042) – After numerous autoclave cycles, a brown or gray discoloration can accumulate on the polished ends of the light guide. This is a residue that forms during the autoclave drying cycle. It is recommended that the Optics Maintenance Kit be used every 50 autoclave cycles. Following the instructions included in the kit helps promote optimum curing performance.

To remove composite that adheres to the tip of the light guide:

- Soak gauze in alcohol.
- Wrap the gauze around a razor blade.
- Rub the tip of the light guide with the gauze-protected razor blade until the composite is removed.

Note: Never use a blunt instrument to remove composite because it will scratch the glass.

- IV. Handpiece and Battery Charger – Plastic parts can be cleaned following these guidelines:
 - a) All plastic parts may be wiped with CaviWipes™ following the direction on the package. A quaternary ammonium compound product like CaviCide™/CaviWipes™ or similar is recommended (containing 20% alcohol or less).
 - b) Do not use: Denatured alcohol, Lysol®, phenol, ammonia complex or iodine complex solutions.
 - c) Wipe, do not spray, solution onto plastic parts. Prevent liquids from entering openings on unit.

CAUTION: Do not immerse the plastic parts in solutions. Use of solutions other than those recommended may damage plastic parts and will void product warranty. This unit is not autoclaveable.

INFECTION CONTROL MEASURES

CAUTION: To prevent cross-contamination, a disposable plastic sleeve must be used over the Demi^{Plus} with each use.

WARRANTY

Kerr hereby warrants that for a period of three years (one year for batteries) from the date of purchase, this instrument shall be free from defects in material and workmanship and will perform satisfactorily under normal use and service. To register your warranty, complete the online warranty registration within 30 days of purchase at www.kerrdental.com/warranty. If multiple units have been purchased, fill out one registration form per unit. Please have your dealer invoice available. You can also register your warranty by calling 1.800.KERR.123 (1.800.537.7123).

THE WARRANTY STATED HEREIN IS THE SOLE WARRANTY APPLICABLE TO KERR PRODUCTS. KERR EXPRESSLY DISCLAIMS ANY AND ALL OTHER WARRANTIES EXPRESSED OR IMPLIED, INCLUDING WARRANTIES OR MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. KERR'S LIABILITY WITH RESPECT TO ITS PRODUCTS IS EXPRESSLY LIMITED TO THE REMEDIES SET FORTH ABOVE. THE REMEDIES ARE THE BUYER'S EXCLUSIVE REMEDIES. KERR SHALL UNDER NO CIRCUMSTANCES BE LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES. THE WARRANTY DOES NOT APPLY TO ANY DAMAGE OR DEFECT CAUSED BY THE USER'S FAILURE TO FOLLOW THE FOLLOWING INSTRUCTIONS FOR USE OR THE OTHER USER ERROR IN OPERATING OR MAINTAINING THE INSTRUMENT. ALTHOUGH KERR OFFERS SERVICES TO REPAIR INSTRUMENTS AND PROVIDES REPLACEMENT PARTS AFTER THE EXPIRATION OF THE WARRANTY, KERR ONLY GUARANTEES THE AVAILABILITY OF SUCH SERVICES AND REPLACEMENT PARTS FOR TWO YEARS AFTER THE INSTRUMENT MODEL HAS BEEN DISCONTINUED.

ACCESSORIES AND REPLACEMENT PARTS

| Part No. | Fiber Optic Light Guides |
|----------|--|
| 21020 | 8 mm Curved Turbo Light Guide |
| 921551 | 11 mm Curved Turbo Light Guide |
| 20812 | 13 mm Curved Light Guide |
| Part No. | Eye Protection |
| 20816 | Protective Light Shield |
| Part No. | Accessories |
| 21042 | Optics Maintenance Kit |
| 20399 | Disposable Hardness Disk, (pack of 6) |
| 910726 | Handheld LED radiometer |
| Part No. | Replacement Parts |
| 921918-1 | Replacement Kit, Battery, Demi ^{Plus} |
| 921919-1 | Replacement Kit, Charger, Demi ^{Plus} |
| 921920-1 | Replacement Kit, Wand, Demi ^{Plus} |
| 921656 | Plug-in Universal Power Supply |

The Demi^{Plus} is intended for use in the electromagnetic environment specified below. The customer or user of the Demi^{Plus} shall insure it is used in such an environment.

Guidance and manufacturer's declaration – electromagnetic emissions

| EMISSIONS TEST | COMPLIANCE LEVEL | ELECTROMAGNETIC ENVIRONMENT -GUIDANCE |
|---|------------------|---|
| RD emissions CISPR 11 | Group I | The Demi ^{Plus} uses RF energy only for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment. |
| RD emissions CISPR 11 | Class B | The Demi ^{Plus} is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power distribution network that supplies buildings or for domestic purposes. |
| Harmonic emissions | Class A | |
| Voltage fluctuations/flicker emissions IEC 60601-3-3 | Complies | |

Guidance and manufacturer's declaration – electromagnetic immunity

| IMMUNITY TEST | COMPLIANCE LEVEL | ELECTROMAGNETIC ENVIRONMENT -GUIDANCE |
|---|--|--|
| Electrostatic discharge (ESD) IEC 61000-4-2 | +/-4kV contact +/-8kV air | Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Electrical fast transient/burst IEC 61000-4-4 | +/-2kV for power supply lines +/-1kV for input output lines | Mains power quality should be that of a typical commercial or hospital environment. |
| Surge IEC 61000-4-5 | +/-1kV for differential mode +/-1kV for common mode | Mains power quality should be that of a typical commercial or hospital environment. |
| Voltage dips, short interruptions and Voltage variations on power supply input lines IEC 61000-4-11 | 30% <i>UT</i> for 0.5 cycle 60% <i>UT</i> for 5 cycles <5% <i>UT</i> for 5 sec | Mains power quality should be that of a typical commercial or hospital environment. Note: <i>UT</i> is the AC mains voltage prior to the application of the test level. |
| Power frequency (50/60Hz) magnetic field IEC 61000-4-8 | 3 A/m | Power frequency magnetic field should be at levels characteristic of a typical commercial or hospital environment. |
| Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3 | 10Vrms 150 KHz to 80 MHz 10V/m 80MHz to 2.5GHz | Portable and mobile RF communications equipment should be used no closer to any part of the Demi ^{Plus} including the power supply and cable than the recommended separation distances calculated from the equation applicable to the frequency of the transmitter. Field strength from fixed RF transmitters, as determined by an electromagnetic site survey, ¹ should be less than the compliance level in each frequency range ² . Interference may occur in the vicinity of equipment marked with the following symbol:  |

Note 1: At 80Mhz and 800Mhz the higher range applies

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

¹ Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radio, amateur radio, AM and FM broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. In the measured field strength in the location in which the Demi^{Plus} is used exceeds the applicable RF compliance levels above, the Demi^{Plus} should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Demi^{Plus}.

² Over the frequency range of 150Khz to 80Mhz, field strengths should be less than 10V/m.

The Demi^{Plus} is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the Demi^{Plus} can help prevent electromagnetic distance between portable and mobile RF communications equipment (transmitters) and the Demi^{Plus} as recommended below according to the maximum power output of the communications equipment.

Recommended separation distances between portable and mobile RF communications equipment and the Demi^{Plus}

| Radiated Maximum output power of transmitter W | Separation distance according to frequency m | | |
|---|---|--|---|
| | 150 KHz to 80Mhz $d = 0.4\sqrt{P}$ | 80 MHz to 800 MHz $d = 0.4\sqrt{P}$ | 800 MHz to 2.5 GHz $d = 0.7\sqrt{P}$ |
| 0.01 | 0.1 | 0.1 | 0.1 |
| 0.1 | 0.2 | 0.2 | 0.3 |
| 1 | 0.4 | 0.4 | 0.7 |
| 10 | 1.3 | 1.5 | 2.2 |
| 100 | 4.0 | 4.0 | 7.0 |

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80Mhz and 800Mhz the separation distance for the higher frequency range applies

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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Kerr Italia S.r.l.
Made in U.S.A.