TM

# PREMISE<sup>TM</sup>

Tri-Modal Composite Restorative

Instructions For Use

**PREMISE****Tri-Modal Composite Restorative**

PREMISE is a light-cured, resin-based composite dental restorative designed for direct placement. By utilizing three types of fillers (tri-modal), PREMISE is designed to offer high polishability, high mechanical strength, and decreased polymerization shrinkage. It is indicated for all restorations.

There is no need to refrigerate PREMISE, and it is recommended that the composite be stored at ambient temperature. The expiration date is marked on syringes and on Unidose® refill boxes. Do not use after expiration date.

It is also recommended that a reliable, high-power, visible light source (such as the Optilux 501, L.E.Demetron I, L.E.Demetron II, or Demi) be used to cure this product. Please note that PREMISE uses the photo-initiator camphorquinone which activates in the approximate 465 nm range.

Every effort should be made to protect the composite from visible light prior to use.

**PLEASE NOTE THE FOLLOWING:**

You have purchased this composite in either a bulk syringe or Unidose delivery. Either delivery has certain handling and storage requirements that if followed will provide excellent performance for the user over the life of the material. Please follow the appropriate recommendations indicated below for the delivery choice that you have made.

Syringe Delivery: To ensure optimal performance of syringes during use: After dispensing the composite, relieve pressure by turning the syringe handle back two turns. This will prevent material build-up in the cap and provide reliable composite use after storage.

Unidose Delivery: Unidose tips are designed for SINGLE PATIENT USE ONLY. Do not re-cap and/or re-use the Unidose tip once material has been dispensed for that patient.

**PRIOR TO PLACEMENT -- RECOMMENDATIONS ON PROPER BONDING**

- Isolate tooth from adjacent teeth. Use clear, thin strips. Rubber dam is highly recommended.
- Moist dentin surfaces (rather than dry, desiccated dentin) are shown to actually enhance the bond strengths of today's modern bonding agents (such as OptiBond Solo Plus, OptiBond FL, or OptiBond All-In-One). NOTE: "Moist" does not infer blood or saliva contamination.
- Placement sequences such as those outlined in this instruction booklet should be followed closely.
- Take a moment to identify the surface to which you are bonding.
- Etching dentin is recommended.
- Please take care to ensure that your air line is free of oil and other contaminants.

**OPTIBOND SOLO PLUS and OPTIBOND ALL-IN-ONE**

If you have purchased PREMISE in kit form, you will notice that OptiBond Solo Plus is included. This is the recommended adhesive for use with all Kerr composites. OptiBond Solo Plus is 15% filled with the same 0.4-micron barium glass found in the Kerr Point 4 composite. This unique technology provides a high level of protection against microleakage, while sustaining high-bond strengths to a variety of surfaces. The filler not only reinforces the hybrid zone, but penetrates the dentin tubules as well, creating a true "structural bond." With OptiBond Solo Plus, adhesion promoters are carried in an ethanol solvent, thus greatly reducing the tedious need for multiple coats and constant re-application commonly found with acetone adhesives. OptiBond Solo Plus is indicated for both direct and indirect (in conjunction with a luting agent) bonding applications.

OptiBond Solo Plus is available in two delivery systems, Bottle and Unidose. Kerr's Unidose delivery adhesive system was the first of its kind, designed as a progressive step toward ensuring fresh material for each bonding procedure as well as the ultimate in infection control. It is protected under U.S. Patent 5,860,806.

OptiBond Solo Plus is a "total-etch" adhesive. OptiBond All-In-One is a single-component, "self-etch" adhesive and may also be used with PREMISE.

#### **DIRECT BONDING PREMISE TO ENAMEL/DENTIN WITH OPTIBOND SOLO PLUS -- TOTAL-ETCH TECHNIQUE**

1. Isolate and dry.
2. Place Kerr Gel Etchant\* on enamel and dentin for 15 seconds.
3. Rinse thoroughly, ensuring that all etch is removed.
4. Dry lightly, do not desiccate.
5. If using Solo Plus in bottle delivery, shake bottle briefly before dispensing. If using Solo Plus in Unidose delivery, twist capsule at middle, and discard stem side.
6. Apply OptiBond Solo Plus to enamel/dentin surface with light brushing motion for 15 seconds.
7. LIGHTLY air thin for maximum 3 seconds.
8. Light cure for 20 seconds.\*\*
9. Proceed with placement of PREMISE composite.

Note: Kerr Disposable Mixing Wells, Gel Etchant tips, and Applicator Tips are designed for single use only.

#### **DIRECT BONDING PREMISE TO ENAMEL/DENTIN WITH ALL-IN-ONE -- SELF-ETCH TECHNIQUE**

1. Prepare cavity. Wash thoroughly with water spray and air dry (Do not desiccate).  
Note: For optimal adhesion, always prep (bevel) all enamel and/or margins prior to application of OptiBond All-In-One adhesive.
2. Shake OptiBond All-In-One adhesive bottle briefly.  
Note: Upon first use, shake OptiBond All-In-One adhesive bottle vigorously for 10 seconds.
  - Bottle Delivery: Dispense 2-3 drops of OptiBond All-In-One adhesive into a clean well. Replace cap immediately after dispensing.
  - Unidose™ Delivery (INTENDED FOR SINGLE PATIENT USE ONLY): Open the Unidose device and insert the applicator brush into the container to saturate the tip.
3. Using the disposable applicator brush, apply a generous amount of OptiBond All-In-One adhesive to the enamel/dentin surface. Scrub the surface with a brushing motion for 20 seconds.
4. Apply a second application of OptiBond All-In-One adhesive with a brushing motion for 20 seconds.
5. Dry the adhesive with gentle air first and then medium air for at least 5 seconds with oil-free air.
6. Light cure for 10 seconds.\*
7. Apply composite according to manufacturer's instructions for use.

#### **PLACEMENT OF PREMISE COMPOSITE**

PREMISE is available in a wide variety of shades. For details on all available shades and their various uses, please refer to "Explanation of PREMISE Shading System" below for more information.

Due to the decreased polymerization shrinkage of PREMISE, increments of 2.5 mm may be used. Please remember to cure each layer completely with a high quality curing unit such as the Optilux 501, L.E.Demetron I, L.E.Demetron II, or Demi.

1. Select the desired shade(s).
  - Syringe Delivery: After extruding material out on pad, make sure to back off the syringe at least two full turns to relieve pressure on the material.
  - Unidose Delivery: Insert tip and properly align within dispenser. Place Unidose tip at the deepest portion of the preparation.
- Caution: Extrude slowly with even pressure. Increments should be no more than 2.5 mm at a time.
2. After placing an increment, stroke the composite to ensure marginal adaptation.
3. Light cure each increment for 40 seconds.\*\*
4. When multiple surfaces are available, cure each surface for this recommended time.

## **EXPLANATION OF PREMISE SHADING SYSTEM**

The body and opaque shades of PREMISE are based on the VITA shading system. There are also 4 unique translucent and 2 extra light shades to assist you in obtaining your esthetic and restorative goals. Further explanation of these specialty shades is provided below.

### *Translucent Shades*

The translucent shades are designed to help create highly esthetic restorations with relative ease. Designed to be placed as a final layer in very thin depths, these specialty shades give the dentist everything he/she needs to virtually replace tooth structure with lifelike results. Our goal is to give you the tools to build true depth and beauty to every restoration.

Each translucent shade has a specific contrast ratio and slight color. Please refer to the table below for the approximate values and description:

<b>Shade</b>	<b>Description</b>	<b>Contrast Ratio</b>	<b>Translucency %</b>
Amber	Yellow/Red	45-55	45-55
Grey	Grey	45-55	45-55
Clear	Untinted	45-55	45-55
Super Clear	Untinted	20-30	75-80

### *Extra Light Shades (XL1, XL2)*

The extra light shades can be useful to match bleached as well as deciduous teeth. Both of them are lighter than the lightest VITA shade. In order to become familiar with these shades, some experimentation outside of the mouth is recommended.

The Starter Kit and the Masters Kit contain a shade guide to aid you in the shade taking process. The tabs in the shade guide are designed to be an accurate indicator of what each shade looks like after it is cured.

### *Packable Shades*

There are three packable shades available for your use when performing posterior composites. These shades, coded to the VITA shading system, have a higher viscosity than the other PREMISE shades. While they incorporate the same tri-modal technology as the other shades, the overall filler loading of the packable shades is higher. As a result, the shrinkage values of the packable material are slightly reduced relative to the other shades. Therefore, the recommended increment of the packable shades is 3mm, which should save the dentist time during the restorative process.

The following is a listing of the shades available with the PREMISE Composite System:

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B1 B2 B3 B4

C1 C2 C3 C4

D2 D3 D4

Opaque Shades: A2 A3 A3.5 A4

B1 B2 C2

D2

Translucent Shades: Amber

Grey

Clear

Super Clear

Extra Light Shades: XL1

XL2

(XL2 is whiter than XL1)

Packable Shades: A2 A3 A3.5 A4

Unidose Composite Gun Sterilization guideline: Steam Autoclave: 30 minutes @ 270°F.

## **OPTIGUARD**

OptiGuard is a composite sealant that can be used as a final overlay on top of all composite restorations. Modern high-speed finishing materials may create tiny cracks that can lead to premature composite wear, as well as continued propagation below the surface. Applying OptiGuard on the composite surface may increase the restoration's chances for extended wear and durability.

### ***DIRECTIONS FOR USE – OPTIGUARD***

1. Rinse final restoration and dry
2. Using Kerr Gel Etchant,\* etch the restorative surface and approximately 2mm of enamel beyond the margin for 15 seconds.
3. Rinse tooth thoroughly, and dry with oil-free air.
4. After dispensing a small amount of OptiGuard in mixing well, apply with Kerr Applicator to etched surface in a thin layer. Blow lightly to even the application on the surface.
5. Light cure the surface for 20 seconds.\*\*
6. Check occlusion when necessary.

**\* CAUTION:** Contains phosphoric acid. Avoid contact with skin, eyes, and soft tissue. In case of contact with skin or eyes, flush immediately with water. Get medical attention for eyes. Do not take internally.

**\*\* Recommended Cure Times:** Demi, 5 seconds, L.E.Demetron II, 5 seconds; L.E.Demetron I, 10 seconds; or Optilux 501 in Boost mode, 10 seconds, Ramp Mode, 20 seconds, or Regular Mode, 20 seconds. For all other lights, see manufacturer's recommendation.

**CAUTION:** Uncured methacrylate resin may cause contact dermatitis and damage the pulp. Avoid contact with skin, eyes and soft tissue. Wash thoroughly with water after contact.

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## **LIMITED WARRANTY LIMITATION OF KERR'S LIABILITY**

Kerr's technical advice, whether verbal or in writing, is designed to assist dentists in using Kerr's product. In the event of a defect in material or workmanship, Kerr's liability is limited, at Kerr's option, to replacement of the defective product or part thereof, or reimbursement of the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Kerr. In no event shall Kerr be liable for any indirect, incidental, or consequential damages.

EXCEPT AS EXPRESSLY PROVIDED ABOVE, THERE ARE NO WARRANTIES, BY KERR, EXPRESS OR IMPLIED, INCLUDING WARRANTIES WITH RESPECT TO DESCRIPTION, QUALITY, OR FITNESS FOR A PARTICULAR PURPOSE.

Manufactured By:  
**Kerr Corporation**  
1717 West Collins Avenue  
Orange, CA 92867  
TEL. (714) 516-7400  
FAX (714) 516-7635  
International TEL. (800) 537-7123  
International FAX (800) 537-7345  
Customer Service (800) KERR-123

European Union Representative:  
**Kerr Italia S.r.l.**  
Via Passanti, 332  
I-84018 Scafati (SA) Italy  
+39-081-850-8311

Kerr Distribution Facilities:  
**KerrHawe SA**  
Via Strecce 4  
CH-6934 Bioggio  
Switzerland  
+41-91-610-0505

**Kerr Australia Pty. Ltd.**  
Unit 10  
112-118 Talavera Road  
North Ryde 2113  
New South Wales, Australia  
+61-29-870-7500

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**Sybron Dental Specialties Japan, Inc.**