ENGLISH
Revolution Formula 2
 Flowable, light-cure hybrid resin restorative

Revolution Formula 2 is a flowable, light-cure hybrid resin restorative suitable for Class III, Class IV, and Class V restorations. Additional functions include: repair of enamel defects, repair of porcelain restorations, minor occlusal build-ups in non-stress bearing areas, pit and fissure sealant, cement for ceramic/composite veneers, incisal abrasions, and core build-ups. This unique composite which incorporates BIS-GMA chemistry along with a proprietary glass filler (approximately 60%) yields a flowable non-slumping restorative material. Revolution's combination of flowability and direct-application system simplifies material placement and minimizes finishing, which results in a consistently superior restoration.

PRIOR TO PLACEMENT – RECOMMENDATIONS ON PROPER POSITIONING

- Isolate tooth from adjacent teeth. Use clear, thin strips. Rubber dam is highly recommended.
- Moist dentin surfaces (rather than dry, desiccated dentin) have been shown to enhance the bond strengths of today's modern bonding agents (such as OptiBond All-In-One, OptiBond S, or OptiBond FL). NOTE: "Moist" does not infer blood or saliva contamination.
- Placement sequences such as those outlined in this instruction booklet should be followed closely.
- Take a moment to identify the surface to which you are bonding.
- Please take care to ensure that your air line is free of oil and other contaminants.

Directions for Use

1. Isolate and dry preparation.
 - a) Apply a self-etch bonding system such as OptiBond All-In-One.
 - b) If performing the total-etch technique, etch enamel and/or dentin for 15 seconds with Kerr Gel Etchant (37.5% phosphoric acid gel*). Rinse thoroughly, and lightly air dry. Apply a total-etch bonding system such as OptiBond S.
2. Select the appropriate shade. Dispense a layer of Revolution Formula 2 directly onto the preparation along the axial and gingival margins. Use a brushing action to allow the material to flow into the marginal interfaces and recesses of the preparation. Incremental layers should be no more than 2 mm. Light-cure each increment.
4. Build contours and ridges to desired shape by flowing the composite into the preparation with a brush-like motion.
5. Finish and polish using conventional techniques.
6. For a glaze finish, apply one or two drops of Revolution Formula 2, brush over the surface, and cure for 20 seconds**.

Shading System

Revolution Formula 2 is available in 16 Vita shades and 5 special shades as follows:

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Notes:

1. Replace syringe cap after each use to prevent the resin from hardening in the syringe.
2. Revolution Formula 2 dispensing tips are intended for single patient use only to prevent cross-contamination between patients.

Storage Conditions: Store Revolution Formula 2 at ambient temperature.

***CAUTION:** Contains phosphoric acid. Avoid contact with skin, eyes, and soft tissue. In case of contact with skin or eyes, flush immediately with water, get medical attention for eyes. Do not take internally.

**** Recommended Cure Times when using Kerr Demetron Lights:**

Demi/Demi Plus	5 seconds
L.E. Demetron II	5 seconds
L.E. Demetron I	10 seconds
Optilux 501	Boost Mode – 10 seconds Ramp Mode – 20 seconds Regular Mode – 20 seconds

For all other lights, see manufacturer's recommendation.

CAUTION: Uncured methacrylate resin may cause contact dermatitis and damage the pulp. Avoid contact with skin, eyes, and soft tissue. Wash thoroughly with water after contact.

Limited Warranty-Limitation of Kerr's liability

Kerr's technical advice, whether verbal or in writing, is designed to assist dentists in using Kerr's product. Such advice does not expand Kerr's limited warranty or relieve the dentist of testing Kerr's products to determine their suitability for the intended uses and procedures. The dentist assumes all risk and liability for damages arising out of the improper use of Kerr's product. In the event of a defect in material or workmanship, Kerr's liability is limited, at Kerr's option, to replacement of the defective product or part thereof, or reimbursement of the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Kerr. In no event shall Kerr be liable for any indirect, incidental or consequential damages.

EXCEPT AS EXPRESSLY PROVIDED ABOVE, THERE ARE NO WARRANTIES, BY KERR, EXPRESS OR IMPLIED, INCLUDING WARRANTIES WITH RESPECT TO DESCRIPTION, QUALITY OR FITNESS FOR A PARTICULAR PURPOSE.

Composite Material

Revolution Formula 2
 Flowable Light Cure Composite

Kerr Corporation
 1717 West Collins Avenue
 Orange, CA 92867
 U.S.A.
 TEL: (800) KERR 123
 TEL: +1-714-516-7400
 intl.kerrdental.com

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