User Manual Soft-Core® Endodontic Heater and Obturators



oft-Core

POSTED, 01 Jul 2010, G-DWG-REL

SybronEndo
A Division of Sybron Dental Specialties

³ 2: Heating of the obturator

he selected obturator in the Soft-Core heater and activate the heater, re information on how to use the heater, please see section 1 of this



3: Application of sealer

the obturator is heating, mix any heat resistant sealer with long g time (a non-eugenol sealer is recommended). Make sure that it canal is completely dry, before applying sealer. Place a very thin

ary. Utilize a paper point or a Soft-Core Size Verifier y a thin layer of sealer on the canal walls. When ing more than one canal, apply sealer in all at once. er on the wall of the canal. Too much sealer is neither desirable nor

Immediately insert the obturator into the canal to the working distance, without twisting the handle, and while using a firm and steady pressure.

Important: The obturator should be inserted to the full working length.

excess central core and gutta percha STEP 5: Removal of the handle,

removed once the gutta percha has hardened. The handle, excess core and gutta percha may be

1. For canals less than 18mm: Twist the handle and plastic core to one side to break off the excess.

2. For all Canals: Stabilize the plastic core small inverted cone bur. with a finger or cotton pliers, and cut away the handle and excess plastic core with a

movement of the core to minimize the chance of disturbing the apical seal.

Trim away the extra gutta percha using an endo excavator or similar instrument Important: Be careful to prevent vertical

THIN LAYER OF SEALER -

obturator into the canal STEP 4: Insertion of the

is ready for use, the obturator may be removed from the heater. When the oven indicates that the obturator







ction 1


nical Specifications

he Soft-Core Heater is intended to be used in heating Soft-Core obturators. The heater should not be used for any other purpose.

ut Power	100-240V AC, 50/60 hz
ver Supply	9V DC, 1.7A
ver Consumption (120V)	43.2 w
." Temperature	110 (degree) C
with	All Soft-Core Obturators
rranty	1 year



Temperature Status



- 1. Connect the power supply to the oven's DC inlet and

NOTE: When the oven is used for the first time, it may have a slight scent during the first use. This will not affect heater performance.

heating cycle will begin.

Product Precautions Read this entire manual before use

- Do not allow any persons to use this equipment This product should only be used in accordance with other purpose, will void the warranty and may cause directions in this manual. Use of this product for any
- with-out trained supervision.
- Do not reach inside the oven chamber while unit is
- plugged in or until after unit has cooled.

 Remove any obturator or gutta percha that is in contact with heater chamber walls or heating elements.

 Disinfect heater and sterilize obturator rest after each
- beating chamber.
 Do not allow obturators to warm in oven for more than Do not immerse unit in water or spray any liquid into procedure.
- Do not leave plugged in unit unattended.
- Do not attempt to service unit. Heaters should only be service by an authorized Soft-Core distributor.

Pre-heating the Heater


- an electrical outlet.
- When the unit is plugged in the oven temperature status will be red and timer light will be green.
 Once the heater is pre-heated, it will chime and the
- time should be no more than 15 minutes. temperature status light will turn green. Pre-heating

Stand-By Feature
To place the heater in stand-by mode: hold down the timer/stand-by button for 3 seconds. The heater will then go on stand-by and the heating elements will turn off. To re-activate the heater: if the heater is in stand-by, hold re-activate the heater. down the timer/stand-by button for 3 seconds and the

POSTED, 01 Jul 2010, G-DWG-REL

Soft-Core Endodontic Obturator

Core Endodontic Obturator is a core of biocompatible plastic, sized ing to ISO standards sizes 20 to 60, coated with thermoplastic gutta. The obturator sizes correspond to the ISO standards of root canal single obturator is all that is needed to totally obturate a root canal. sult will always be a perfectly placed root canal filling with a tight seal.



1. Plastic handle with a metal insertion pin

A Soft-Core Endodontic Obturator consists of three parts:

The injection molded plastic handle is color-coded according to the size of the obturator. The stainless steel insertion pin is 9 mm long, and is permanently attached to the plastic handle. The hollow portion of the plastic core surrounds 6 mm of the pin.

2. Outer layer of thermoplastic Natural GP™ gutta

make the gutta percha soft and adhesive. temperatures above approx. 100°C. When cooled it returns to a firm rigid state. Another heating will once again that it becomes soft and highly adhesive when heated to The gutta percha is reversible thermoplastic. This means

3. Plastic core with endo stop

The tapered plastic core is sufficiently flexible to negotiate curved canals.

A rubber endo stop is provided to indicate the working distance measurement. The plastic core and gutta percha are radiopaque. The Soft-Core plastic core is slightly oval-shaped. This secures adequate backflow of the heated gutta percha. The length of the core is 24 mm. The coronal 6 mm of the core is hollow to accommodate the metal insertion pin.