| | | |
|------|------------------------------|-----|
| (DE) | Gebrauchsanweisung | 3 |
| (GB) | Instructions for use | 11 |
| (FR) | Mode d'emploi | 18 |
| (ES) | Instrucciones de uso | 26 |
| (IT) | Istruzioni per l'uso | 34 |
| (PT) | Instruções de uso | 42 |
| (NL) | Gebruiksaanwijzing | 50 |
| (SE) | Bruksanvisning | 58 |
| (DK) | Brugervejledning | 66 |
| (NO) | Bruksanvisning | 74 |
| (FI) | Käyttöohjeet | 82 |
| (GR) | Οδηγίες χρήσης | 90 |
| (CZ) | Návod k použití | 98 |
| (HU) | Használati utasítás | 106 |
| (LV) | Lietošanas instrukcija | 114 |
| (LT) | Naudojimo instrukcija | 122 |
| (PL) | Instrukcja obsługi | 130 |

Produktbeschreibung

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Instructions for use 

Product description

iBond Total Etch is a light-cured single-component adhesive for use in adhesive restorative dentistry. iBond Total Etch has been developed for the adhesive bonding of plastic filling materials (e.g. composites, compomer, Polyglas[®]) to dental hard substances and laboratory-fabricated restorations (e.g. ceramics).

With iBond Total Etch, priming, bonding and desensitizing can be done in one step. iBond Total Etch is an ethanol based solution of light-activated, adhesive resins. Before using iBond Total Etch, the dental hard substance is conditioned with an etching gel (etch&rinse).

Composition

iBond Total Etch contains:

- Methacrylate
- Ethanol
- Fillers
- Photoinitiators
- Glutaraldehyde

Indications

- Adhesive fixing of direct composite, Polyglas® and compomer restorations
- Adhesive fixing of indirect laboratory-fabricated ceramic, Polyglas® and composite restorations (inlays, onlays, veneers and crowns)
- Treatment of hypersensitive tooth regions

Read the instructions for use carefully before use. Observe the safety instructions when using iBond Total Etch.

Application

1. Adhesive fixing of direct composite, Polyglas® and compomer restorations

1.1 Preparation

- The use of a rubber dam is recommended.
- Clean tooth with a toothpaste free from oil and fluoride.
- Prepare cavity as specified by the rules of adhesive filling treatment.
- Rinse prepared site with water and blow-dry.
- For deep cavities, we recommend protecting the dentine near the pulp with an appropriate base filling (e.g. calcium hydroxide preparation and glass ionomer cement base filling). Cover the dentine near the pulp only. Leave the remainder of the cavity free for adhesion of iBond Total Etch.

1.2 Conditioning

Prior to application of iBond® Total Etch to the cavity preparation, pre-treatment of the tooth structure with an adequate etching gel is mandatory (e.g. iBond Etch 35 Gel/Gluma Etch 35 Gel, iBond Etch 20 Gel/Gluma Etch 20 Gel by Kulzer).

1.3 Application of iBond Total Etch

- Place iBond Total Etch in a Dappen glass, pick it up with an applicator tip or a soft disposable brush and apply it generously to the complete surface of the cavity.
- Allow iBond Total Etch to set for 15 seconds.
- iBond Total Etch contains ethanol as a solvent. Carefully evaporate solvent and residual moisture under a gentle air stream until no movement of liquid can be detected.

An excessively strong air flow at the start of blow-drying will dilute the bonding and may result in inadequate adhesion.

- The surface must be visibly shiny after application of iBond Total Etch and after evaporation of the solvent. Make certain that the complete surface of the cavity is fully covered. If the cavity does not appear to be completely shiny all over, apply iBond Total Etch a second time as described above. *Make sure that no excess remains in corners of the cavity.* Polymerise iBOND® Total Etch for 20 s with a halogen or LED polymerisation lamp. Use a Kulzer Translux® polymerisation lamp or a polymerisation lamp with comparable intensity (min. 400–500 mW/cm²).

Insufficient lamp intensity will result in inadequate adhesion. Polymerisation lamps should be tested with reliable testing devices at regular intervals. The light emitting window should be placed as close as possible to the plastic surface during polymerisation.

- Then immediately apply and process the filling material in the cavity as directed by the manufacturer.

2. Adhesive fixing of indirect laboratory-fabricated ceramic, Polyglas® and composite restorations (inlays, onlays, veneers and crowns)

2.1 Preparation of the restoration

Prepare restoration surface as directed by the manufacturer's instructions for use.

2.2 Conditioning

Proceed as described in 1.

2.3 Application of iBond Total Etch

Make sure that no excess material is left in the corners of the cavity and that the iBond Total Etch coating is carefully blow-dried before polymerisation.

Otherwise proceed as described in 1.

- Polymerise iBOND® Total Etch for 20 s with a halogen or LED polymerisation lamp (see instructions above in 1.3).
- Apply resin cement (self-curing or dual-curing resin cement) as directed by the manufacturer's instructions for use.

3. Treatment of hypersensitive tooth regions and cleaning the tooth

- Carefully rinse the surface with water and blow-dry.
-

3.1 Conditioning hypersensitive tooth regions

Proceed as described in 1.

3.2 Application of iBond Total Etch

Proceed as described in 1.

- Polymerise iBOND® Total Etch for 20 s with a halogen or LED polymerisation lamp (see instructions above in 1.3).
- Carefully remove the oxygen inhibition layer with an alcohol-soaked pellet.

If the iBond Total Etch desensitisation effect is not sufficient, apply again as described above, cure and carefully remove the inhibition layer with an alcohol-soaked pellet.

Safety instructions

- iBond Total Etch contains ethanol. Ethanol is very flammable. Do not inhale vapour.
- iBond Total Etch contains glutaraldehyde. If it is inhaled or swallowed, it may damage health. Seek medical attention immediately if swallowed.
- iBond Total Etch contains methacrylate and glutaraldehyde. Contact with skin or mucous membranes may cause irritation or sensitisation. Prevent contact with skin, mucous membranes and eyes by taking precautions such as using a dental dam (patient), wearing rubber gloves (dentist) and wearing safety glasses (patient, dentist). Rinse eyes with water for several minutes after contact. Seek ophthalmological medical attention if pain persists. Wash skin with water and soap after contact.

Special instructions:

- For dental use only.
- Keep out of reach of children
- Materials containing eugenol may affect the polymerisation of iBond Total Etch
- Container must be at room temperature before use
- Do not use iBond Total Etch after the date of expiry

Side effects

Allergies to the product or its contents may occur in isolated cases. If an allergy is suspected, request a list of its contents from the manufacturer. In the event of a skin reaction or an allergy, do not use this product.

Storage

Store at temperatures of 4–25 °C (40–77 °F). The product can be stored in refrigerator. Use the adhesive at room temperature only. Do not expose to direct sunlight. The shelf life will be reduced by unsuitable storage conditions and the function of the product may be affected. Close the container tightly immediately after use.

Our instructions for use that are provided orally, in writing and by experiments are offered to the best of our knowledge. They are nevertheless non-binding, including in regard to third party trademark rights, and do not release the user from appropriately testing our delivered products for their suitability for the intended procedures and purposes. The products are used, applied and processed out beyond our control, and you are therefore exclusively responsible. We of course guarantee the highest quality of our products in accordance with our general terms of sales and delivery.

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Kulzer GmbH
Kulzer, LLC