## RC-Prep® Root Canal Preparation Cream

**Intended Use:**

*Premier® RC-Prep® Root Canal Preparation Cream is intended to help remove calcifications and lubricate the canal to permit more efficient instrumentation. RC Prep cleans and lubricates the canal so the instruments can cut freely without binding.*

**Description of Device and Contents**

RC-Prep Root Canal Preparation Cream is a water based paste, white in color, formulated to soften, lubricate, whiten and deodorize the dental canals. RC-Prep Root Canal Preparation Cream can be dispensed by the end user using a multitude of different methods.

RC-PREP IS DISPENSED IN:

- Jars of: 18g and 227g
- One box of: 2 prefilled syringes -9g each
- Dispensing Pump 18g
- MicroDose pipettes 0.2mL Qty: 100
- 3mL syringes Qty: 5

Note with Syringe use:

- Disposable tips are designed for single patient use only.
- Use a protective plastic sheath barrier on the syringe in order to mitigate potentia cross-contamination between patients.

**Contraindications:**

Do not use RC-Prep in cases of:

Known allergy to EDTA or any other contents of RC-Prep

Immediately wash skin with soap and water after contact or flush eyes if contact occur

**PRECAUTIONS:**

- For Jars and Pumps, dispense onto small clean surface, i.e. mixing pad to prevent cross-contamination
- Syringe tips shall be wiped down with disinfecting wipes prior to placing into a dispensing syringe.
- Use a protective plastic sheath barrier on the syringe in order to prevent potential cross-contamination between patients.
- Microdose pipettes are designed for single patient use only. Do not re-use
- Use rubber dam prior to use to prevent contact with intraoral tissues
- Always irrigate the canals and chamber with sodium hypochlorite solution until all evidence of bubbling has ceased.
- It is recommended to wear protective clothing when using this product.

**WARNINGS:**

- Do NOT** seal RC-Prep into the pulp chamber to open calcified canals. Residual organic tissues might react with the peroxide and cause pressure and discomfort.
- Do NOT** leave on tissues extracoronally or it may cause irritation
- Corrective dental procedures may be necessary if not all excess RC-Prep is removed.
- Pipettes and syringe tips are single use only.
- For use by a trained dental professional only
- Do not use after the shelf life period indicated on product labeling/packaging

**Adverse Reactions:**

Adverse reactions may include but are not limited to:

- Allergic reaction to RC-Prep
- Tissue damage from the prolonged exposure extracoronally or on intraoral tissues to RC-Prep
- Pain, discomfort, and/or abnormal sensations due to the presence of RC-Prep
- Infection from multi-use of pipettes and syringe tips

**Materials:**

Ethylene diamine Tetraacetic Acid (EDTA) and Urea Peroxide in a water soluble glycol base.

**PRIOR TO USE**

- For Jars**, dispense onto small clean surface, i.e. mixing pad to prevent cross-contamination
- For Pumps**, prime pump prior to use. Slowly press the tube dispenser tab several times to remove air or snap off bottom cap and gently push pump up. Replace bottom cap and begin use.
- For syringes**, disinfect syringe using an intermediate level disinfection process (liquid contact) as recommended by the Centers for Disease Control and use a protective plastic sheath barrier on the syringe. Ensure a new tip is placed on the dispensing end for each use.
- For Microdose**, microdose pipettes shall be wiped down with disinfecting wipes prior to placing into a dispensing syringe.

**DIRECTIONS FOR USE**

*(Please follow carefully)*

- Isolate tooth to be treated (e.g., rubber dam isolation)
- Open the pulp chamber to provide adequate access to the root canal(s).
- Remove chamber contents with sterile burs and excavators and flush with sodium hypochlorite solution.
- Fill pulp chamber with RC-Prep. **CAUTION:** if back pressure or resistance is felt, stop immediately as to not express material out the apex of the tooth.
- The sodium hypochlorite solution will react with the peroxide and the bubbling will help open the orifices of the canal. A fine instrument, reamer or file, is then introduced into the canal and using a slight pumping or mixing motion will enhance the bubbling and float out debris that is aspirated away.
- The first instrument should not be placed more than half the way down the canal. This action is then followed by gentle irrigation with sodium hypochlorite solution.
- The procedure is repeated adding RC-Prep to the pulp chamber.
- This time, the instrument is carried to approximately 3/4 the distance of the canal.
- When a vital pulp is present, it is not unusual to see the entire pulp float to the surface or be removed with the file.
- After gentle irrigation, the instrument is again inserted to approximately one millimeter short of the radiographic apex and an x-ray film taken to confirm the measurement.
- In teeth with two or more roots, the efferescant activity will tend to cleanse the pulp chamber and open the orifices of the canals, making it easier to find the entrance and easier to insert the root canal instruments.
- RC-Prep is used only with the first two or three sizes of instruments during debriement of the canal.
- There is usually enough residual activity of the RC-Prep to lift out the remaining pulp tissue and debris.
- Always irrigate the canals and chamber with sodium hypochlorite solution until all evidence **DO NOT** seal RC-Prep into the pulp chamber to open calcified canals. Residual organic tissues might react with the peroxide and cause pressure and discomfort.
- After preparation, thoroughly flush and evacuate canal using H2O to remove all irrigants.
- Dry canals with absorbent points and seal canals and chamber with medicaments of choice.

**Storage:**


Store in a cool place or at room temperature. Storage temperature is acceptable from 8°C/46°F to 25°C/77°F. Avoid excess heat. Tightly recap after use.

**Shelf Life:**

RC-Prep has a shelf life of 2 years from the date of manufacturing. The shelf life was determined based on the stability of the formulation of RC-Prep.

- Do not use after expiration date. See expiry date on the label.

**For Further Information:** Visit www.premierdentalco.com or call Premier at 610-239-6000.

Symbols	Symbol Meaning
	Expiry date
	Manufacturer
	Batch number
	Catalog number
	Consult IFU
	Caution: See instructions for use
	The upper and lower limits of temperature of use, storage and transportation
	Irritant (Skin and Eye) Skin Sensitizer. Acute Toxicity
<b>Rx ONLY</b>	*CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician

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