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Praxis system for finishing and polishing of composite restorations is composed by abrasive discs and pop-on mandrel, which provides agility at side inversions, as well as in disc replacements. Available in 3/8" and 1/2" sizes and four grits, codified by colors.

Use: Install the finishing/polishing device on the mandrel. Make sure that the disc is fully seated on the mandrel. Insert the mandrel shank into the chuck of the contra-angle handpiece. Perform a test simulating the use to verify the product's effective rotation. For better results, use the four grits, successively from the coarsest to the finest. On product must be applied dry, with low to medium speed and with intermittent touches. Discs of coarse and medium grits must be applied at 10,000 rpm rotation speed. Discs of fine and x-fine grits must be applied at 30,000 rpm rotation speed.

Warnings, precautions and restrictions: Professional use only. The discs are for single use only. Its reuse may cause cross contamination. The use of the mandrel or discs with similar products is of the user's sole responsibility. To avoid an accidental swallowing of this product, it is necessary the application of a total isolation technique of the operator's field. If there is proven allergy to any of its components, the product must not be used. Wear protective goggles to shield against particles (Dentists and patients). Wear a surgical mask to avoid inhalation of particles. Use of a rubber dam is highly recommended to avoid the potential ingestion of dental materials or other tissues. Rotating devices must be handled with care and control to avoid injury to the patient. Contact with soft tissues (skin, gingiva, mucosa) may result in injury to the tissue. The use of applied suction techniques applied during the procedure help limit the problem of particles discharging while using rotary instrumentation.

Connection with counter angle: Praxis mandrel is manufactured in compliance with specifications from the rule ISO 1797.

Cleaning and sterilization: Mandrel - Cleaning: Carry out the cleaning with the aid of water and enzymatic detergent 4 enzymes (Amylase, Protease, Lipase and Carbohydrase) and with pure pH range between 6.5 - 7.5. Prepare the dilution solution according to the manufacturer's instructions. The water should be potable. For best effectiveness, it is recommended that the temperature be between 40°C and 55°C. Submerge the product and leave in contact for 5 minutes with the solution. Remove parts and rinse under running water for at least 10 seconds. Dry the piece with absorbent paper or disposable wipes. Sterilization: The mandrel should be autoclaved (126°C, 16min), after each patient. Reprocessing Limits: Repeated processing has minimal effect on the mandrel due to its composition. This is a reusable device. The manufacturer guarantees up to 240 reuses. Inspection: Carefully inspect the mandrel to ensure that all debris has been removed. Observe damage / breakage of the product that would prevent proper operation and do not use in this case.

Composition: Discs: Aluminum oxide, polyethylene terephthalate, synthetic rubber resin, metal eyelets, wax-based pigments. Mandrel: Stainless steel.

Storage: Store below 30°C (86°F), in a dry place, protected from light.

Burgundy: Coarse | Pink: Medium
Light Pink: Fine | White: X-Fine

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EC REP Obelis S.A
TV DENTAL LTD