# Resinlay Pattern Photo Gel

PL

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ISO 13485 BPF

En

## Light-curing resin in gel for casting and carving patterns.

Resinlay Pattern Photo is a light-curing resin developed for general impression, which can be cured using any standard light-curing units. Its gel consistency makes it easy to apply, with a direct application using an applicator tip attached to the syringe. In addition, the resin has thixotropic properties, allowing the gel to be applied on a precise location. Besides that, due to its formulation, the resin allows full polymerization and burn, even of thicker layers.

**Indication:** This product is indicated for use by both, the dentist and the dental laboratory technician. Among its main indications are the manufacture of: Copings, Abutments, Pontics and bridges, Attachments, Veneers, Implant bars. Additionally, the product can be used for splinting procedure of transfer molding in cases of multiple implants.

**Use:** Prior to application, the area shall be isolated with the product Isolard TDV or glycerin. Remove the syringe cap and attach the applicator tip. Apply the resin directly on the area to be worked on. It is suggested to the application be performed by layers of up to 2mm thickness. In case of application of a very thick layer, it is recommended to apply thinner overlapped intermediate layers. The resin can be applied directly on the area to be worked on. If a thicker layer is to be applied, do not build it at once, build intermediate overlapped layers with smaller thickness. Before building the next layer, light cure the surface for 20 seconds, without removing the oxygen-inhibition layer. After the use, the cap should be replaced in the syringe, preventing the resin from polymerizing. It is recommended to remove the residual material from the applicator tip with a jet of compressed air, ensuring that the resin does not polymerize and obstructs the resin flow in further applications. Perform the finishing of the piece. The oxygen-inhibition layer can be removed by the finishing.

**Warnings, precautions and restrictions:** Professional use. The product must not be used in patients who show allergic reactions to any of its components. To prevent the product from polymerizing, close the syringe immediately after use. The packaging should be protected from heat and direct exposure to sunlight. The Material Safety Data Sheet for this product can be obtained by telephone. Disposal: Waste must be disposed of in an industrial or sanitary landfill, in compliance with local legislation in force. Do not dispose of waste in sewers, waterways or household waste bin. The product is a medical device. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**Sterilization:** The applicator tips are autoclavable (126°C, 147kPa, 16min). Before sterilization, execute the cleaning with the aid of water and enzymatic detergent. The sterilization can cause slight change in the coloring, but does not affect its functionality.

**Storage:** To be stored in a cool dry place, protected from sunlight. Do not store above 30°C.

**Composition:** Methacrylate and activators.

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**MD** Medical Device  
Dispositif Médical

**ISO 13485** **GMP**  
Good Manufacturing Practices

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**EC REP** Obelis S.A