

Quick-Stat™ - Hemostatic Agents

Descriptions:

Quick-Stat™ Liquid is an aqueous 15.5% ferric sulfate solution used for immediate and profound control of bleeding.

Quick-Stat™ Gel is a 20% ferric sulfate gel used for immediate and profound control of bleeding. It works on both hard and soft tissue.

Quick-Stat™ Free is a 25% Aluminum Chloride gel. It is ideal for cases in the esthetic zone as it does stain.

Indications for Use:

Quick-Stat™ Gel, Quick-Stat™ Liquid, and Quick-Stat™ Free is intended for sulcus retraction prior to impression making and to control bleeding and gingival oozing in restorative and operative dentistry used with gingival retraction cord. Quick-Stat™ Gel, Quick-Stat™ Liquid, and Quick-Stat™ Free facilitate the insertion of the cord into the sulcus.

Directions for Use - Preparation:

1.) For 1.2ml syringe

- Remove luer cap.
- Securely attach applicator tip (Stat-Flo™ tip).
- Verify flow prior to intraoral use.

2.) For 30mL Prefilled Syringe

- Remove luer cap from 30mL prefilled syringe.
- Attach one end of the docking port to the 30mL prefilled syringe and attach the other end of the docking port to a 1.2mL syringe.
- Slowly depress 30mL prefilled syringe plunger while guiding 1.2ml syringe plunger to desired fill.
- Separate syringes from the docking port and re-cap 30mL prefilled syringe.
- Securely attach applicator tip (Stat-Flo™ tip) to 1.2mL syringe.
- Verify flow prior to intraoral use.

3.) For 30ml bottle (Quick-Stat™ Liquid only)

- Place a small amount of Quick-Stat™ Liquid into a dappen dish or similar receptacle to soak cord prior to placement.

NOTE: For bottles, do not place syringes, cord, or instruments into bottle. Recap after dispensing liquid.

NOTE: To avoid cross contamination, use new syringes and tips for additional volumes.

Directions for Use - Impressions:

1.) Follow "Preparation" steps above.

2.) Slowly expel solution while simultaneously rubbing firmly against bleeding tissue.

Quick-Stat™ - Hemostatic Agents

- 3.) Continue expelling material and rubbing solution into the tissue to stop bleeding and remove coagulum buildup. Recommended contact time is 1-3 minutes.
- 4.) As coagulum appears, gently rinse away with air/water and use suction to remove. Hemostasis is achieved when new coagulum stops forming.
- 5.) When hemostasis is achieved, firmly rinse with air/water spray to ensure hemostasis is maintained.
- 6.) If bleeding continues, repeat steps 3-5 above.
- 7.) Displace tissue by packing appropriate size Stat-Pak™ Retraction cord or preferred retraction cord into sulcus using Vista-Pak™ Packers. For optimum displacement and hemostasis, place a small amount of Quick-Stat™ Gel, Quick-Stat™ Liquid, or Quick-Stat™ Free in a dish and soak cord prior to packing.
- 8.) For best results, clean preparation with chlorhexidine prior to the final impression.
- 9.) Wait 1-3 minutes before removing cord.
- 10.) Remove cord, rinse with a firm air/water spray, check for hemostasis, and make impression.

Directions for Use – Direct and Indirect Bonded Restorations:

- 1.) If adequate space is available to retract tissue during restorative procedure, follow “Impression” steps 1-7 above.
- 2.) Thoroughly rinse preparation and surrounding tissue using a firm air/water spray. Check to ensure area is free from blood. If bleeding reoccurs, apply more hemostatic agent to sulcus.
- 3.) Scrub preparation surface with pumice or chlorhexidine.
- 4.) If no residual bleeding occurs, wait 1-3 minutes and remove cord. If using an unsoaked cord for tissue displacement, cord may be left in place during restorative procedure to protect soft tissue or may be removed.
- 5.) Rinse again with a firm air/water spray and dry.
- 6.) For direct restorations, apply bonding agent and restorative per manufacturer’s instructions.
- 7.) For indirect restorations, proceed with cementing the final restoration.

NOTE: For definitive bonding and luting of restorations, thoroughly scouring and cleaning the preparations are crucial to avoid microleakage and decreased bond strength.

NOTE: Ideally, the preparation should be cleaned with a product such as Clean & Boost™ or CHX Restorative prior to bonding. If not thoroughly cleaned, residual hemostatic agent or coagulum may contaminate dentin and/or enamel potentially resulting in reduced bond strength and microleakage.

Precautions and warnings:

- For professional use only.
- Hemostatics can prevent quality adhesion and lead to staining microleakage under “direct bonded” restorations. They must be thoroughly cleaned with an air/water spray and pumice or chlorhexidine.

Quick-Stat™ - Hemostatic Agents

- For definitive bonding and luting of restorations, thoroughly scouring and cleaning the preparations are crucial to avoid microleakage and decreased bond strength.
- Hemostatics must be washed from the preparation site with an air/water spray to avoid staining.
- When using self-etch bonding agents, the tooth/preparation surface must be scrubbed with pumice or chlorhexidine and thoroughly washed. This is not necessary when a phosphoric etch bonding system is being used.
- To avoid blue discoloration, do NOT mix with epinephrine or any epinephrine-impregnated cord.
- Always thoroughly wash hemostatic from tissue before using other chemicals on the site.
- Do not use Quick-Stat™ Free on patients with known allergies to aluminum chloride or chemical sensitivities. Do not use Quick-Stat™ Gel or Quick-Stat™ Liquid on patients with known allergies to ferric sulfate or chemical sensitivities.
- Soft tissue discoloration is temporary.
- Prefilled syringes may be used several times as long as they are properly sealed with a syringe cap and disinfected after each usage.
- Do not allow product to be ingested.
- Use protective clothing and eyewear (patient and user) when loading and handling Quick-Stat™ Gel, Quick-Stat™ Liquid, or Quick-Stat™ Free.
- Avoid contact with eyes and skin.
- Dispose of used tips and empty syringes properly. Please refer to individual state laws to dispose of correctly.
- Use only recommended tips. To avoid cross contamination, do not re-use tips.
- Keep out of reach of children.
- Verify flow of all syringes prior to applying intraorally. If resistance is met, replace tip and re-check.

NOTE: For SDS please visit vistaapex.com.



Rx ONLY

REF

504600, 504640, 504650, 504660, 502804, 502800, 503550, 503555, 317120, 504155

 **Inter-Med, Inc.**
2200 South St.
Racine, WI 53404 USA
Tel: +1-262-636-9755

Made in USA
504600-I-V-ENG (7)

EC REP

CE
0459

Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands
Tel: +31(0)70 345 8570
Fax: +31(0)70 346 7299