Block Out Gel LC

Instructions for use

Block Out Gel LC is a light-curing methacrylate-based block-out material. Its blue contrast shade facilitates application on a model.

Not for intra-oral use!

Before application of **Block Out Gel LC**, plaster models should be left to dry for a minimum of 2 h to enhance adhesion. Drying can be speeded up with a warm airjet (hairdryer) or by placing it in a dryer.

The block-out material can be applied directly from the syringe.

To polymerize the material place the model in a light-box for 2 min or light-cure each blocked-out area for 10 - 30 s with a halogen polymerization device.

After light-curing remove the inhibition layer, best by wiping it off with a tissue soaked in alcohol.

Close syringe tightly after use and store in a dark and cool area.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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VOCO GmbH