Grandio Seal

in accordance with DIN EN ISO 6874

Carefully read instructions prior to use

Product description:

Grandio® Seal is a light-curing, highly filled, fluoride-containing nanohybrid fissure sealant with outstanding mechanical properties.

Grandio® Seal contains 70% w/w inorganic fillers in a methacrylate matrix (Bis-GMA, TEGDMA). The sealant is light-curing.

Indications:

- sealing/filling of pits and fissures
- sealing/facing of damaged enamel surfaces
- covering of caries predilection sites during orthodontic treatment
- sealing of deciduous teeth

Contraindications:

Grandio® Seal contains BIS-GMA, TEGDMA and BHT. Do not use in cases of known allergies to these ingredients.

Application:

Fissures should be sealed as early as possible, as soon as the eruption of the tooth allows treatment and moisture control is possible. After isolating the operating site, thoroughly clean the surfaces to be sealed (rotating brush/cleaning paste free of fluorides), spray and dry.

In case of incipient carious lesions or unclear diagnosis, extend the fissure with a fine instrument.

Etching

Apply etching gel, work in with a short-haired disposable brush and let it act for 15 - 30 sec according to the condition of the enamel. Thoroughly spray and dry. The etched enamel should have an opaque and matt-white appearance (contamination should be avoided).

Follow the instructions for use of the etching gel.

Application

Remove the protective cap from the syringe. Screw a new cannula onto the **Grandio® Seal** syringe. In order to ensure the function of the non-dripping syringe do not pull back the plunger during or after use.

Grandio® Seal flows homogeneously from the cannula upon medium pressure and thus is very easy to dispense accurately. Check the cannula for clogging if stronger pressure becomes necessary.

Apply **Grandio® Seal** from the margin into the fissure or cavity and allow it to penetrate for 15 - 20 sec. Remove any air bubbles with the fine cannula.

Light-curing

Polymerize **Grandio® Seal** for 20 sec per surface and check the sealing. Conventional polymerization devices are suitable for light-curing this material. The light output should be at least 500 mW/cm² on halogen polymerization devices and 300 mW/cm² on LED devices. Re-seal undercuts, remove excess cautiously (e.g. finishing diamond), polish if necessary (rubber polisher). Then apply a fluoride varnish on the treated and neighboring teeth.

Storage:

Remove the cannula after use and carefully close the syringe immediately. Store **Grandio® Seal** at temperatures between 39°F - 73°F (4°C - 23°C). Avoid unnecessary exposure to light.

Information/precautionary measures:

- An operating site which cannot be kept sufficiently dry involves the risk of a (partial) premature loss of the sealing. In such cases, the application of a fluoride varnish or sealing with a glass ionomer cement is recommended until perfect sealing is possible.
- If the etched surface is contaminated (e.g. by saliva), the adhesion of **Grandio® Seal** may be reduced. In case of contamination clean, etch again for approx. 10 - 20 sec, spray and dry.
- Adhesion can be increased by using a bonding agent after etching (e.g. in case of difficult moisture control/extended sealing). Follow the bonding agent manufacturer's instructions for use. Then apply **Grandio® Seal**.
- Only expose **Grandio® Seal** to the operating light for a short time as otherwise polymerization will start (viscosity increases).
- Check the sulcus after application and remove residues of material.
- Avoid contact with the mucosa as far as possible. Do not use cleaning or drying agents containing phenolic substances (e.g. eugenol, thymol) as they inhibit complete polymerization.

- Carious lesions (e.g. around fillings) can be stopped by appropriate sealing if renewal of the filling is not desired or not possible. For this purpose, clean the areas affected, prepare as minimally as possible and etch in the enamel area. In such cases, sealing of the complete filling is recommended.
- In the case of orthodontic treatment in areas of reduced oral hygiene, sealing of these areas can provide good and permanent protection. Fluoridation should be continued in all cases to protect the unsealed tooth areas from caries.
- Do not use **Grandio® Seal** after the expiration date.

PRODUCT ORDERING INFORMATION:

Syringes 5 x 2.0gm, Voccid 5ml, application cannulae type 45	REF 1060
Syringes 2 x 2.0gm, application cannulae type 45	REF 1061
Application cannulae type 45, 100 pieces	REF 2147

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. VOCO recognizes its responsibility to replace products if proven to be defective. VOCO does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

Descriptions and data constitute no warranty of attributes and are not binding.

CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.

No person is authorized to provide any information which deviates from the information provided the instructions for use.

For questions or comments, please call 1-888-658-2584.

Keep this material out of reach of children.

For dental use only.

An explanation of the symbols used in labeling can be found at www.voco.dental/symbols

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VOCO GmbH Phone +49 (4721) 719-0