

# **TECHNICAL BULLETIN**

VIRACLEAN® October 2021 1



HOSPITAL GRADE DISINFECTANT AUST L 69000

## **BENEFITS:**

- Is non flammable
- Is suitable for use on healthcare surfaces
- Does not contain chlorine or oxidising biocides
- Will not cause staining on fabrics or carpets
- Is pH neutral
- Non corrosive to metals
- Has a mild lemon fragrance

VIRACLEAN® is the result of years of intensive research and development. A considerable investment in both time and money has been made in the interests of improving the safety and efficacy of hospital grade disinfectants. VIRACLEAN® can be used for both cleaning and disinfecting, demonstrating beneficial versatility for healthcare and associated environs.

## MICROBIOLOGICALLY TESTED:

VIRACLEAN® is proven to kill:

- Acinetobacter
- Candida albicans
- Coronaviruses including SARS- CoV-2 (COVID-19)
- Enterococcus faecalis (VRE)
- Escherichia coli (E coli)
- Hepatitis B group virus
- Herpes Simplex virus
- Influenza virus
- Klebsiella pneumoniae (CPE/CRE)
- Proteus vulgaris
- Pseudomonas aeruginosa
- Salmonella choleraesuis
- Staphylococcus aureus (MRSA or Golden Staph)

Independent test reports available on request.

# **PRODUCT DESCRIPTION:**

VIRACLEAN® is a pink liquid with a lemon fragrance. It is easy to use and is conveniently packaged in squeeze or spray bottles, 5 litre containers, and 15 litre containers. VIRACLEAN® is formulated to replace most cleaning products and hospital grade disinfectants. VIRACLEAN® can be used for cleaning and disinfecting. This product is not intended for use on therapeutic devices. INSTRUMAX® PINK is available as an instrument disinfectant low level for use in disinfecting non critical medical devices.

ABN No. 24 000 906 678 PO Box 1076, North Sydney NSW 2059 Product Support Hotline: 1800 833 566 www.whiteley.com.au



Protecting lives and livelihoods since 1933



# **TECHNICAL BULLETIN**

VIRACLEAN® October 2021 2

#### CHEMICAL COMPOSITION AND PROPERTIES:

Surfactant Blend:	Blended
Active Ingredient:	Benzalkonium Chloride 0.42% w/v
Colour :	Light Pink
pH:	Neutral (7)
Specific Gravity:	0.99
Dangerous Goods Rating:	Not applicable
Shelf Life:	18 months
Safety Equipment:	Use of gloves and eye protection recommended

### **DIRECTIONS FOR USE:**

The use of gloves is recommended whenever using disinfectants.

To remove organic soiling (blood, sweat, sputum etc) from a surface, it is recommended to apply VIRACLEAN® liberally around the soiled area and work towards the centre. Soak up using fresh, disposable cloth. Once clean, to disinfect the surface, apply undiluted and leave to dry for 10 minutes at room temperature. If required, soak up using fresh, disposable towelling. Used towelling should then be disposed to a suitable receptacle.

#### MATERIALS COMPATIBILITY:

VIRACLEAN® has been used on a vast range of materials over several decades and has excellent general materials compatibility. However it has not been tested on every type of plastic or material. Should you notice degradation, cracking or a change in the colour of materials please immediately discontinue use of VIRACLEAN® on that material.

### EDUCATIONAL VIDEO AVAILABLE:

An educational video is available on 'How to clean and disinfect blood spills using VIRACLEAN®'. To watch the educational video visit www.whiteley.com.au/education.

ENVIRONMENT	PRODUCT USE
Nursery/Childcare:	Cleaning cribs, cots, change tables.
CSSD	Bench tops, stainless steel surfaces, sinks, trolleys, tap fittings.
Bathroom (toilet/shower)	Sinks, tile surfaces, and all hard surfaces.
Intensive Care Units	Cleaning and disinfection of all hard surfaces.
Manual Cleaning of Theatre Equipment	VIRACLEAN <sup>®</sup> can be used to clean and disinfect equipment and cables and will not adversely affect surfaces.
Operating Theatres	Cleaning and disinfecting of tables in theatre.
Blood Spills (excludes medical devices)	VIRACLEAN <sup>®</sup> is ideal for cleaning and disinfecting blood spills on all hard surfaces.
Telephones	VIRACLEAN <sup>®</sup> is ideal for cleaning and disinfecting phones where large numbers of bacteria are found.
Dental & GP Surgeries (excludes medical devices)	VIRACLEAN® can be used in Dental & GP Surgeries for cleaning and disinfecting surfaces including chairs, stools, tables, benches and light fittings.
Laboratories & Pathology Clinics	VIRACLEAN® is excellent for wiping down benches, disinfection of white rooms, cleaning and disinfection of fume cupboards, etc.

ABN No. 24 000 906 678 PO Box 1076, North Sydney NSW 2059 Product Support Hotline: 1800 833 566 www.whiteley.com.au



Protecting lives and livelihoods since 1933



VIRACLEAN® October 2021 3

SAFETY DATA: Information is listed in the Safety Data Sheet.

PACK SIZE: 12 x 500mL Squeeze

12 x 500mL Surface Spray 2 x 5 Litre 15 Litre

#### HEALTH AND SAFETY STATEMENT:

There is no way of predicting an individual's response to cleaning or disinfecting products. All cleaning and disinfecting products have a potential to cause skin and other irritation in susceptible individuals. Direct physical contact with the skin should be avoided. If irritation does occur when using this product or any other cleaning/disinfecting solution then the individual should be removed from direct contact with the solution until a satisfactory safe working environment for that individual is determined.

### **ENVIRONMENTAL STATEMENT:**

Whiteley manufactures products that comply with Australian and International environmental policy. Whiteley aims to ensure that every product is 100% biodegradable, whilst achieving optimal performance and benefit to customers. All Whiteley products are supplied in 100% recyclable containers. Whiteley and its staff are fully committed to ensuring that our environment and those international environs in which our products are used are enhanced by the use of technology, and that these environments are preserved for the enjoyment of future generations.

© 2021 This document remains the property of Whiteley Corporation Pty Ltd.

Alterations are not permitted without prior written authorisation from Whiteley Corporation Pty Ltd.

ABN No. 24 000 906 678 PO Box 1076, North Sydney NSW 2059 Product Support Hotline: 1800 833 566 www.whiteley.com.au



Protecting lives and livelihoods since 1933