# LuxaCrown

## Inhalt | Contents

Gebrauchsinformation - DEUTSCH	4
Instructions for use - ENGLISH	11
Mode d'emploi - FRANÇAIS	18
Istruzioni d'uso - ITALIANO	25
Instrucciones de empleo - ESPAÑOL	32
Instruções de uso - PORTUGUÊS	39
Gebruiksaanwijzing - NEDERLANDS	46
Brugsanvisning - DANSK	53
Användarinstruktioner - SVENSKA	60
Instrukcja użycia - POLSKI	67
Информация по применению - РУССКИЙ	74
Kullanım kılavuzu - TÜRKÇE	81
ข้อแนะนำการใช้ - ภาษาไทย	88
Petunjuk penggunaan - BAHASA INDONESIA	95

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## Product description

---

LuxaCrown is a self-curing composite for the fabrication of semi-permanent crowns and bridges that are cemented with a permanent or temporary luting cement. The material is mixed automatically.

## Indications

---

Crown restorations:

- for restoration of the anatomical form in order to provide durable protection for the remaining tooth;
- in order to restore chewing function;
- for esthetic corrections.

Bridge restorations:

- in order to restore chewing function;
- for esthetic corrections.

## Contraindications

---

Do not use the material:

- on preparations with a minimum wall thickness of  $< 1.5$  mm;
- on tooth cores with a height of  $< 3$  mm;
- in cases of allergies to any of the ingredients or, as can occur in rare cases, in the event of contact allergies.

## Side effects

---

There are no known systemic side effects to date.

## Interactions

---

Leftover methacrylate-based materials may alter the setting behavior of the silicone impression and bite registration materials.

- ▶ If a restoration was prepared beforehand, carefully clean the preparation and the surrounding tissue with alcohol-soaked gauze and water spray before taking the impression.

## Basic safety information

---

- For dental use only!
- Keep out of the reach of children!
- Avoid contact with the skin! In the event of accidental skin contact, immediately wash the affected area thoroughly with soap and water.
- Avoid eye contact! In the event of accidental contact with the eyes, immediately rinse thoroughly with plenty of water and consult a physician if necessary.

## Timing

---

0:00–0:40 min*	Place into the mouth
1:45 –2:20 min	Remove from the mouth
5:00 min*	Setting completed Process restoration

\*The indicated times apply at a room temperature of 23°C/73°F and a normal relative humidity of 50%. Higher temperatures shorten, lower temperatures increase these times.

## Recommended use

---

1. Before preparation of crowns or bridges, take an impression using alginate, silicone or a thermoplastic impression material. Alternatively, the temporary can be made using a laboratory-prepared thermoforming matrix. Strictly follow the appropriate manufacturer's instructions.

**Note:** In a silicone impression, undercuts should be compensated for and, if necessary, relief channels should be cut.

2. Dry the prepared teeth and compensate for any undercuts present in the preparation.
3. Lightly coat the prepared teeth, surrounding tissue and any present composite core build-ups with Vaseline or a similar separating agent.

**Note:** The convergence angle must be kept as small as possible in the preparation. Ensure an appropriate connector cross-section with splinted crowns and bridges.

4. Select the desired LuxaCrown shade.
5. Fill the impression with LuxaCrown. Start by applying the material to the occlusal surfaces and fill up working toward the gingival areas. Keep excess material to a minimum.

**Note:** To avoid bubbles, keep the mixing tip immersed in the material throughout application.

6. After a maximum time period of 0:40 minutes from the start of mixing, place the filled impression with moderate pressure on the prepared teeth and fix.
7. Monitor the setting process intraorally (e.g. with a probe), as the curing process is influenced by the temperature in the mouth and as it is recommended to take off the restoration during the elastic phase.

8. Remove the restoration from the mouth during the elastic phase, i.e. in the time period from 1:45 to 2:20 minutes after the start of mixing.
9. Remove the oxygen inhibition layer at the surface of the restoration (e.g. with alcohol).

**Caution!** Grinding dust can get in the eyes and airways when working on or finishing the restoration.

- ▶ Wear a face mask and protective goggles.
  - ▶ Make sure that grinding dust is properly extracted – do not breathe it in.
10. About 5:00 min after the start of mixing, finish the restoration using cross-cut burs, flexible discs, etc.
  11. If necessary, coat the surface of the restoration with a light-curing high-gloss glaze (e.g. Luxatemp-Glaze & Bond) or manually polish it with a rubber polisher or

small goat hair brush and a suitable polishing paste.

12. Thoroughly clean the preparation with water spray, dry it and thoroughly remove the separating medium. In doing so, ensure that the dentine retains a little residual moisture and does not become too dry.
13. Slightly roughen the inner crown surfaces with a bur or sand blaster.
14. Cement the restoration using a suitable permanent or temporary cement (e.g. PermaCem 2.0, TempoCem). Observe the manufacturer's instructions.

**Note:** Cements that contain eugenol may prevent the polymerization of resin-based luting cements during the definitive cementation.

15. Check the occlusion and adjust accordingly.

## Individualization of and repairs to restorations

---

**Note:** Repairs can only be performed in cases where small defects are present that do not adversely affect stability.

**Caution!** Grinding dust can get in the eyes and airways when working on or finishing the restoration.


- ▶ Wear a face mask and protective goggles.
- ▶ Make sure that grinding dust is properly extracted – do not breathe it in.

1. Individualization of or repairs to in site restorations: Slightly roughen the restoration with a bur or sand blaster at the corresponding site. Afterwards apply a bonding agent (e.g. Luxatemp-Glaze & Bond) in accordance with the manufacturer's instructions.
2. Apply LuxaCrown to the prepared sites.
3. Finish the corresponding site in the usual manner after 5 min.

**Note:** In order to accelerate polymerization, the finished restoration can be placed in hot water at a temperature of 50°C for a few minutes.

Alternatively, the restoration can also be individualized or repaired using the flowable and light-curing composite EcuSphere-Flow. Observe the manufacturer's instructions.

## Using the Automix cartridge

### Inserting the cartridge

1. Push up the lever [A] at the back of the Automix-Dispenser and pull back the slide [B] all the way.

2. Lift up the plastic lever [C], insert the cartridge and lock it using the plastic lever [C].

**Note:** Make sure that the notches on the cartridge line up with those on the Automix-Dispenser.

### Attaching a mixing tip

1. Remove the cartridge cap [D] or used mixing tip after turning it 90° counterclockwise, and discard it.
2. Attach a new mixing tip [E].

**Note:** In order to achieve an optimum mixing result, DMG recommends the use of the mixing tips available from DMG. All mixing tips are also available in refill packs.

**Note:** Make sure that the notches on the mixing tip line up with those on the cartridge.

3. Twist the mixing tip clockwise through 90° to lock it.

## Applying material

- ▶ When using a cartridge for the first time, squeeze out an approximately pea-sized amount of material and discard it.
- ▶ Mix the material in the mixing tip using the lever [F] on the Automix-Dispenser and apply directly.

**Note:** After application, leave the used mixing tip on the cartridge as a seal!

## Removing the cartridge

1. Push up the lever [A] at the back of the Automix-Dispenser and pull back the slide [B] all the way.
2. Lift up the plastic lever [C] and remove the cartridge.

## Disinfection of the Automix-Dispenser

The Automix-Dispenser can be reused after disinfection, although it should be replaced if there are visible signs of wear.

Disinfection can be carried out using a commercially available dipping disinfectant. DMG recommends using only disinfectant solutions listed by the RKI (Robert Koch Institute). The Automix-Dispenser is also autoclavable.

## Composition

Glass filler material in a matrix of multifunctional methacrylates; catalysts, stabilizers, and additives. Does not contain methyl methacrylate. Filler content: 46 wt.% = 26 vol.% (0.02 to 1.5  $\mu\text{m}$ )

## Storage and shelf life

- ▶ Store at room temperature (15–25°C/59–77°F) in a dry place.
- ▶ Keep away from direct sunlight.
- ▶ Do not use after the expiration date!

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**DMG** Chemisch-Pharmazeutische Fabrik GmbH