
	Ø	%	L	Stop	Ring	RPM	Recommended Torque (N.cm)
Revo-S+ SC1	25	6	21	●	○	250 - 400	1.8 - 2.5
Revo-S+ SC2	25	4	21 - 25 - 31	●	●	250 - 400	1.8 - 2.5
Revo-S+ SU	25	6	21 - 25 - 31	●	●	250 - 400	1.8 - 2.5

	Ø	%	L	Stop	Ring	RPM	Recommended Torque (N.cm)
Revo-S+ SC1	25	6	21	●	○	250 - 400	1.8 - 2.5
Revo-S+ SC2	25	4	21 - 25 - 31	●	●	250 - 400	1.8 - 2.5
Revo-S+ SU	25	6	21 - 25 - 31	●	●	250 - 400	1.8 - 2.5

1. Indication

Revo-S+: Coronary flaring and root canals shaping during the non-surgical endodontic treatment.

For use by dental professionals only.

2. Contraindications

There is no contraindication to the use of Revo-S+ for endodontically treating a tooth by orthograde way.

3. Complications

In cases of complex canal anatomy, per-operative risks (instrumental breakage, ledge, stripping, zipping, false path, perforation, etc.) may occur and lead to a risk of infectious processes.

4. Warnings/precautions

- Keep presentUDI and packaging information until the last use.
- Instruments supplied non-sterile: Follow the reprocessing instructions (§8) before use.
- Use in accordance with the recommended technique (§7).
- Be sure to observe the maximum recommended number of uses.
- Do not use for endodontic retreatment.
- Use in continuous rotation at the recommended speed.
- The decision to use the endodontic instrument must be referred to a clinical case expert, particularly if a root canal anatomy is considered to be too complex.
- In cases of complex anatomy, the maximum use torque may be different from the manufacturer's recommended value.
- Respect the recommendations for patients identified with a risk of infectious endocarditis.
- Contains nickel and titanium: should not be used in persons with a known allergy to these metals.
- Respect good dental practices, in particular the use of a dental dam and gloves.
- Check the condition of the instrument before use between each canal. If the instrument is damaged or shows signs of wear, do not use it.
- Do not use the instrument if there is any doubt about compliance with the storage conditions.
- Do not use the instrument if there is any doubt about the identification of the product.
- Year of manufacture: see labelling.
- Inform the manufacturer and the regulatory authority in the country of any serious incident related to the instrument.

5. Clinical claims

Clinical performance in normal conditions:

- Shaping ability (debris removal and cutting efficiency)
- Respect of canal anatomy
- Geometry adapted to irrigation
- Maintains the apical foramen in its original position

Safety, in normal conditions:

- Limits the risk of fracture
- Limits the screwing effect
- Limits the risk of shaping aberration

6. Characteristics

Type 2 rotating instruments under EN ISO 3630-1:2019 standard:

- Maximum recommended torque (not required): see table 'Recommended Torque'.

Material of the operative part: Nitinol.

Use in connection with an endodontic contra-angle according to EN ISO 1797:2017 (Type 1).

Number of uses: maximum of 5 canals recommended if the instrument is not visibly damaged.

Medical device class according to Directive 93/42/EEC and MDR 2017/745: IIa.

7. Protocol

Operating procedure:

Where unconstrained penetration of a manual file 15 is possible up to apical constriction upon initial penetration, using Revo-S+ SC1 is unnecessary and Revo-S+ SC2 should be used directly.

- Revo-S+ SC1 is used with a "wave" motion using three successive movements, from coronary to apical, centered in the canal, allowing progression of a few millimeters or until a constraint is felt. At each end of progression, a series of brushing movements can be performed. The initial canal portions - cervical and middle - are thus successively shaped using a "crown-down" approach.
- Revo-S+ SC2 and Revo-S+ SU are used with a so-called "wave" motion using three successive movements, from coronary to apical, centered in the canal, allowing progression of a few millimeters. When the instrument reaches the apical limit of preparation, a series of selective ascent brushing movements should be performed to finalise cleaning and shaping.

Irrigate abundantly between each step and each passage of the instrument.

8. Reprocessing instructions

General recommendations	- For all metal devices, the use of anticorrosion disinfecting and cleaning agents is recommended. - For your own safety, please wear personal protective equipment (gloves, glasses and mask). - Do not use cleaning or disinfection agents containing phenol, aldehyde and alkaline composition. - Always follow the instructions for use provided by the manufacturer of the products.
Limitations on reprocessing	- Due to the product design and the materials used, no definite limit to the maximum number of performable processing cycles can be specified. The service life of medical devices is determined by their function and careful handling. Multiple use disinfection and re-sterilisation cycles may lead to increased risk of file separation. - The user must ensure that the processing method used, including resources, materials and personnel, is appropriate and meets the applicable requirements. - The state of the art and national laws require that validated procedures be followed.
Material needed	- Gloves, masks, gown as recommended by the manufacturer of cleaning agent and detergent - Tap or deionized water - Disinfectant (neodisher® Septo Active) - Detergent (neodisher® MediZym) - Small soft brushes - Container - Ultrasonic tub or washer-disinfector - Class B sterilisation apparatus Remark: All material used should be cleaned and replaced regularly. Identify material used for each step of the process (initial treatment, cleaning or rinsing).

1	Initial treatment Place used products in a container or a wipe with tap water at 20-40 °C and 1.0 % neodisher® Septo Active for 5 to 15 min. Rinse the products with tap water at 20-40 °C for 1 min. If there is waiting time before the next step, make sure the device stays moist by placing it in a wet wipe. Do not exceed 1 hour waiting time. Remarks: - Do not use fixing agents or hot water (> 40 °C), as this causes fixation of residues and can impair successful cleaning. - Follow instructions and respect concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or other defects on devices).
2	Preparation before cleaning If the devices have visible impurities, manually pre-clean is recommended by brushing under tap water at 20-40°C for at least 1 min with a soft brush until all impurities have been removed. Remark: Follow instructions and respect concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or other defects on devices).
3	Visual inspection Inspect the used products and discard damaged products (broken, untwisted or abnormally bent).

manual	4	Cleaning Insert the products into an ultrasound apparatus beaker. Run ultrasound apparatus for 10-30 min with tap water and 0.5-2.0% neodisher® MediZym. Remarks: Follow instructions, observe water quality, concentrations and cleaning time specified by the manufacturer of the cleaning solution.
	5	Rinsing Rinse the products with tap water at 20-40°C for 1 min. Remark: It is recommended to use deionized water.
	6	Drying Dry the products with compressed air until products are visibly dry.
automatic	4	Cleaning/Rinsing/Drying Place the instruments in the tray of the washer/disinfector's sliding trolley. Run cleaning cycle with 0.2-1.0% neodisher® MediZym. Perform drying. Remarks: - Disinfection (thermal or chemical-thermal) is not required as the products are sterilised after cleaning. - Follow the instructions and concentrations specified by the manufacturer of the detergent solution. - Follow the instructions of the washer-disinfector and check the success criteria have been met after each cycle as specified by the manufacturer. - The final rinse step should be with deionized water. For other steps, follow the water quality defined by the manufacturer. - Use only approved washers-disinfectors that comply with EN ISO 15883, and are maintained and validated regularly.
	5	
	6	
	7	Visual inspection Inspect the used products. Repeat steps 4, 5 and 6 if the product is not visibly clean or discard any damaged products.
	8	Packaging Place the instruments in a paper-plastic pouch for steam sterilisation in compliance with ISO 11607 and EN 868 standards. Remarks: - For sharp devices that are not contained within a box, silicon tubes should be placed around the devices to prevent the packaging from being pierced. - Seal the pouches as recommended by the pouch manufacturer. If a thermo-sealer is used, the process must be validated and the thermo-sealer must be calibrated and qualified.
	9	Sterilisation Sterilise the products using steam: - Apparatus: class B - Minimum temperature: 132°C - Minimum time: 3 min - Absolute pressure: 2.2 bar - Minimum drying time: 20 min Check the physico-chemical indicators and cycle parameters. The temperature settings at 134°C and duration of 18 min are mandatory for prion inactivation according to French regulations. Remarks: - When sterilising multiple instruments in one autoclave cycle, ensure that the steriliser's maximum load is not exceeded. - Place the pouches in the steam steriliser according to the recommendations specified by the steriliser manufacturer. - Use only a pre-vacuum air removal steam steriliser that meets the requirements of EN 13060 (class B, small steriliser) and EN 285 (full-size steriliser), with saturated steam.
	10	Storage Store the products in a dry, clean and dust-free environment at the temperature specified by the paper-plastic pouch by the steam steriliser manufacturer. Remark: Check the packaging and the medical devices before using them (packaging integrity, no humidity and use-by date). In case of damage, a complete rerun should be performed.

9. Storage and transport conditions

There are no special storage and transport conditions for non-sterile devices.


10. Disposal

After use, instruments must be placed in a secure container, used to collect cutting or sticking instruments (like needles or disposable bistouries) as per good dentistry practice.

11. Symbols

	NiTi raw material		Root canal preparation
	Quantity		Assortment
 CLASSICS	'CLASSICS' handle		Medical device
	Sterilizable in a steam sterilizer (autoclave) at the temperature specified		

Year of CE marking: 2021
 Date of IFU revision : 2021-06-18
 Reference: 30005409-C

	Ø	%	L	Stop	Ring	RPM	Recommended Torque (N.cm)
Revo-S+ SC1	25	6	21	●	○	250 - 400	1,8 - 2,5
Revo-S+ SC2	25	4	21 - 25 - 31	●	●	250 - 400	1,8 - 2,5
Revo-S+ SU	25	6	21 - 25 - 31	●	●	250 - 400	1,8 - 2,5

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