3M Filtek[™]

One Bulk Fill Restorative

You choose it for speed. Now it's more aesthetic, too.



You'll like that it's simple. They'll be happy that it looks good, too.

3M[™] Filtek[™] One Bulk Fill Restorative

It's true. Now bulk fill placement can be more aesthetic.

Posterior restorations are challenging—which is why bulk fill composites that simplify and speed up the procedure have become so popular.

However, until now, the faster, easier placement of bulk fill has been at the sacrifice of aesthetics. So, as a global leader in aesthetic and restorative technologies, we went back to our experts for a solution you're sure to love:

3M[™] Filtek[™] One Bulk Fill Restorative.

Designed for the posterior, its unique optical properties and improved opacity give you the simplicity of one-step placement—up to 5 mm—*without compromising aesthetic results*.

Now you and your patients have something to smile about.

Fast and easy one-step placement

M Filtek

- Increased opacity for improved aesthetics
- Stress relief to enable up to 5mm depth of cure
- Excellent adaptation
- No need for layering or expensive dispensing devices
- Excellent handling and sculptability
- 3M's true nanotechnology provides excellent wear resistance and brilliant polish retention
- High radiopacity
- Unique tip design for easier access to deep cavities
- Available in 5 shades

Up to 5 mm

How 3M Science enables the simplicity of one-step bulk placement with improved aesthetics.

The technologies behind 3M[™] Filtek[™] One Bulk Fill Restorative

Smart Contrast Ratio Management

Increases opacity without sacrificing depth of cure


Further development of 3M nanotechnology allows us to manage opacity and depth of cure in a smarter way. With our unique Smart Contrast Ratio Management technology, we've increased opacity without reducing depth of cure.

In its cured state, 3M[™] Filtek[™] One Bulk Fill Restorative has a higher opacity than 3M[™] Filtek[™] Bulk Fill Posterior Restorative, resulting in improved aesthetics.

Did you know?

Contrast Ratio is the measurement of the translucency or opacity of the material.

Contrast Ratio




*The contrast ratio is the average of all shades.

50K magnification

Innovative AUDMA and AFM monomers


Reduces shrinkage and relieves stress to allow bulk fill placement up to 5 mm.



3M's TRUE nanofiller technology

Provides excellent wear resistance and brilliant polish retention.

5K magnification of nanocluster



Now that you've repaired the tooth, help them protect the restoration at home.

3M[™] Clinpro[™] Tooth Crème 0.21% Sodium Fluoride **Anti-Cavity Toothpaste**



Contains 0.21% NaF (950 ppm)

Non-prescription-strength formula 3M[™] Clinpro[™] Tooth Crème 0.21% Sodium Fluoride Anti-Cavity Toothpaste is the answer. In addition to a lower level of fluoride, it contains functionalised Tri-Calcium Phosphate (fTCP), available exclusively from 3M-so it releases calcium and phosphate at the tooth surface, which helps to maintain healthy teeth.

Protects against cavities and helps reverse white spots
Protects against cavities and helps reverse white spots



Proven through 5 years of clinical use

Cleans teeth with low abrasion

Replaces the patient's current toothpaste in the toothbrushing step, making compliance easy

Ordering Information



3M™ Filtek™ One Bulk Fill Restorative				
Refill	Item #	ltem #		
Shade	Syringe	Capsule		
A1	4866A1	4867A1		
A2	4866A2	4867A2		
A3	4866A3	4867A3		
B1	4866B1	4867B1		
C2	4866C2	4867C2		

SYRINGE SHADE REFILLS include: 1 - 4g syringe; Technique Guide; Instructions for Use.

CAPSULE SHADE REFILLS include: 20 - 0.2g capsules; Technique Guide; Instructions for Use.



3M[™] Clinpro[™] Tooth Crème 0.21% Sodium Fluoride Anti-Cavity Toothpaste

Item # Product Information

12216 Vanilla Mint 113g tube

www.3m.com.au/FiltekOne www.3m.co.nz/FiltekOne



3M Oral Care 3M Australia Pty Limited Building A, 1 Rivett Rd North Ryde NSW 2113 Ph: 1300 363 454 3M.com.au/dental

3M New Zealand Limited 94 Apollo Drive, Rosedale Auckland 0632 Ph: 0800 80 81 82 3M.co.nz/dental

Claims supported by 3M data on file. Please refer to instructions for use. 3M, ESPE, Clinpro and Filtek are trademarks of 3M. Information correct at time of printing. © 3M 2017. All rights reserved.