

# MICRO-TOUCH<sup>®</sup> DENTA-GLOVE<sup>®</sup> Green Neoprene

DENTA-GLOVE<sup>®</sup> Series

## Neoprene, powder-free dental exam gloves with ergonomic design for excellent flexibility and grip

This non-sterile, neoprene dental exam glove is designed to deliver the ultimate fit, comfort, tactile sensitivity, flexibility and grip. Featuring our ERGOFORM<sup>™</sup> ergonomic design technology, this glove helps reduce muscle effort that can cause stress on joints, ligaments and tendons.

### KEY FEATURES AND BENEFITS:

- Unique formulation for extra sensitivity with flexibility and comfort
- These neoprene dental gloves are strong and durable to avoid snagging and tearing through long procedures and donning
- Low to no odor
- Tested Safe for Oral Examination, ensuring that the glove is safe for patient contact

### Recommended For

- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Enhanced grip in wet environments
- Extra comfort during long, sustained procedures
- Protection from Type I latex allergy in HCW's or patients
- Oral examinations

### FEATURED TECHNOLOGIES



## MICRO-TOUCH<sup>®</sup> Everyday, Everywhere Protection

MICRO-TOUCH<sup>®</sup> premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

# MICRO-TOUCH<sup>®</sup> DENTA-GLOVE<sup>®</sup> Green Neoprene

D E N T A - G L O V E <sup>®</sup> S e r i e s

## PRODUCT INFORMATION

Material	Neoprene (Polychloroprene)	Not Made From Natural Rubber Latex	Yes
Color	Green	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.08 / 3.1
Finger Thickness (mm/mil)	0.11 / 4.3	Allergy Prevention	Latex (Type I)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	No
Sterile	No	Antistatic	Not Tested
Glove Length (mm/inches)	240 / 9.5	Product Segmentation	DENTA-GLOVE <sup>®</sup> Series

## ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Ansell Product Code	3130060XS	31300600S	31300600M	31300600L	3130060XL

## PACKAGING AND STORAGE

Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case *90 gloves per dispenser (Size XL) *900 gloves per case (Size XL)
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

## STANDARDS AND CERTIFICATIONS

Category III, EN ISO 374-1:2016, EN ISO 374-5:2016, EN 420:2003 + A1:2009, FDA21 CFR 177-2600, ISO 10993-10:2002, ISO 13485, ISO 9001

## Sensitisers

Zinc Dibutylthiocarbamate (ZDBC)

Only a very small number of users may be sensitive to this ingredient(s) and hence may develop irritant and/or allergic contact reactions.

Contact your Ansell representative for ordering or more information.

North America  
US Tel: 800 952 9916  
CA Tel: 1-844-494-7854  
Email:  
insidesalesus@ansell.com

Latin America & Caribbean  
Tel: +52(442) 296 20 50  
Email: cslac@ansell.com  
Brazil Tel: +55-11-3356-3100  
Brazil Email:  
luvas.medicas@ansell.com

Europe, Middle East & Africa  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu

Asia Pacific  
Tel: +603 8310 6688  
Email:  
apac.medical@ansell.com

Australia & New Zealand  
Tel: +61 3 9270 7270  
Email:  
protection@ap.ansell.com

Ansell, <sup>®</sup> and <sup>™</sup> are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending:  
[www.ansell.com/patentmarking](http://www.ansell.com/patentmarking).

© 2020 Ansell Limited. All Rights Reserved.