

For the daily and weekly cleaning and maintenance of dental aspiration and suction systems.

Liquid concentrate for simple preparation of solution for daily and weekly flushing of

- Aspiration systems
- Suction System
- Spittoons

Advantages:

- $\cdot \ \text{Three actions} \\$
 - decontaminates, sanitises and deodorises
- · Alkaline cleaner for removal of daily and weekly deposits
- · Simple and safe application
- · Foam free formula
- · Economical concentrate for preparation of:
 - ~ 2% dilution for daily cleaning and deodorising
 - ~ 6% dilution for weekly maintenance and decontaminates
- · Aldehyde free, Phenol free, Chlorine free







Dentalife Australia Pty Ltd Factory 9 /505 Maroondah Highway Ringwood, Victoria 3134 Australia Phone (03) 9879 1226 www.dentalife.com.au