Dysilk®

Black braided multifilament silk

Natural Non-absorbable

Suggested procedures

- Dental
- Oral maxillo facial
- Ophthalmic
- Subcuticular
- General / Intestinal
- Plastic

Features

- Excellent handling and tying properties
- Very low memory prevents knotting and tangling during suture removal
- Specially treated braid allows ease of handling
- Braid made of a core of twisted silk around which a cover is plaited
- Uniform thread diameter and surface smoothness for good knotting characteristics
- Doesn't become hard or difficult to handle
- Ties down with a high degree of knot security
- Available in sterile pre-cut lengths of various sizes
- Sutures made with 300 series stainless steel
- Up to five-year shelf life



Manufacturing Australian sutures for the world since 1974

Commonly used Dysilk[®] product codes and specifications.

Contact us for any other suture specification.

		Thread									
Needle*	USP Metric	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5		
3/8 Circle Reverse C	Cutting										
10.5mm 🔻 🧹	45cm	S601E10									
	45cm	S602	S502	S402							
12mm	76cm		S1502								
13mm 🔻 💟	45cm	S6040		S4040							
▼\ /	45cm				S307	S207	S07				
28mm	76cm				S1307	S1207	S107				
• /	45cm					S2031					
35mm	76cm				S13031	S12031					
40mm	76cm					S2013	S1013	S1113			
45mm	50cm					S2010	S010				
3/8 Circle Reverse P	Premium C	utting		1		1	1	1	1		
11mm	/ 45cm	DS602D									
13mm 🛛 💟	45cm		S5049								
	45cm	S604	S504	S404	S304	S204					
16mm	76cm		S1504	S1404	S1304	S1204					
	45cm	S605	S505	S405	S305						
19mm	76cm		S1505	S1405	S1305						
	45cm			S406	S306	S206					
24mm	76cm			S1406	S1306	S1206					
26mm	45cm				S3047	S2047					
	45cm				S309	S209	S09		S29		
30mm	76cm				S1309	S1209	S109				
1/2 Circle Reverse C	Cutting										
13mm 🔻 💟	45cm			S40#7							
16mm	45cm			S4015							
18mm	45cm			S4014	S3014						
25mm	76cm					S12034					
30mm	76cm				S13038	S12038	S1038				
	45cm		1				S012				
40mm	76cm						S1012	S1112			

* Needle images and sizes are approximations only.



Commonly used Dysilk[®] product codes and specifications.

Contact us for any other suture specification.

	Thread								
Needle*	USP Metric	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	
Straight Premium Cutting									
60mm	100cm				S3030	S2030	S030		
3/8 Circle Conventional Cuttin	g		1	1	1		1	1	
16mm	45cm 76cm		S5037 S15037	S14037					
3/8 Circle Fineline® (Conventional Premium) Cutting									
16mm	45cm	S6025	S5025	S4025	S3025				
20mm	45cm			S4026	S3026				
25mm	45cm		S5027	S4027					
3/8 Circle Round Bodied Taper									
10mm • 🗸	76cm		S5052	S4052					
13mm • 🗸	76cm	S6053	S5053						
16mm	76cm		S5054	S4054	S3054				
1/2 Circle Round Bodied Tape	-		1	1	1		1	1	
16mm	76cm			S4060	S3060	S2060			
20mm	76cm		S5065	S4065	S3065	S2065			
25mm	76cm			S4068	S3068	S2068	S068	S168	
30mm	76cm				S3069	S2069	S069	S169	
35mm	76cm					S2070	S070	S170	
1/2 Circle Round Bodied A-CU	TE [®] (Cu	tting Tip) Taper	1	I	I	I	1	
25mm	76cm					S2032			
Sterile Pre-Cut Lengths - Non	-Needle								
N/A	1x100cm			SA4040	SA3040	SA2040	SA040	SA140	
N/A	1x180cm			SA4072	SA3072	SA2072	SA072	SA172	
N/A	5x60cm			SE4024	SE3024	SE2024	SE024	SE124	
N/A	13x60cm			SM4024	SM3024	SM2024	SM024	SM124	
N/A	6x45cm			SF4018	SF3018	SF2018	SF018	SF118	
N/A	10x76cm			SJ4030	SJ3030	SJ2030	SJ030	SJ130	
N/A	12x76cm					SL2030			

* Needle images and sizes are approximations only.



Proudly Australian owned and operated.


Dynek is a leading, family owned manufacturer of sutures, committed to technical and surgical excellence. Established in 1974, Dynek has a proud history of innovation and quality as the only Australian manufacturer of sutures, with a singular focus of excellent patient outcomes.



dynek.com



1300 449 262 www.adamdental.com.au



- Always refer to the instructions for use when using Dynek sutures.
- Dysilk[®] is included on the ARTG (AU), Medsafe (NZ), and CE Mark (EU).
- Knot tying requires the standard surgical technique of flat and square ties with additional throws as indicated by surgical circumstance.
- A minimum order quantity may apply to all products.

QA: BE082/1 06/2022