

GenENDO Handfiles

A selection of the most popular handfile sizes to get to **the essentials of ENDO**.



GenENDO K-Files



Tried-and-tested range



GenENDO K-Reamers



Sustainable and reusable instruments⁽¹⁾



GenENDO H-Files



Manufacture 'Origine France Garantie' certified

Best-sellers

→ Pack of 6 instruments

→ Selection of 6 files:
> Ø 015 to 040

→ Available lengths:
> 21 mm
> 25 mm
> 31 mm

→ Diameters from 008 to 050⁽²⁾

→ Conicity of 2%

(1) See instructions for use



(2) You can find the list of references and instructions for use on our product page or on [genendo.com](https://www.genendo.com)



Micro-Mega SA
12 rue du Tunnel
25000 BESANCON / France
T +33 3 81 54 42 42
F +33 3 81 54 42 30
commercial@micro-mega.com

GenENDO K-Files, GenENDO K-Reamers, GenENDO H-Files
Advertising | Year of CE-marking: 2020 | Medical device class I according to regulation 2017/745
Medical device for dental care, meant for professional dental use only | See product labelling and instructions for use available on website: [genendo.com](https://www.genendo.com)

80000190-B 02.21