

Demi[™]Plus

LED Light Curing System

The best just got better.

Demi[™] Plus, powered by Demetron[™] Technology, is the next generation in premium curing lights. It combines sophisticated engineering with innovative design advancements that enhance performance and deliver consistent cures while allowing you to focus on what's truly important—your practice.



Demi Plus Technology

True 5-second cure



cure for composite shades delivers a true 5-second Periodic Level Shifting (PLS) Consistent output of

Improved cooling system



comfortable for the patient that allows the unit to and improved cooling system Fan-free for quiet operation touch for the clinician, and continue to feel cool to the



single charge. 10-second cures with a technology can yield 300 The lithium-ion battery

Advanced battery technology

Lightweight, ergonomic design



can reduce hand and arm Multiple holding options

& Rockwell **Hardness Depth of Cure**

	Depth of	Depth of Cure (mm)		Rockwell Hardness	Hardness	Ji
	5 Seconds	5 Seconds 10 Seconds	Top	5 Seconds Bottom	Top	10 Seconds op Bottom
Harmonize A2E	4.08	5.07	79.6	72.8	82.0	79.6
Harmonize C4E	3.50	4.12	78.0	63.5	82.0	75.4
Herculite Ultra A2E	4.81	5.74	78.8	72.7	80.2	79.5
Herculite Ultra C4E	3.13	3.67	75.5	49.6	79.1	70.8
Simile A2E	4.13	4.77	83.6	79.0	83.4	81.0
Herculite Ultra Flow A2	4.26	4.88	69.9	66.4	71.6	70.8
Herculite Ultra Flow A4	2.98	3.42	69.7	60.4	72.1	69.7

24-hour testing conducted by Kerr R&D. Data available upon request

Bulk Fill Products

Bottom /Top	Depth of Cure	Products
95.6%	5mm	SonicFill A3
99.6%	5mm	Filtek One Bulk Fill A3
97.9%	4mm	Tetric Evocream Bulk Fill IVA
98.7%	4mm	SDR Flow+
99.7%	4mm	GrandioSo X-tra A3
98.2%	4mm	Estelite Bulk Fill Flow A3

Pulsing Power of PLS

80.72	80.58	76.22	Rockwell 10 seconds
TPHLV	TPHST=CeramX	Tetric Evoceram	
87.28	79.53	79.57	Rockwell 10 seconds
Filtek Supreme Ultra Body	Herculite Ultra	Harmonize	

Demi with PLS 5-Second Exposure

		C	UTPU	mW/cr	m2
		1100	1200	1300	1400
	0				
	_	F		\Rightarrow	
TIME	2	F		\Rightarrow	
TIME IN SECONDS	ω	F		\Rightarrow	
01	4	F		\Rightarrow	
	ιπ	L		\Rightarrow	
		L			

Advanced Features

by Demetron Technology. A Kerr curing light system powered



Longer Warranty

Our three-year warranty ensures you get the most value out of your premium investment.

Demetron Light Guide



921676 | 4 mm | Curved Turbo

For preventing overlapping on the gingiva in Class V.



952213 | 8 mm | Curved Turbo+

Increases efficiency of light output and reduces overall curing time.



21020 | 8 mm | Curved Turbo

Increases efficiency of light output and reduces overall curing time.



921746 | 8 mm | Extended Turbo+

Extended length improves access to the posterior.



921551 | 11 mm | Curved Turbo

Increases efficiency of light output and reduces overall curing time.

Light guides are compatible with most Demetron Light Curing Systems, including Demi, LED II, and OptiLux 501.

Ordering Information

910860-1 Demi Plus LED Light Curing System (120 V)

Contains: Handpiece, 8 mm Turbo Extended Light Guide, Lithium Ion Battery Pack, Battery Charger, Power Supply,

Protective Light Shield.

921918-1 Demi Plus Battery Pack921919-1 Demi Plus Battery Charger

921920-1 Demi Plus Handpiece Replacement952269 Demi and Demi Plus Power Supply

910726 L.E.D. Radiometer



