

Maxcem Elite[™] and Maxcem Elite[™] Chroma Universal Resin Cement

One cement to meet all your needs.

Universal Strength and Simplicity.

Simplify your inventory requirements and keep using your preferred cementation technique, with Maxcem Elite[™] and Maxcem Elite[™] Chroma Universal Resin Cement. The only cement you need to meet all of your indirect restoration needs, delivering universal versatility, higher bond strengths, predictable handling, and exceptional ease of use.



Features & Benefits



Best of Both Worlds

Enjoy the simplicity of a self-adhesive cement and the confidence of an adhesive resin cement when used with a universal bonding agent. Maxcem Elite™ and Maxcem Elite™ Chroma Universal Resin Cement provide excellent cementation for indirect restorations in self-cure, light-cure, or dual-cure modes.



Universal Compatibility

Simplify with one cement for your indirect restoration needs. Maxcem Elite[®] and Maxcem Elite[®] Chroma Universal Resin Cement are compatible with all dental substrates and indirect restorations including anterior or posterior ceramics, PFMs and metal restorations, and CAD/CAM materials.



Boost Your Bond Strength

A self-etch/self-adhesive cement, Maxcem Elite[™] Universal Resin Cement offers one of the most secure holds compared to other leading self-adhesive cements on the market in both dry and moist environments.



Smart and Simple

One-peel" removal simplifies your workflow even more. In addition, Maxcem Elite" Chroma Universal Resin Cement dispenses pink before fading at the gel state — indicating the optimal time to clean up excess cement. With no hand mixing, no refrigeration required and no guesswork, you'll enjoy easier and consistent clean up every time.

ONE CEMENT TO MEET ALL YOUR NEEDS

Use Maxcem Elite™ Universal Resin Cement as a self-adhesive or pair with any universal bonding agent, such as OptiBond™ Universal on preparations that need more retention.

Long-Lasting Esthetics

Maxcem Elite[™] and Maxcem Elite[™] Chroma Universal Resin Cement offer long-lasting color stability and exceptional esthetics. With a proprietary and proven redox initiator system also used in NX3 Nexus[™], any color shifts over time are virtually imperceptible.




The Ultimate Power Combo

For restorations that require even more strength, Maxcem Elite[™] Universal Resin Cement paired with OptiBond[™] Universal delivers excellent bond compatibility, with a wide range of dental substrates for all indirect procedures. This power combo ensures high strength for more confident results.


Because there's no need to light-cure when combining with OptiBond™ Universal, you'll also be able to simplify your protocol while boosting your bond strength.




Maxcem Elite™ vs. Other Self-Adhesive Cements



Bond Compatibility of Maxcem Elite™ with Various Adhesives



Maxcem Elite™/OptiBond™ Universal vs. Other Cement/Adhesive Systems



Data on file. Maxcem Elite" Universal Resin Cement is a trademark of Kerr Corporation. All other trademarks are properties of their respective owners.

Case Studies

Images courtesy of Dr. Matthew Miller

Failing Direct Composite Veneers

Teeth 7-11 (4 years old)



Teeth 6-11 Cemented with Maxcem Elite™ Chromal (Total Etch Technique with OptiBond™ Universal)



Before



After: Notice the cement color change in the old veneers



Veneers Removed



Porcelain Veneers

Veneers Immediately After Cementation with Maxcem Elite™ Chroma (Total Etch Technique with OptiBond[™] Universal)



Veneers at 1 Year Follow-Up



Before



Crowns 6 Month Follow-Up (Maxcem Elite Chroma™ Self-Adhesive Technique)



Crowns 1.5 Year Follow-Up (Maxcem Elite Chroma™ Self-Adhesive Technique)



Before



Day of Placement



Ordering Information

Maxcem Elite™ Universal Resin Cement 34417 Intro Kit 34418 Bulk Pack 33872 Clear Refill (2 dual syringes) White Refill (2 dual syringes) 33873 33874 White Opaque Refill (2 dual syringes) 33875 Yellow Refill (2 dual syringes) 32862 Automix Regular Tips 4:1 (pkg of 50) 33879 Automix Wide Tips 4:1 (pkg of 10) 33880 Curved Intraoral Tips (pkg of 10) Curved Root-Canal Tips (pkg of 10) 33881

Maxcem Elite™ Chroma Universal Resin Cement

36299 Bulk Pack 36300 Clear Refill (2 dual syringes) 36301 White Refill (2 dual syringes) 36302 Yellow Refill (2 dual syringes) 36303 Pink Automix Regular Tips (pkg of 50)

OptiBond™ Universal

36517 Bottle Kit

36518 Unidose Refills (pkg of 100)

36519 Bottle Refill (5ml)

26909 Gel Etchant Dispensing Tips (pkg of 30)

31297 Gel Etchant Refill Pack

Visit Go.Kerrdental.com/Maxcem for more information.

Data on file, Maxcem Elite Universal Resin Cement is a trademark of Kerr Corporation. All other trademarks are properties of their respective owners.

