

Truly Universal Adhesive Resin Cement



Your practice is our inspiration...

NX3 Perfection in permanent cements



NX3 Defining a new generation in permanent cements

an innovative chemistry for unmatched esthetics applications, including veneers. adhesion and great versatility. Delivery systems time. NX3 cement can be used for all indirect tions and a light-cure cement for unlimited work include an automix syringe for dual-cure indica-NX3 is a permanent resin cement system with

Research

matrix, NX3 is the first truly color-stable adhebility for both total-etch and self-etch adhesives sive resin cement with enhanced bond compati initiate polymerization and curing, resulting in a benzoyl peroxide and tertiary amine pair to free redox initiator system and optimized resin with self-etch adhesives. Now with an aminediscoloration over time due to incompatibility Most self-cure or dual-cure resin cements use

Features

CAM blocks, ceramic, porcelain, Superior bonding. Bonds to all substrates. resin and metal. Excellent adhesion to dentin, enamel, CAD/

both dual-cure and light-cure cements. True color stability. Long-term esthetics for compatibility. No dual-cure activator required Self-etch or total-etch. Bonding protocol

Shear Bond Strength (MPa) to Dentin – Self-Cured SELF-CURE MODE NO ACTIVATOR

Ordering

33642B NX3 Nexus[®] Third Generation Intro Kit w/risk-free trial offer

Technique Guide; Directions for Use Gel (3 g) of each shade: clear, white and each shade: clear, white and yellow; 1 Syringe Try-In clear, white and yellow; 1 Syringe Light-Cure (1.8 g) of Contains 1 Syringe Dual-Cure (5 g) of each shade: yellow; 1 Bottle Silane Primer (5 ml); 30 Automix Tips;

33682 NX3 Light-Cure Kit

10 Gel Etchant Syringe Tips; Directions for Use 10 Applicator Tips; 1 Syringe Gel Etchant (3 g): Unidose Pouches; 1 Bottle Silane Primer (5 ml) each shade: clear and white; 10 OptiBond Solo Plus shade: clear and white; 1 Syringe Try-In Gel (3 g) of Contains 1 Syringe Light-Cure (1.8 g) of each

NX3 Automix Dual-Cure Syringe Refills

33644	33643
White	Clear

i Yellow	VVIIIUG
W	'n

33646	33645
Bleach	Yellow

33647	
White	
Opaque	

NX3 Light-Cure Syringe Refills

33649	33648
White	Clear

- 33650 Yellow
- 33651 Bleach
- 33652 White Opaque


NX3 Try-In Gel Syringe Refills

- 33657 33656 White Clear
- 33658 Yellow
- White Opaque

33659

- 33660 Bleach

- 33655 50 Automix Tips



Trademarks are property of their respective owners

Nexus® Third Generation

Variolink® II

Calibra®

RelyX[™] ARC

Bond strength study conducted by Dr. James Dunn of Loma Linda University

NX3

37°C water 28 weeks in After

Initial

Color Stability of Dual-Cure Resin Cements – Clear Shade



©2008 Kerr Corporation Lit. No. 34000

Kerr Corporation 1717 West Collins Avenue, Orange, CA 92867 800.KERR.123