

OptiBond[®] Universal

Single-Component Universal Adhesive



ONE STEP. ONE BRAND. ANY PROCEDURE.


OptiBond[™] Universal is here. Our one-bottle, one-coat bonding agent combines the reliability of patented OptiBond GPDM technology with the ease-of-use of a Universal solution. Exceptional bond strength to all surfaces and substrates, in any etching technique.

Efficacy and flexibility coming together at last.

THE GOLD STANDARD IN ADHESION.




Kerr's proven **OptiBond™ GPDM Monomer** ensures more effective and reliable adhesion than other adhesive monomers to both the tooth structure and the resin material. It contains both acidic phosphate groups to effectively anchor the adhesive to the hydrophilic surface of the tooth and more methacrylate functional groups than other adhesive monomers (i.e. MDP) for a more effective interface with the hydrophobic resin. This creates a better compatibility with both substrates and an overall higher bond strength.



Shear Bond Strengths of


Shear Bond Strengths of Universal Adhesive to Enamel



Kerr OptiBond[™] Universal yields outstanding bond strength on both dentin and enamel, creating maximum flexibility with outstanding results under different etching techniques.

SIMPLICITY WITH NO COMPROMISE.


OptiBond[™] Universal technology delivers excellent penetration into dentinal tubules, ensuring extraordinary bond strength and protection against microleakage and post-op sensitivity, with an application time of 30 seconds. Its unique nanoetching capability delivers more effective enamel etching than other single-component adhesives, creating a deeper etched surface for higher micro-mechanical retention, and therefore higher bond strength.




Uncut enamel 50,000x - OptiBond Universa




Uncut enamel 50,000x - leading competitor




Dentine tubules 2,000x – OptiBond Universal




Dentine tubules 2,000x - leading competitor





Shear Bond Strengths, MPa

Using Kerr OptiBond[™] Universal/NX3[™] Cement together provides enhanced bond compatibility and therefore outstanding indirect bond strength to dentin.



The combination of Kerr OptiBond[™] Universal/NX3[™] Cement yields outstanding indirect bond strength to various dental substrates, without a separate ceramic or metal primer.

OptiBond Universal Single-Component Universal Adhesive

ALL ETCHING TECHNIQUES IN ONE BOTTLE.

Reconstructing the tooth structure involves different substrates with different proven technology for different substrates and techniques, from Total Etch to Self-

Keri OptiBond Universa

See How Harmonize[™] and OptiBond[™] Universal Deliver Esthetic Restorative Solutions



Case 1 - Dr. Scott Coleman completes a class IV in tooth #8 to restore the incisal, distal corner. A selective etch technique was used with OptiBond™ Universal and cured with the Demi™ Ultra curing light. A layer of Harmonize[™] dentin shade A3 was applied lingually and cured, followed by Harmonize[™] enamel shade A2 with additional curing.



Case 2 - Dr. Scott Coleman restores multiple teeth with Harmonize[™], using a total etch technique with OptiBond[™] Universal followed by Harmonize[™] dentin shade A1 followed by enamel B1 and a final layer of Harmonize[™] Translucent Blue.

"In a survey of 100 dentists, 100% reported no post-operative sensitivity upon follow-up.*"



Case 3 - Dr. Drew Ballard starts with a dark dentin shade followed by an enamel shade using Harmonize[™] and OptiBond Universal.





ORDERING INFORMATION

36517	OptiBond [™] Universal Bottle Kit
36518	OptiBond [™] Universal Unidose (100 Pack)
36519	OptiBond [™] Universal Bottle Refill
26909	Gel Etchant Dispensing Tips (pkg of 30)
31297	Gel Etchant Refill Pack
36440	Harmonize [™] A1E Unidose
36441	Harmonize [™] A2E Unidose
36442	Harmonize [™] A3E Unidose
36432	Harmonize [™] A1D Unidose
36433	Harmonize™ A2D Unidose
36434	Harmonize™ A3D Unidose
36414	Harmonize [™] Translucent Clear Syringe

Visit KerrDental.com for additional shades and delivery system.

Kavo Kerr

PREFERRED LOYALTY PROGRAM BUY NOW. EARN REWARD POINTS. SIGN UP AT KAVOKERRPREFERRED.COM

kerrdental.com/optibonduniversal *Internal Data on File MKT-17-0478

