



## Soft-Core™

Carrier-based obturation — finally a better alternative.



## 5 key reasons to obturate with Soft-Core.™




### INJECTION MOLDED

Unlike other brands, Soft-Core carriers are injection molded for consistency.



### BENDABLE HANDLE

The bendable handle allows for easier posterior access. The removable handle and hollow core creates a pilot hole for a post drill.




### THINNER CORE

The core is thinner than other leading brands, allowing greater flexibility to navigate curved canals.



### DENSER FILLS

A thinner core also means more gutta percha for denser fills, up to a .10 taper.



### UNIVERSALLY COMPATIBLE

Soft-Core's one universal size taper fits most common root canal preparations.

# 1

# 2

# 3

# 4

# 5

Soft-Core™ from Kerr Endodontics, includes carriers, verifiers and an oven for a complete carrier-based obturation system. Available either in the classic Soft-Core line, or branded to match TF Adaptive, Twisted File, and K3XF file systems.



The Soft-Core heater holds up to 4 obturators at a time.


Bendable, easy-to-remove, twist-off handle

Removable handle creates a pilot hole for a post drill

Thinner core for greater flexibility to navigate curved canals

More gutta percha fills up to .10 taper

One size taper fits most common preps



### A better solution for carrier-based obturation

The Soft-Core system was designed with the clinician in mind. With tip sizes from #20 through #60, matching verifiers and the Soft-Core oven, 3-dimensional warm vertical obturation is at your fingertips.

### Simplicity


With a universal one-size taper, Soft-Core obturators have been designed to accommodate most common root canal preparations. Available in the convenient 6-pack, or economy size 36 pack, as well as branded TFA, TF, and K3XF to match your NIT file system of choice.

### Efficiency

The twist-off handle is easy to remove and leaves behind a pilot hole for when a post is to be used. Because the handle is flexible, it can be pre-bent to allow for greater posterior access.

### Predictability

Unlike other leading carrier brands, Soft-Core carriers are injection molded in order to obtain consistent, accurate tolerances required in Root Canal Treatment.



**Soft-Core**  
Carrier-Based Obturators

## SOFT-CORE OBTURATORS & VERIFIERS

### SOFT-CORE OBTURATORS

6  
PACK

SIZE	PRODUCT	ITEM
#20	972-3020	Obturator (6pk)
#25	972-3025	Obturator (6pk)
#30	972-3030	Obturator (6pk)
#35	972-3035	Obturator (6pk)
#40	972-3040	Obturator (6pk)
#45	972-3045	Obturator (6pk)
#50	972-3050	Obturator (6pk)
#55	972-3055	Obturator (6pk)
#60	972-3060	Obturator (6pk)

### SOFT-CORE OBTURATORS

36  
PACK

SIZE	PRODUCT	ITEM
#20	972-3120	Obturator Economy Pk (36pk)
#25	972-3125	Obturator Economy Pk (36pk)
#30	972-3130	Obturator Economy Pk (36pk)
#35	972-3135	Obturator Economy Pk (36pk)
#40	972-3140	Obturator Economy Pk (36pk)
#45	972-3145	Obturator Economy Pk (36pk)
#50	972-3150	Obturator Economy Pk (36pk)
#55	972-3155	Obturator Economy Pk (36pk)
#60	972-3160	Obturator Economy Pk (36pk)

### SOFT-CORE VERIFIERS

6  
PACK

SIZE	PRODUCT	ITEM
#20	972-3220	Size Verifier #20 (6pk)
#25	972-3225	Size Verifier #25 (6pk)
#30	972-3230	Size Verifier #30 (6pk)
#35	972-3235	Size Verifier #35 (6pk)
#40	972-3240	Size Verifier #40 (6pk)
#45	972-3245	Size Verifier #45 (6pk)
#50	972-3250	Size Verifier #50 (6pk)
#55	972-3255	Size Verifier #55 (6pk)
#60	972-3260	Size Verifier #60 (6pk)

## TF™ ADAPTIVE OBTURATORS & VERIFIERS

### TF™ ADAPTIVE OBTURATORS

6  
PACK

SIZE	PRODUCT
SM2	815-1550
SM3	815-1551
ML1	815-1554
ML2	815-1555
ML3	815-1556

### TF™ ADAPTIVE OBTURATORS

36  
PACK

SIZE	PRODUCT
SM2	815-1560
SM3	815-1561
ML1	815-1562
ML2	815-1563
ML3	815-1564

### TF™ ADAPTIVE VERIFIERS

6  
PACK

SIZE	PRODUCT
SM2	815-1552
SM3	815-1553
ML1	815-1557
ML2	815-1558
ML3	815-1559

## K3/K3XF OBTURATORS & VERIFIERS

### K3/K3XF OBTURATORS

6  
PACK

SIZE	PRODUCT
#60	971-3060
#55	971-3055
#50	971-3050
#45	971-3045
#40	971-3040
#35	971-3035
#30	971-3030
#25	971-3025
#20	971-3020

### K3/K3XF VERIFIERS

6  
PACK

SIZE	PRODUCT
#60	971-3160
#55	971-3155
#50	971-3150
#45	971-3145
#40	971-3140
#35	971-3135
#30	971-3130
#25	971-3125
#20	971-3120

## TF™ OBTURATORS & VERIFIERS

### TF™ OBTURATORS

6  
PACK

SIZE	PRODUCT
#25	822-3025
#30	822-3030
#35	822-3035
#40	822-3040
#45	822-3045
#50	822-3050

### TF™ OBTURATORS

36  
PACK

SIZE	PRODUCT
#25	822-3625
#30	822-3630
#35	822-3635
#40	822-3640
#45	822-3645
#50	822-3650

### TF™ VERIFIERS

6  
PACK

SIZE	PRODUCT
#25	822-3125
#30	822-3130
#35	822-3135
#40	822-3140
#45	822-3145
#50	822-3150

## ALSO RECOMMENDED:

ITEM	PRODUCT
Soft-Core Oven	972-3002
LA Axxess Stainless Steel Bur (6pk)	815-1401
Absorbent Points (assorted cell pack)	16221
Sealapex Xpress (2 - 10.5g syringes)	33639
PeerlessPost Intro-Kit	977-1040

## SOFT-CORE INTRO-KIT PRODUCT #: 972-3000

1 SOFT-CORE OVEN



+

VERIFIERS (6PK)  
2 PACKS EACH: #20, #25, #30, #35



+

OBTURATORS (6PK)  
3 PACKS EACH: #20, #25, #30, #35



+

#1 LA AXCESS STAINLESS  
STEEL BURS (6PK) 1 PACK



For more information visit [kerrdental.com](http://kerrdental.com) or call us at 800-KERR-123