

Invisible below  
the temporary  
restoration

Visible during  
removal of  
excess  
material

# Bifix® Temp

DUAL-CURING, COMPOSITE-BASED TEMPORARY  
LUTING MATERIAL



# Bifix® Temp

## EASE OF USE AND OUTSTANDING AESTHETICS COMBINED IN A SINGLE PRODUCT

With Bifix Temp the perfect aesthetic result is already achieved when luting temporary restorations. The translucent, tooth-like Universal shade is ideal for highly aesthetic or very thin and thus transparent temporary restorations, does not shine through them and therefore enhances the natural appearance of the restoration. What's more, it is also extremely easy to use.

### Indication


Non-permanent luting of temporary restorations

### Advantages


- Perfect fit of restoration – film thickness less than 5 µm does not affect the occlusion
- Flexible workflow – adjustable time settings for polymerisation due to dual-curing
- Ease of use – in daily routine proven QuickMix syringe yields perfect components self-mixture
- Excess material is at a stretch removable after a brief period of light-curing
- Cost saving – minimal material loss by shortest mixing tip



### Influence of the luting material on Structur 3



Opaque temporary luting materials can have a major effect on the aesthetics of temporary restorations being fitted. In the anterior region in particular this is a highly undesirable factor.



Source: Results of survey among testing dentists

### Presentation

- REF 1055 QuickMix syringe 5 ml universal, accessories
- REF 2261 Dispenser – QuickMix
- REF 2188 Mixing tips type 9, 50 pcs.

VOCO GmbH  
 Anton-Flettner-Straße 1-3  
 27472 Cuxhaven  
 Germany  
 Freecall: 00 800 44 444 555  
 Fax: +49 (0) 4721-719-140  
 info@voco.com  
 www.voco.dental

Available from:

