



## Prophylaxis Set

## Prophylaxis Set

PROTECTION. CARE. CONSERVATION.

Keeping teeth healthy for a long time is the common goal towards which you and your patients work together with us. For this to succeed, VOCO products make an important contribution: With instantly effective **fluoride varnishes** and **desensitising agents** suitable for a great variety of indications. With **fissure sealants** offering high filler content. With **cleaning pastes** in various grades. Our **prophylaxis set** offers you the opportunity to try a number of different products from the oral care section. If you like these products, then we recommend placing a repeat order straight away. The REF numbers of the respective products can be found in our catalogues. Simply enter these in the order form and either fax the form directly to us or hand it to your VOCO dental consultant for processing.

## CleanJoy®

Fluoride-containing tooth cleaning and polishing paste in three cleaning grades

- Stable, homogenous consistency
- Free of parabens
- No splashing (2,000 - 3,000 r / min)
- Traffic light coding: simple to use without confusion
- RDA values: coarse (red) = 195, medium (yellow) = 127, fine (green) = 16\*
- Can be used with any conventional instrument
- Contains xylitol and fluoride (700 ppm)
- Pleasant light mint flavour



## VOCO Profluorid® Varnish

Fluoride-containing dental desensitising varnish (5 % NaF)

- Easy to apply on wet surfaces
- Esthetic, tooth-shaded varnish
- Fast desensitisation and fluoride release (5 % NaF  $\cong$  22,600 ppm fluoride)
- Available in six pleasant flavours\*: melon, mint, cherry, caramel, bubble gum and cola lime
- Economical application in thin layers, high yield



## Remin Pro®

### Protective dental care with fluoride and hydroxy apatite

#### Remin Pro®

- Contains fluoride (1,450 ppm), hydroxy apatite and xylitol
- Provides a balanced oral flora and thus protects the teeth against harmful acid attacks
- Strengthens the teeth after whitening and professional cleaning
- Cariostatic properties from the contained xylitol
- Available in three aromatic flavours\* (melon, mint and strawberry)
- Can be used daily and especially well-suited for use at home



## Grandio® Seal

### Light-curing nano-hybrid fissure sealant

- Use of nano-fillers for optimal flow behaviour
- With 70 % w/w highest filler content in its class
- Outstanding physical properties in abrasion, strength and shrinkage values
- Excellent handling
- Optimal wetting properties
- Perfect marginal adaptation



\*set does not include all product presentations

