

SAFETY DATA SHEET

Section 1: IDENTIFICATION

V-WIPES

Recommended Use: Disinfectant & Cleaner.

Product Code: 210561 (Carton of 4 Canisters), 210579 (Carton of 12 Flat Packs)

Whiteley Medical

A division of Whiteley Corporation Pty Ltd (A.C.N. 000 906 678)
Postal Address: P. O. Box 1076 North Sydney NSW 2059
Telephone Number: (02) 9929 9155 Facsimile: (02) 9929 9077

Web: www.whiteley.com.au

Emergency Telephone Number: Poisons Information Centre (National) 131126

Distributed in New Zealand by: Whiteley Corporation NZ Limited

Address: PO Box 22-519 Otahuhu, Auckland, New Zealand 1640
Telephone Number: 092737313, 0800 257 352
Emergency Telephone Number: New Zealand Poisons Centre 0800 POISON or 0800 764 766

Section 2: HAZARDS

GHS Classification

Not a dangerous substance or mixture according to GHS or HSNO criteria.

Signal Word None

Hazard Statements

None

Precautionary Statements

None

Response Statements

None

Storage Statements

None

Disposal Statements

None

Issue No. 1 GHS Page 1 of 5

SAFETY DATA SHEET V-WIPES

Section 3: COMPOSITION INFORMATION Ingredient CAS No **Proportion** 64-17-5 <10% Ethanol To 100% Ingredients deemed not to be hazardous Not applicable Section 4: FIRST AID Eye (Contact) Hold eyelids apart and flush the eye continuously with running water. Skin (Contact) Remove contaminated clothing and flush skin and hair with running water. Inhalation(Breathing) Remove to fresh air. Ingestion (Swallowing) DO NOT induce vomiting. Give water to drink if conscious. For advice, contact a Poisons Information Centre (Phone 131126) or a doctor. **Advice to Doctor** Treat symptomatically for neutral detergent. **First Aid Facilities** Ensure an eye wash is available and ready for use. **Additional Information** No aggravated medical conditions are known to be caused by exposure to this product. Section 5: FIREFIGHTING MEASURE Suitable Extinguishing Use water fog (or if unavailable fine water spray), dry chemical, carbon dioxide Media or alcohol stable foam. **Hazards From** Carbon dioxide and carbon monoxide may be produced in the case of fire or **Combustion Products** during thermal decomposition. **Precautions For Fire** Keep containers cool by spraying with water. Fighters and Special **Protective Equipment Additional Information** Hazchem Code – Not applicable.

| Section 6: ACCIDENTAL RELEASE MEASURES | | |
|--|---|--|
| Emergency Procedure | SAA/SNZ HB76: Dangerous Goods – Initial Emergency Response Guide – Not applicable. | |
| Spills / Clean up | Disposal to infectious waste if contaminated recommended. If not contaminated then disposal to normal waste acceptable. Refer to State/Territory Land Waste Management Authorities if applicable. Containers are recyclable and can be disposed of by a licensed waste contractor. Containers can be disposed of to general waste or rinsed thoroughly and recycled. | |
| Section 7: HANDLING AND STORAGE | | |
| Precautions for Safe Handling | Contact Whiteley Corporation sales representative for advice when using this product for any application other than that outlined on the label or technical bulletin. Any non-intended or non-authorised use of this product may result in personal injury or damage to equipment. Store product in original container. Wash thoroughly after handling product. | |
| Conditions for Safe Storage | Store in a cool, dry, well ventilated area. Keep container tightly sealed. Store below 30°C. | |

Issue No. 1 GHS Page 2 of 5 SAFETY DATA SHEET V-WIPES

Section 8: EXPOSURE CONTROL/PERSONAL PROTECTION

National Exposure Standards – Source: National Exposure Standards for Atmospheric Contaminants in the Occupational Environment [NOHSC:1003] & New Zealand Workplace Exposure Standards 2002.

IngredientCAS NoES-TWAES-STELEthanol64-17-51000ppm-

Biological Limit Values No data available.

Engineering Controls Not required.

Personal Protective Equipment

Eye/face protection - Safety glasses / face shield / chemical resistant

goggles should be worn to prevent eye contact.

Skin protection – Use Nitrile gloves or similar to prevent skin contact. Respiratory protection – Respirator is not usually necessary but if

required use a half face filter respirator suitable for organic vapours.

| Section 9: PHYSICAL AND CHEMICAL PROPERTIES | | | |
|---|-------------------------------|--------------------------|-------------------|
| Appearance | White and pink apertured wipe | Odour | Citrus |
| Odour Threshold | No data available | рН | 6 – 8 |
| Freezing Point | No data available | Initial Boiling Point | No data available |
| Boiling Range | No data available | Flash point | No data available |
| Evaporation Rate | No data available | Flammability | No data available |
| Upper flammability limit | No data available | Lower flammability | No data available |
| Vapour Pressure | No data available | Vapour Density | No data available |
| Relative Density | No data available | Solubility | No data available |
| Partition Coefficient: n-octanol/water | No data available | Autoignition temperature | No data available |
| Decomposition temperature | No data available | Viscosity | No data available |

| | Section 10: STABILITY AND REACTIVITY |
|----------------------------------|--|
| Chemical Stability | Product is stable and will not undergo any hazardous reactions under normal conditions of use and storage. |
| Conditions to avoid | Avoid high temperatures (store below 30°C). |
| Incompatible materials | None known. Do not mix with other chemicals. |
| Hazardous decomposition products | None known. |
| Hazardous reactions | None known. |

Issue No. 1 GHS Page 3 of 5

SAFETY DATA SHEET V-WIPES

Section 11: TOXICOLOGICAL INFORMATION No data available. **Acute toxicity** Skin corrosion/irritation May cause irritation. Serious eye damage/ May cause irritation and reddening. irritation Respiratory or skin No data available. sensitisation **Germ Cell Mutagenicity** No data available. Carcinogenicity No components of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC. No data available. Reproductive toxicity **Specific Target Organ** No data available. **Toxicity - Single Exposure Specific Target Organ** No data available. **Toxicity – Repeated Exposure Aspiration Hazard** No data available. Section 12: ECOLOGICAL INFORMATION **Ecotoxicity** No data available. Persistence and No data available. degradability Bioaccumulative potential No data available. Mobility in soil No data available. Other adverse effects No data available. Section 13: DISPOSAL CONSIDERATIONS **Disposal method** Disposal to infectious waste if contaminated recommended. If not contaminated then disposal to normal waste acceptable. Refer to State/Territory Land Waste Management Authorities if applicable. Containers are recyclable and can be disposed of by a licensed waste contractor. Containers can be disposed of to general waste or rinsed thoroughly and recycled.

Issue No. 1 GHS Page 4 of 5

None.

Special precautions

SAFETY DATA SHEET V-WIPES

Section 14: TRANSPORT INFORMATION

Not classified as Dangerous Goods by the criteria of the Australian Dangerous Goods Code (ADG Code).

UN Number

UN Proper Shipping Name
Class and subsidiary risk
Packing Group

Special precautions for
Hazchem Code

Not applicable.
Not applicable.
Not applicable.
Not applicable.
Not applicable.

Section 15: REGULATORY INFORMATION

Poisons Schedule (SUSDP): Not applicable.

HSNO Approval Code: Not applicable.

All ingredients are listed in the Australia Inventory of Chemical Substances (AICS).

This document has been produced in accordance with the requirements of the Globally Harmonised

System of Classification and Labelling.

Section 16: OTHER INFORMATION

Date of 12th September 2016 **preparation:**

© This document remains property of Whiteley Corporation Pty Ltd.

Alterations are not permitted without prior written authorisation from Whiteley

Corporation Pty Ltd.

Issue No. 1 GHS Page 5 of 5