



**Water Hardness  
Test Kit**  
Product Code:  
303004



**HACH Chloride  
Test Strips (40 strips)**  
Product Code: 303022  
\*Suitable for testing Chloride levels  
between 30-600ppm.



**Conductivity  
Meter**  
Product Code:  
303025

## Endoscopy and CSSD Water Quality Requirements for Compliance AS 5369:2023

As a safety requirement, healthcare organisations are required to regularly test and monitor their water, in accordance with AS 5369:2023.

This includes monthly water hardness, conductivity, chloride testing, and annual endotoxin testing.

Whiteley offers water testing kits and endotoxin compliance testing, to help your endoscopy and sterilising departments achieve compliance.

Whiteley also offers practical compliance support. We can train your staff and provide ongoing assistance, to make sure you follow industry best practice.

### Water testing for Endoscopy & CSS

#### Annual endotoxin testing for compliance with AS 5369:2023

<b>Final Rinse Water Endotoxin Testing*</b>	Endotoxin: AERs 30 EU/mL, Washer Disinfectors 0.25 EU/mL	Annual Endotoxin Testing Fee Product Code: 303024
---	---	---

For testing items such as conductivity, chloride and hardness, these are to be done in accordance with the AER manufacturers recommendations.

### Water testing for Endoscopy & CSSD

As above with the following recommendations:

AS 5369:2023 WATER QUALITY USED FOR REPROCESSING OF REUSABLE MEDICAL DEVICES		
Substance	Water sample	Maximum concentration levels
<b>Water Hardness</b>	Sample Location	150mg/L
<b>Chloride**</b>	Supply Water	120mg/L
<b>Conductivity</b>	Final Rinse Water	<30µs/cm

\*Microbiological Compliance testing is also required monthly and should be conducted with your normal microbiology test provider.

\*\*Contact your Whiteley™ Sales Representative for further information when testing below 20mg/L.