



P R O M O T I O N

---

OCT - DEC 2023

# PRIME LINE



AeroAcoustic blade



Glass rod fiber optics



Reverse flow protection



Automatic pressure regulation<sup>1</sup>



Refined ceramic bearings



Scratch resistance  
Iconic Prime body

25  
W

Powerful & Precise  
Power Head - 25W

55  
dB

Incredibly Quiet  
Only 55 dB noise



Available for  
KaVo\*, NSK\*, W&H\*  
and Sirona\* Couplings



## Coupling Connection

Kavo  
NSK  
W&H  
Sirona

## Power Head 25 Watt

HP21KL  
HP21NL  
HP21WL  
HP21SL

## Small Head 18 Watt

HP22KL  
HP22NL  
HP22WL  
- N/A -

**\$1,683**

**BUY 2  
GET 1  
FREE**

**2 YEARS  
UNCONDITIONAL  
WARRANTY\***

\*Excluding physical damage on Prime Line handpieces

## COUPLINGS



- ✓ MK-dent & Kavo with LED & water
- QC6016KW

**\$495**



- ✓ NSK connection with LED
- QC6016NT

**\$385**



- ✓ W&H connection with fibre optic
- QC5016W

**\$385**



- ✓ Sirona connection with LED & water
- QC6016SW

**\$495**

- ✓ Non-optic
- QC4014K

**\$296**

- ✓ Non-optic
- QC4014NT

**\$296**

- ✓ Non-optic
- QC4014W

**\$296**

- ✓ Non-optic

**-NA-**

1) Only available for Kavo coupling 2) FairCare warranty applicable on Prime Line high speeds sold after 1st Oct 2018

Buy 3 **GET 15% OFF**  
Buy 4 **GET 20% OFF**  
Mix & Match

**NEW RANGE!**



2 Years FairCare Warranty  
Unconditional warranty  
except physical damage



SAFE-TECH  
PUSH  
BUTTON

⚙️ 1:5 Red ●

⌚ Max. 200,000 rpm

○ Internal water

🔧 For FG-burs  
(Ø 1.6 mm)

• LP15L (F/Optic)

**\$1,826**

**BUY 3**  
**\$1,552.10<sup>EA</sup>**

**BUY 4**  
**\$1,460.80<sup>EA</sup>**



⚙️ 1:1 Blue ●

⌚ Max. 40,000 rpm

○ Internal water

🔧 For CA-burs  
(Ø 2.35 mm)

• LP11L (F/Optic)

**\$1,287**

**BUY 3**  
**\$1,093.95<sup>EA</sup>**

**BUY 4**  
**\$1,029.60<sup>EA</sup>**



⚙️ 5.4:1 Green ●

⌚ Max. 7,407 rpm

○ Internal water

🔧 For CA-burs  
(Ø 2.35 mm)

• LP54L (F/Optic)

**\$1,320**

**BUY 3**  
**\$1,122.00<sup>EA</sup>**

**BUY 4**  
**\$1,056.00<sup>EA</sup>**



⚙️ 1:1 Blue ●

⌚ Max. 40,000 rpm

○ Internal water

🔧 For HP-burs  
(Ø 2.35 mm)

• LP01L (F/Optic)

**\$1,045**

**BUY 3**  
**\$888.25<sup>EA</sup>**

**BUY 4**  
**\$836.00<sup>EA</sup>**

AIR MOTORS ♡



- ✓ Internal water
- ✓ LED

• AM1006

**\$990**



- ✓ Internal water

• AM1004

**\$748**



- ✓ External water

• AM0014

**\$440**



• AM0004

**\$440**

**NEW**

KAVO COMPATIBLE HEADS ♡



- Kavo 68LH compatible
- Blue (1:1)
- LCH11

**\$440**



- Kavo 67LH compatible
- Green (2:1)
- LCH21

**\$440**



- Kavo 31LR compatible
- Green (2:1)
- LCH21P

**\$440**

## ECO LINE TURBINES



**Now with 1 Year FairCare Warranty<sup>1</sup>**  
Unconditional warranty except physical damage



- 62 dB Operating Noise
- AeroPerformance blade
- Refined ceramic bearings
- Glass rod fiber optics
- Smooth Eco body



Available for  
KaVo\*, NSK\*, W&H\*  
and Sirona\* Couplings



**BUY 3  
or MORE**  
**\$759**  
each

3 x Eco Line HS  
with F/O for any  
coupling connection

**SAVINGS FROM**  
**\$392**

Coupling Connection	Power Head 24W (F/Optic)	Small Head 16W (F/Optic)	Power Head 24W (N/Optic)	Small Head 16W (N/Optic)
KaVo	HE21KL	HE22KL	HE21K	HE22K
NSK	HE21NL	HE22NL	HE21N	HE22N
W&H	HE21WL	HE22WL	HE21W	HE22W
Sirona	HE21SL	- N/A -	- N/A -	- N/A -

**\$890 (FIBRE OPTIC)**

**\$649 (NON OPTIC)**

Buy 3  
**GET 10% OFF!**

## ECO LINE POWER TURBINE



60 dB Operating Noise

AeroPerformance blade

Refined ceramic bearings

Glass rod fiber optics

Smooth Eco body

**18 MONTH  
WARRANTY\***

• HE20KL **\$1,265** **BUY 3 \$1,089<sup>EA</sup>**

Available for only KaVo coupling

## ECO LINE FIXED CONNECTION TURBINES

- Direct 4 hole connection
- Refine ceramic bearings
- Push Button System
- HE17 - Standard head with 4 hole water spray
- HE18 - Small head with 3 hole water spray



**BUY 3  
GET  
10% OFF**  
**\$485<sup>10</sup>**  
EACH

**\$539**

1) FairCare warranty (except physical damage) applicable on Eco Line sold after 1st Jan 2019 2) Eco Line Power Head 24 Watt

**Buy 3 GET 15% OFF**  
**Buy 4 GET 20% OFF**  
 Mix & Match



**FairCare Warranty<sup>1</sup>**  
 1 year unconditional warranty  
 except physical damage



- ⚙️ 1:5 Red ●
- Internal water
- 🔧 For FG-burs  
(Ø 1.6 mm)

● **LE15L (F/Optic)**  
**\$1,485**

**BUY 3**  
 \$1,262.25<sup>EA</sup>

**BUY 4**  
 \$1,188.00<sup>EA</sup>

● **LE15 (Non-optic)**  
**\$1,287**

**BUY 3**  
 \$1,093.95<sup>EA</sup>

**BUY 4**  
 \$1,029.60<sup>EA</sup>



- ⚙️ 1:1 Blue ●
- Internal water
- 🔧 For CA-burs  
(Ø 2.35 mm)

● **LE11L (F/Optic)**  
**\$990**

**BUY 3**  
 \$841.50<sup>EA</sup>

**BUY 4**  
 \$792.00<sup>EA</sup>

● **LE11 (Non-optic)**  
**\$770**

**BUY 3**  
 \$654.50<sup>EA</sup>

**BUY 4**  
 \$616.00<sup>EA</sup>



- ⚙️ 5.4:1 Green ●
- Internal water
- 🔧 Screw-in & Snap-on  
Prophy Cups

● **LE54LP (F/Optic)**  
**\$1,210**

**BUY 3**  
 \$1,028.50<sup>EA</sup>

**BUY 4**  
 \$968.00<sup>EA</sup>

● **LE54P (Non-optic)**  
**\$935**

**BUY 3**  
 \$794.75<sup>EA</sup>

**BUY 4**  
 \$748.00<sup>EA</sup>



- ⚙️ 1:1 Blue ●
- Internal water
- 🔧 For HP-burs  
(Ø 2.35 mm)

● **LE01L (F/Optic)**  
**\$935**

**BUY 3**  
 \$794.75<sup>EA</sup>

**BUY 4**  
 \$748.00<sup>EA</sup>

● **LE01 (Non-optic)**  
**\$748**

**BUY 3**  
 \$635.80<sup>EA</sup>

**BUY 4**  
 \$598.40<sup>EA</sup>



**Buy More  
Save More**

- ⚙️ 1:1 Blue ●
- Push button

● **LB11/11**  
**\$418**

**BUY 3**  
 \$339.00<sup>EA</sup>

**BUY 4**  
 \$319.00<sup>EA</sup>



- ⚙️ 1:1 Blue ●
- Latch Type

● **LB11/11C**  
**\$418**

**BUY 3**  
 \$339.00<sup>EA</sup>

**BUY 4**  
 \$319.00<sup>EA</sup>



- ⚙️ 4:1 Green ●
- Snap-on Cups

● **LB41/11P**  
**\$462**

**BUY 3**  
 \$396.00<sup>EA</sup>

**BUY 4**  
 \$385.00<sup>EA</sup>



- ⚙️ 1:1 Blue ●

● **LB02**  
**\$418**

**BUY 3**  
 \$339.00<sup>EA</sup>

**BUY 4**  
 \$319.00<sup>EA</sup>

1) FairCare warranty (except physical damage) applicable on Eco Line sold after 1st Jan 2019

Buy 3 GET 1 FREE



EMS Non-optic  
Type EN-041

- SC21E

\$440

BUY 3  
GET 1  
FREE



EMS Piezon Non-optic  
Type EN-061

- SC21E2

\$539

BUY 3  
GET 1  
FREE



EMS Piezon LED  
Type EN-060A

- SC21EL

\$715

BUY 3  
GET 1  
FREE



KaVo Piezo LED  
Type 1.007.3995

- SC21KL\*

\$715

BUY 3  
GET 1  
FREE

NEW



Acteon-Satelec Non-optic  
Type F 12 281

- SC21S

\$440

BUY 3  
GET 1  
FREE



Acteon-Satelec LED  
Type F 12 609

- SC21SL

\$715

BUY 3  
GET 1  
FREE

\*Only for use with EMS compatible scaler tips

## SCALER TIPS

Buy 3 GET 1 FREE  
While stocks last

- ✓ High-quality Diamond-Like-Carbon coating
- ✓ Extremely hard and wear-resistant surface

### • SCALING DLC COATED Only available for Satelec fitting



### • PERIO DLC COATED Only available for Satelec fitting




### • ENDO DLC



### • SCALING STAINLESS STEEL Only available for Satelec fitting



### • PERIO STAINLESS STEEL Available for EMS & Satelec fittings



# Universal Torque Wrench

## TORQ CONTROL®

Precision instrument with a Torque Adjustment System for tightening of Prosthetic components at specific Torque, for any implant system.



### Automatic Declutching when Selected Torque is reached

- Compliance with Protocols for screw tightening
- Limits the risk of the screw fracturing or loosening



### Universal connection of Prosthetic mandrels

- 7 different Torque levels (from 10 to 35 N.cm)
- Can be used with ALL Implant Systems



### Micro-Head and Ergonomic Handle

- Easy access in restricted area with Excellent Visibility



### Easy and Quick to use

- Reduced motion compared with using a ratchet wrench while tightening



**\$1,210**


**Lightweight  
ONLY 135 gms**

Order code: 15501

## 3 Easy Steps



1. Please select the recommended Torque



2. Start rotating



3. Screw automatically stops when selected torque is reached



**Made in  
France**

[www.anthogyr.com](http://www.anthogyr.com)

**Anthogyr**  
A Straumann Group Brand

## Buy 1 Polisher

Get 1 Bottle of Propy Powder FREE

Get 50% DISCOUNT ON 1 ADDITIONAL TIP

• PR1011

**\$1,375**

- Ergonomically shaped, compact and light-weight air-driven polishing device

- 360° rotatable outer handle due to unique RotaGrip® for simple and fatigue-proof operation



- Propy Powder (300 gms) \$ 60.50 each

- Lemon - PP1012



### Additional Propy tip



• SP3011

**\$374**

### Coupling connection order codes:

Kavo - PR1011K, NSK - PR1011N, W&H - PR1011W, Sirona - PR1011S

## ✓ LEDS

## Buy 4 GET 10% OFF

- ✓ Quality illumination, No Glare
- ✓ 5500K Natural Daylight Colour
- ✓ Clear Vision & Sharper Contrast
- ✓ Longer Durability & Economical



**\$88**

- BU8012 - KaVo & MK-dent coupling
- BU8012KM - KaVo & MK-dent motors
- BU8012WH - W&H handpieces
- BU8012WLT - W&H airmotors LT25
- BU8012WEA - W&H motors 40LT, 50 LT & A51LT
- BU8012SG - Sirona coupling and motors
- BU8012N - NSK coupling
- BU8012NT - NSK motors TiM40 & Ti205L
- BU8012NM - NSK coupling KCL,SCL,PTL LED

## ✓ PREMIUM LUBRICATION OIL



- ✓ Fully Synthetic
- ✓ For any brand of handpieces and motors

- LU1011 (500 ml)

**\$68.20**



- ✓ Fully Synthetic
- ✓ For KaVo Quattrocare 2 lubrication system

- LU1022 (500 ml)

**\$68.20**



- ✓ For W&H Assistina lubrication system

- LU1031 (500 ml)

**\$88**

MANUFACTURER



**MK-DENT**

MK-dent GmbH

Marie-Curie-Str. 2

22941 Bargteheide

GERMANY

[www.mk-dent.com](http://www.mk-dent.com)

YOUR SPECIALIST DEALER

**ADAM DENTAL**  
supplies

• All catalogue prices are inclusive of GST • Promotion valid until 31st Dec 2023 or while stocks last

KaVo, W&H, Sirona, NSK, Bien-Air, EMS and Acteon-Satelec are registered trademarks. MK-dent GmbH has no economic connection with the aforementioned companies.