




PROMOTION
JAN - MAR 2024



AeroAcoustic blade



Glass rod fiber optics



Reverse flow protection



Automatic pressure regulation¹



Refined ceramic bearings



Scratch resistance
Iconic Prime body

25
W

Powerful & Precise
Power Head - 25W

55
dB

Incredibly Quiet
Only 55 dB noise



Available for
KaVo*, NSK*, W&H*
and Sirona* Couplings

Coupling
Connection

Kavo
NSK
W&H
Sirona

Power Head
25 Watt

HP21KL
HP21NL
HP21WL
HP21SL

Small Head
18 Watt

HP22KL
HP22NL
HP22WL
- N/A -

\$1,683

Only
\$1,183
with
Trade-In

Save \$500 with
Trade-In handpiece of

- ✓ Any brand
- ✓ Any type
- ✓ Any condition

**2 YEARS
UNCONDITIONAL
WARRANTY***

Excluding physical damage on Prime Line handpieces

COUPLINGS



- ✓ MK-dent & Kavo
with LED & water
- QC6016KW

\$495



- ✓ NSK connection
with LED
- QC6016NT

\$385



- ✓ W&H connection
with fibre optic
- QC5016W

\$385

- ✓ Sirona connection
with LED & water
- QC6016SW

\$495

- ✓ Non-optic
- QC4014K

\$296

- ✓ Non-optic
- QC4014NT

\$296

- ✓ Non-optic
- QC4014W

\$296

- ✓ Non-optic

-NA-

1) Only available for Kavo coupling 2) FairCare warranty applicable on Prime Line high speeds sold after 1st Oct 2018

PRIME LINE CONTRA - ANGLES

Buy 3 GET 15% OFF
Buy 4 GET 20% OFF
 Mix & Match

NEW RANGE!



2 Years FairCare Warranty
 Unconditional warranty
 except physical damage



SAFE-TECH
PUSH
BUTTON

1:5 Red ●

- Max. 200,000 rpm
- Internal water
- For FG-burs (Ø 1.6 mm)

• LP15L (F/Optic)

\$1,826

BUY 3
\$1,552.10^{EA}

BUY 4
\$1,460.80^{EA}



1:1 Blue ●

- Max. 40,000 rpm
- Internal water
- For CA-burs (Ø 2.35 mm)

• LP11L (F/Optic)

\$1,287

BUY 3
\$1,093.95^{EA}

BUY 4
\$1,029.60^{EA}



5.4:1 Green ●

- Max. 7,407 rpm
- Internal water
- For CA-burs (Ø 2.35 mm)

• LP54L (F/Optic)

\$1,320

BUY 3
\$1,122.00^{EA}

BUY 4
\$1,056.00^{EA}



1:1 Blue ●

- Max. 40,000 rpm
- Internal water
- For HP-burs (Ø 2.35 mm)

• LP01L (F/Optic)

\$1,045

BUY 3
\$888.25^{EA}

BUY 4
\$836.00^{EA}

AIR MOTORS



- Internal water
- LED

• AM1006

\$990

NEW
 Shorter
 and
 Lighter



- Internal water

• AM1004

\$748



- External water

• AM0014

\$473



• AM0004

\$473

NEW

KAVO COMPATIBLE HEADS



- Kavo 68LH compatible
- Blue (1:1)
- LCH11

\$440



- Kavo 67LH compatible
- Green (2:1)
- LCH21

\$440



- Kavo 31LR compatible
- Green (2:1)
- LCH21P

\$440

ECO LINE TURBINES



- 62 dB Operating Noise
- AeroPerformance blade
- Refined ceramic bearings
- Glass rod fiber optics
- Smooth Eco body

Available for
KaVo*, NSK*, W&H*
and Sirona* Couplings



Now with 1 Year FairCare Warranty¹
Unconditional warranty except physical damage

**BUY 3
or MORE
\$759
each**



**3 x Eco Line HS
with F/O for any
coupling connection**

**SAVINGS FROM
\$392**

Coupling Connection	Power Head 24W (F/Optic)	Small Head 16W (F/Optic)	Power Head 24W (N/Optic)	Small Head 16W (N/Optic)
KaVo	HE21KL	HE22KL	HE21K	HE22K
NSK	HE21NL	HE22NL	HE21N	HE22N
W&H	HE21WL	HE22WL	HE21W	HE22W
Sirona	HE21SL	- N/A -	- N/A -	- N/A -

\$890 (FIBRE OPTIC)

\$649 (NON OPTIC)

**Buy 3
GET 10% OFF!**

ECO LINE POWER TURBINE



- 28 Watt Power
- 60 dB Operating Noise
- AeroPerformance blade
- Glass rod fiber optics
- Refined ceramic bearings
- Smooth Eco body

• HE20KL
\$1,265 **BUY 3
\$1,089^{EA}**

Available for only KaVo coupling

ECO LINE FIXED CONNECTION TURBINES

- Direct 4 hole connection
- HE17 - Standard head with 4 hole water spray
- Refine ceramic bearings
- HE18 - Small head with 3 hole water spray
- Push Button System



1) FairCare warranty (except physical damage) applicable on Eco Line sold after 1st Jan 2019 2) Eco Line Power Head 24 Watt

ECO LINE CONTRA - ANGLES ♥

Buy 3 GET 15% OFF
Mix & Match



FairCare Warranty¹
1 year unconditional warranty
except physical damage



- ⚙️ 1:5 Red ●
- Internal water
- For FG-burs
(Ø 1.6 mm)

● LE15L (F/Optic)
\$1,485 **BUY 3**
\$1,262.25^{EA}



- ⚙️ 1:1 Blue ●
- Internal water
- For CA-burs
(Ø 2.35 mm)

● LE11L (F/Optic)
\$990 **BUY 3**
\$841.50^{EA}



- ⚙️ 5.4:1 Green ●
- Internal water
- For Screw-in and
Snap-on Prophylaxis Cups

● LE54LP (F/Optic)
\$1,210 **BUY 3**
\$1,028.50^{EA}



- ⚙️ 1:1 Blue ●
- Internal water
- For HP-burs
(Ø 2.35 mm)

● LE01L (F/Optic)
\$935 **BUY 3**
\$794.75^{EA}

● LE15 (Non-optic)
\$1,287 **BUY 3**
\$1,093.95^{EA}

● LE11 (Non-optic)
\$770 **BUY 3**
\$654.50^{EA}

● LE54P (Non-optic)
\$935 **BUY 3**
\$794.75^{EA}

● LE01 (Non-optic)
\$748 **BUY 3**
\$635.80^{EA}

BASIC LINE CONTRA - ANGLES ♥



**Buy More
Save More**



- ⚙️ 1:1 Blue ●
- Push button
- LB11/11

\$418 **BUY 3**
\$349.00^{EA}

- ⚙️ 1:1 Blue ●
- Latch Type
- LB11/11C

\$418 **BUY 3**
\$349^{EA}



- ⚙️ 4:1 Green ●
- Snap-on Cups
- LB41/11P

\$462 **BUY 3**
\$396^{EA}



- ⚙️ 1:1 Blue ●
- LB02

\$418 **BUY 3**
\$349^{EA}

1) FairCare warranty (except physical damage) applicable on Eco Line sold after 1st Jan 2019

ULTRASONIC SCALER

Buy 3 GET 15% OFF!



EMS Non-optic
Type EN-041

• SC21E

\$462

BUY 3
\$392.70^{EA}



EMS Piezon Non-optic
Type EN-061

• SC21E2

\$561

BUY 3
\$476.85^{EA}



EMS Piezon LED
Type EN-060A

• SC21EL

\$748

BUY 3
\$635.80^{EA}



KaVo Piezo LED
Type 1.007.3995

• SC21KL*

\$635.80

BUY 3
\$748^{EA}

NEW



Acteon-Satelec Non-optic
Type F 12 281

• SC21S

\$462

BUY 3
\$462^{EA}



Acteon-Satelec LED
Type F 12 609

• SC21SL

\$748

BUY 3
\$635.80^{EA}

*Only for use with EMS compatible scaler tips

SCALER TIPS

Buy 3 GET 1 FREE
While stocks last

- ✓ High-quality Diamond-Like-Carbon coating
- ✓ Extremely hard and wear-resistant surface

• SCALING DLC COATED Only available for Satelec fitting



ST02SD
SATELEC 1S

\$88



ST03SD
EMS A

\$88



ST04SD
SATELEC 2

\$132

• PERIO DLC COATED Only available for Satelec fitting



ST06SD
SATELEC 10Z

\$132



ST10SD
SATELEC 10P

\$154



ST08SD
SATELEC S07

\$154

• ENDO DLC

• SCALING STAINLESS STEEL Only available for Satelec fitting



ST02S
SATELEC 1S

\$63.80



ST03S
EMS A

\$63.80



ST04S
SATELEC 2

\$88

• PERIO STAINLESS STEEL Available for EMS & Satelec fittings



ST07E/ST07S
EMS DPL3
SATELEC 10X

\$110



ST09E/ST09S
EMS HPL

\$110



ST10E/ST10S
SATELEC 10P

\$110



ST11E/ST11S
EMS RT1

\$110

EMS and Acteon-Satelec are registered trademarks. MK-dent GmbH has no economic connection with the aforementioned companies.

Universal Torque Wrench

TORQ CONTROL®

Precision instrument with a Torque Adjustment System for tightening of Prosthetic components at specific Torque, for any implant system.



Automatic Declutching when Selected Torque is reached

- Compliance with Protocols for screw tightening
- Limits the risk of the screw fracturing or loosening




Universal connection of Prosthetic mandrels

- 7 different Torque levels (from 10 to 35 N.cm)
- Can be used with ALL Implant Systems



Micro-Head and Ergonomic Handle

- Easy access in restricted area with Excellent Visibility



Easy and Quick to use

- Reduced motion compared with using a ratchet wrench while tightening



UNIQUE

\$1,210


**Lightweight
ONLY 135 gms**

Order code: 15501

3 Easy Steps



1. Please select the recommended Torque



2. Start rotating



3. Screw automatically stops when selected torque is reached



**Made in
France**

www.anthogyr.com

Anthogyr
A Straumann Group Brand

▼ PROPHY LINE AIR POLISHING DEVICE

Buy 1 Polisher

Get 1 Bottle of Prophy Powder FREE

Get 50% DISCOUNT ON 1 ADDITIONAL TIP

● PR1011

\$1,375

- Ergonomically shaped, compact and light-weight air-driven polishing device

- 360° rotatable outer handle due to unique RotaGrip® for simple and fatigue-proof operation



- Prophy Powder (300 gms) \$ 60.50 each

- Lemon - PP1012

Additional Prophy tip



● SP3011

\$374

Coupling connection order codes:

Kavo - PR1011K, NSK - PR1011N, W&H - PR1011W, Sirona - PR1011S

▼ LEDS

Buy 4 GET 10% OFF

- ✓ Quality illumination, No Glare
- ✓ 5500K Natural Daylight Colour
- ✓ Clear Vision & Sharper Contrast
- ✓ Longer Durability & Economical



\$88

- BU8012 - KaVo & MK-dent coupling
- BU8012KM - KaVo & MK-dent motors
- BU8012WH - W&H handpieces
- BU8012WLT - W&H airmotors LT25
- BU8012WEA - W&H motors 40LT, 50 LT & A51LT
- BU8012SG - Sirona coupling and motors
- BU8012N - NSK coupling
- BU8012NT - NSK motors TiM40 & Ti205L
- BU8012NM - NSK coupling KCL,SCL,PTL LED

▼ PREMIUM LUBRICATION OIL



- ✓ Fully Synthetic
- ✓ For any brand of handpieces and motors



- LU1011 (500 ml)

\$68.20



- ✓ Fully Synthetic
- ✓ For KaVo Quattrocare 2 lubrication system

- LU1022 (500 ml)

\$68.20



- ✓ For W&H Assistina lubrication system

- LU1031 (500 ml)

\$88

MANUFACTURER



MK-DENT

MK-dent GmbH
Marie-Curie-Str. 2
22941 Bargteheide
GERMANY
www.mk-dent.com

YOUR SPECIALIST DEALER

ADAM DENTAL
supplies

• All catalogue prices are inclusive of GST • Promotion valid until 20th Mar 2024 or while stocks last

KaVo, W&H, Sirona, NSK, Bien-Air, EMS and Acteon-Satelec are registered trademarks. MK-dent GmbH has no economic connection with the aforementioned companies.