


PROTECTING THE OFFICE

Barriers For All Your Needs




SURFACES IN DENTAL PRACTICES ARE AT A HIGHER RISK

than previously thought of spreading infection during dental procedures!



Disease causing Hepatitis C Virus




Disease causing Hepatitis B Virus

Meet the Culprits



Infection causing Methicillin-resistant *Staphylococcus Aureus* (MRSA)




Disease causing Mycobacterium Tuberculosis var: bovis (BCG) (TB)

PLUS, *Candida*, *Clostridium Difficile* Spores, Vancomycin-resistant *Enterococci* (VRE)

3 BARRIERS WERE TESTED

See how effective each barrier was at reducing microbial penetration of methicillin-resistant *Staphylococcus Aureus* (MRSA)



2 Effective Ways

TO REDUCE CROSS CONTAMINATION:




Surface disinfectant wipes. Ideal for disinfecting hard non-porous materials.



Single-use disposable barrier covers. Best protects equipment that comes into patient contact.

Did You Know?

Organisms found on surfaces can survive longer than you may think²:



Not all Barriers ARE CREATED EQUAL


IMPERMEABLE BARRIER: Protects bacteria or infection from passing through.



PERMEABLE BARRIER: Bacteria or viruses can seep through.




6 Considerations FOR USING BARRIER COVERS:



¹ Protective Environmental Surface Barriers in Health Care. The Dental Advisor Number 29 - March 2016
² <http://psmed.net/wp-content/uploads/2015/06/Survival-Times.pdf>

Pinnacle™, Gendex™ & DEXIS™ X-Ray Sensor Sheaths

Disposable sheaths customized to fit various types of sensors.



PSS3100	X-Ray Sensor Sheaths - DEXIS Platinum	100/pk
PSS3310	X-Ray Sensor Sheaths - Schick CDR Size 1	100/pk
PSS3320	X-Ray Sensor Sheaths - Schick CDR Size 2	100/pk
PSS3340	X-Ray Sensor Sheaths - Schick CDR Elite Size 1	100/pk
PSS3350	X-Ray Sensor Sheaths - Schick CDR Elite Size 2,	100/pk
PSS3400	X-Ray Sensor Sheaths - Kodak Size 0	100/pk
PSS3410	X-Ray Sensor Sheaths - Kodak Size 1	100/pk
PSS3420	X-Ray Sensor Sheaths - Kodak Size 2	100/pk
3735	Sensor Cover 1 3/8" x 8" - Gendex	500/pk
3736	Sensor cover 1 5/8" x 8" - DEXIS	500/pk

Pinnacle™ & Allrap™ Chair Sleeve

Tapered design conforms neatly to your chair. It is disposable and reduces the need of chair clean up. Eliminates the need for a headrest cover and with reduced usage of disinfectants it helps extend the life of your chair.



3800	Universal chair cover 27 1/2" W x 24" L	225/pk
3825	Fits A-dec's Cascade & Other Similar Shapes - 32" W x 32" L	200/pk
3850-XL	Chair Sleeve, 29" W x 79" L	125/pk
3870-NS	Pinnacle No-Slip Chair Sleeve, 48" W x 56" L	150/pk
SGHC125	Allrap Chair Sleeve Half	4 boxes/case

Pinnacle™ & Allrap™ Headrest Covers

Disposable plastic cover for headrests.



3511	Pinnacle 11 1/2" W x 9 1/2" L	250/pk
3514	Pinnacle 14" W x 9 1/2" L	250/pk
HRCL250	Allrap 14 1/2" W x 10" L	12 boxes/case
HRCS250	Allrap 11 1/2" W x 10" L	12 boxes/case

Pinnacle™ Cover-All™

This thin sheet of adhesive plastic effectively covers many types of surfaces you might touch during a procedure. It is disposable, easy to apply and remove, and will not leave a residue. The Pinnacle Cover-All Dispenser is a table top dispenser that makes removal of films quick and effortless. Will dispense adhesive films up to 4" wide.



3700-C	Cover All Clear, 4" W x 6" L	1,200/pk
3710-B	Cover-All Blue, 4" W x 6" L	1,200/pk
3715-CN	Cover-All Clear, 2 1/2" W x 6" L	600/pk
AD100C	Cover-All Film Dispenser for 4" W x 6" L film	1

Pinnacle™ Light Sleeve

The Light Sleeve, offers two styles for the different handles, the "L" and the "T". They both allow the user to quickly slip the plastic sleeve over the handle. Its patented design keeps the sleeve from falling off while in use, yet it is easily removed when finished.



3600	T Style Pinnacle Light Sleeve, 5 3/4" W x 4" L	500/pk
3650	L Style Pinnacle Light Sleeve, 11 1/2" W x 8" L	500/pk
TLC500	Allrap Light Sleeve - 5 3/4" W x 4" L	36 boxes/case

Pinnacle™ & Allrap™ Syringe Sleeves



Disposable sleeve made to cover the 3-way syringe, HVE and saliva ejector valves.

Allrap Syringe Sleeve XL

Longer version of syringe sleeve designed to cover more of the dental tubing.

3760	Pinnacle 2 1/2" W x 10" L	500/pkg
3763-WA	Pinnacle with tape strip, 2 1/2" W x 11 1/2" L	500/pk
3767-WO	Pinnacle sleeve with opening	500/pkg
SS500C	Allrap Syringe Sleeve 2 1/2" x 10"	500/pkg
SS500XLC	Allrap Syringe Sleeve 2 1/2" x 20"	500/pkg
SS500	Allrap Syringe Sleeve 2 1/2" x 10"	24 boxes/case
HS500	Allrap Handpiece Sleeve - Universal Cover	60 boxes/case

Pinnacle™ & Allrap™ X-Ray Sleeves



The X-Ray Sleeve's tapered design protects your x-ray from the cross-contamination of surface disinfectants.

Pinnacle X-Ray Sleeve XL

Larger version of X-Ray Sleeve designed to cover larger and wider X-ray heads found on older machines.

3900	Pinnacle 15" W x 26" L	250/pkg
3950-XL	Pinnacle Extra-long, 23" W x 31" L	200/pkg
XRS225C	Allrap Universal 14" W x 24 1/2" L x 13 D"	225/pkg

Pinnacle™ Impression Gun Cover



Offers improved infection control, while keeping the gun clean of impression material. May also fit some large size curing lights.

4600	Impression Gun Cover	400/pkg
------	----------------------	---------

Pinnacle™ & Allrap™ Tray Sleeves



Instruments are placed directly on top of the disposable plastic sleeve, which is made to replace the paper tray cover used on instrument trays. Six custom sizes. Helps minimize staff time spent on cleaning the tray.

3000-B	Pinnacle 10 1/2" W x 14" L	500/pkg
3100-A	Pinnacle 11 5/8" W x 14 1/2" L	500/pkg
3200-E	Pinnacle F 11 1/2" W x 16" L	500/pkg
3300-F	Pinnacle 7 1/2" W x 10 1/2" L	500/pkg
3500-M	Pinnacle 14" W x 19 1/2" L	500/pkg
TS105-14	Allrap 10 1/2" W x 14" L	4 boxes/case
TS115-16	Allrap 11 5/8" W x 16" L	4 boxes/case

Pinnacle™ Side Shields



Optical Side Shields

Enhanced infection control, easy on, easy off design.

4700	Universal size, disposable, clear	250/pkg
------	-----------------------------------	---------

Cure Sleeve™ & Vision Saver™

Dental Cure Sleeve & Disposable Curing Shield

The cure sleeve prevents composite from bonding to end of tip, which decreases curing depth. The Vision Saver protects doctor and staff from the damaging effects of the intense light.

Kerr TotalCare is proud to carry Cure Sleeves. By using our Dental Cure Sleeve, all professionals can safely contribute to infection control and creating a clean environment for patients.




4500	Fits 7 and 8 mm tips	400/pkg
4511-1	Fits 11 mm tips	400/pkg
4513-1	Fits 13 mm tips	400/pkg
DEMI250	Fits DEMETRON®	250/pk
DEMIULTRA-100	Fits Demi™ Plus	100/pk
4575	Vision Saver, Orange	400/pkg

Googles™

Eyeshields

Easier lens attachment. No more fumbling. Larger lens for increased protection against infection, yet not bulky or boxy. Sleek, rounded frame.



GOP10-N	Office pack, Contains: 10 Frames, 20 Clear Lenses, 1 Yikes Neck Strap	100/pkg
GDLR100-N	Clear lenses, refill	25/pkg
GDLR25-N	Clear lenses, refill	
GFA10-N	Assorted frames, pack of 10 will include 2 of 5 colors	
GTDLR25-T	Tinted lenses	25/pkg
YOP5-N	Yikes neck straps	5/pkg

Bib-Eze™

Disposable Bib Clips

Bib-Eze Disposable Bib Holders are lightweight, stretchable paper holders that attach to the dental napkin bib with easy to peel, self-adhesive ends, which can easily be removed after the procedure and thrown away with the dirty napkin.

- Easy to use and helps reduce cross contamination
- Eliminate the cost of cleaning both in time and material.
- Comfortable with no sharp edges or cold metals against patient's neck
- Stretch to fit most patients



35053	Disposable Bib Clip	250/pkg
-------	---------------------	---------

Pinnacle™ & Allrap™ Chair Sleeve Applications



Chair Sleeve 3825



Chair Sleeve XL 3850-XL



No Slip Chair Sleeve 3870-NS



Chair Sleeve 3800



No Slip Chair Sleeve 3870-NS

Pinnacle™ & Allrap™ Tray Sleeve Applications



Tray Sleeve 3100-A



Tray Sleeve 3100-A



Tray Sleeve 3500-M



Tray Sleeve 3500-M



Tray Sleeve 3500-M

Pinnacle™ & Allrap™ X-Ray Sleeve Applications



X-Ray Sleeve 3900



X-Ray Sleeve XL 3950-XL


Pinnacle™ & Allrap™ Headrest Cover Applications



Headrest Cover 3511



Headrest Cover 3511



1-800-KERR-123 | kerrdental.com

MKT-20-0480 ©2020 Kerr Dental

Kerr[™]
TOTALCARE