

NEW:
Optimised
product range

Ivory® Rubber Dam

Your dam for complete moisture control.

Ivory Rubber Dam provides absolute isolation of the treatment area.

- » **100 % latex with no additional fillers:** It is made from natural latex and powder-free, as it has been thoroughly washed in the manufacturing stage to reduce surface proteins.
- » **New latex-free option:** Allows application despite a possible sensitization of users or patients to natural latex.
- » **Strong contrast, greater ease:** The new black rubber dam provides an even stronger contrast in visibility, easing and relaxing the work process.
- » **A hint of mint:** The new mint taste makes the treatment more pleasant for the patients and reduces a possible gag reflex.



Giving a hand to oral health.



KULZER
MITSUI CHEMICALS GROUP

Ivory® Rubber Dam

The comfort of perfect tooth isolation.

The use of rubber dam during restorative therapy and endodontic treatment has proved to make the treatment more successful. Rubber dam ensures a dry and clean working field which is needed for long-term accomplished adhesive restorations. During endodontic procedures rubber dam prevents endodontic failures caused by colonisation of saliva bacteria in the root canal.

Ivory Rubber Dam achieves complete moisture control in the treatment area due to its well-known properties:

- thin, comfortably pliable sheets are powder-free
- 100% latex with no additional fillers
- extremely high tensile strength*
- expansion up to 10 times from original punch hole*
- prevents irritation and damage of soft tissues

The relaunch brings further benefits to ensure a successful treatment:

- stronger contrast with black rubber dam
- latex-free rubber dam in case of sensitization
- pleasant mint taste (with green, black and non-latex rubber dam) for reduced gag irritation

*“Lovely dental dam!
Great elasticity and strength:
It is very easy to place also with
complicated posterior teeth.
Honestly, the color plays its part
and in this case colors of the
dam are great.”*



Art. Code	Colour	Thickness	Mint Taste	Latex Free	Content
Ivory Rubber Dam for adults (Size 6'x6')					
66094052	Blue	Thin			36 sheets
66094053	Blue	Medium			36 sheets
66094054	Blue	Heavy			36 sheets
66094055	Green	Medium	✓		36 sheets
66094056	Black	Medium	✓		36 sheets
66094057	Purple	Medium	✓	✓	15 sheets
Ivory Rubber Dam for children (Size 5'x5')					
66094058	Blue	Medium			52 sheets
66094059	Green	Medium	✓		52 sheets
66094060	Purple	Medium	✓	✓	15 sheets



Ivory Rubber Dam instruments for ideal tooth isolation – Clamps, Punch and Forceps

In addition, the Ivory Rubber Dam system ideally supports these benefits:

- made from high-grade stainless steel
- high performance and durability
- uniquely designed Ivory Clamps meet the most demanding situations
- enables confident fixation of the rubber dam



* Internal data. Data on file.