

Safety Data Sheet

Copyright,2022, 3M Company.All rights reserved. Copying and/or downloading of this information for the purpose of properly utilizing 3M products is allowed provided that: (1) the information is copied in full with no changes unless prior written agreement is obtained from 3M, and (2) neither the copy nor the original is resold or otherwise distributed with the intention of earning a profit thereon.

Document group:	08-7419-8	Version number:	6.00
Issue Date:	16/01/2022	Supersedes date:	25/04/2019

This Safety Data Sheet has been prepared in accordance with the Preparation of Safety Data Sheets for Hazardous Chemicals Code of Practice (Safe Work Australia, December 2011)

SECTION 1: Identification

1.1. Product identifier

3MTM FiltekTM P60 Posterior Restorative Paste

 Product Identification
 Numbers

 70-2010-2605-4
 70-2010-2606-2
 70-2010-2607-0

1.2. Recommended use and restrictions on use

Recommended use Dental Product, Restorative

Restrictions on use For use by dental professionals only.

1.3. Supplier's details

Address:	3M Australia - Building A, 1 Rivett Road, North Ryde NSW 2113
Telephone:	136 136
E Mail:	productinfo.au@mmm.com
Website:	www.3m.com.au

1.4. Emergency telephone number

EMERGENCY: 1800 097 146 (Australia only)

SECTION 2: Hazard identification

This product is classified as a hazardous chemical according to the Model Work Health and Safety Regulations, 2011, in accordance with applicable State and Territory legislation.

Refer to Section 14 of this Safety Data Sheets for product Dangerous Goods Classification.

2.1. Classification of the substance or mixture

Skin Sensitizer: Category 1.

2.2. Label elements

The label elements below were prepared in accordance with the Code of Practice on Preparation of Safety Data Sheets for Hazardous Chemicals (Safe Work Australia, December 2011). This information may be different from the actual product

label.

Signal word

Warning

Symbols Exclamation mark |

Pictograms



Hazard statements H317

May cause an allergic skin reaction.

Precautionary statements

Prevention: P261 P272	Avoid breathing dust/fume/gas/mist/vapours/spray. Contaminated work clothing should not be allowed out of the workplace.
Response: P302 + P352 P333 + P313 P362 + P364	IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse.
Disposal: P501	Dispose of contents/container in accordance with applicable local/regional/national/international regulations.

2.3. Other assigned/identified product hazards

None known.

2.4. Other hazards which do not result in classification

May be harmful if swallowed.

SECTION 3: Composition/information on ingredients

This material is a mixture.

Ingredient	CAS Nbr	% by Weight
Silane Treated Ceramic	444758-98-9	75 - 85
(1-methylethylidene)bis[4,1- phenyleneoxy(2-hydroxy-3,1-propanediyl)] bismethacrylate	1565-94-2	1 - 10
Bisphenol A Polyethylene Glycol Diether Dimethacrylate (BISEMA-6)	41637-38-1	1 - 10
Diurethane Dimethacrylate (UDMA)	72869-86-4	1 - 10
Aluminium oxide	1344-28-1	< 5
Triethylene Glycol Dimethacrylate (TEGDMA)	109-16-0	< 5

SECTION 4: First aid measures

4.1. Description of first aid measures

Inhalation

Remove person to fresh air. If you feel unwell, get medical attention.

Skin contact

Immediately wash with soap and water. Remove contaminated clothing and wash before reuse. If signs/symptoms develop, get medical attention.

Eye contact

Flush with large amounts of water. Remove contact lenses if easy to do. Continue rinsing. If signs/symptoms persist, get medical attention.

If swallowed

Rinse mouth. If you feel unwell, get medical attention.

4.2. Most important symptoms and effects, both acute and delayed

Allergic skin reaction (redness, swelling, blistering, and itching).

4.3. Indication of any immediate medical attention and special treatment required

Not applicable

SECTION 5: Fire-fighting measures

5.1. Suitable extinguishing media

In case of fire: Use a fire fighting agent suitable for ordinary combustible material such as water or foam to extinguish.

5.2. Special hazards arising from the substance or mixture

None inherent in this product.

Hazardous Decomposition or By-Products

<u>Substance</u> Carbon monoxide. Carbon dioxide.

5.3. Special protective actions for fire-fighters

Wear full protective clothing, including helmet, self-contained, positive pressure or pressure demand breathing apparatus, bunker coat and pants, bands around arms, waist and legs, face mask, and protective covering for exposed areas of the head.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

Evacuate area. Ventilate the area with fresh air. For large spill, or spills in confined spaces, provide mechanical ventilation to disperse or exhaust vapors, in accordance with good industrial hygiene practice. Refer to other sections of this SDS for information regarding physical and health hazards, respiratory protection, ventilation, and personal protective equipment.

6.2. Environmental precautions

Avoid release to the environment.

6.3. Methods and material for containment and cleaning up

Collect as much of the spilled material as possible. Place in a closed container approved for transportation by appropriate authorities. Clean up residue. Seal the container. Dispose of collected material as soon as possible in accordance with applicable local/regional/national/international regulations.

<u>Condition</u> During combustion. During combustion.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

A no-touch technique is recommended. If skin contact occurs, wash skin with soap and water. Acrylates may penetrate commonly-used gloves. If product contacts glove, remove and discard glove, wash hands immediately with soap and water and then re-glove. Avoid breathing dust/fume/gas/mist/vapours/spray. Do not eat, drink or smoke when using this product. Wash thoroughly after handling. Contaminated work clothing should not be allowed out of the workplace. Wash contaminated clothing before reuse. Do not get in eyes.

7.2. Conditions for safe storage including any incompatibilities

No special storage requirements.

SECTION 8: Exposure controls/personal protection

8.1 Control parameters

Occupational exposure limits

If a component is disclosed in section 3 but does not appear in the table below, an occupational exposure limit is not available for the component.

Ingredient	CAS Nbr	Agency	Limit type	Additional comments
Aluminium oxide	1344-28-1	Australia OELs	TWA(Inspirable dust)(8	
			hours):10 mg/m3	
Aluminum, insoluble compounds	1344-28-1	ACGIH	TWA(respirable fraction):1	A4: Not class. as human
			mg/m3	carcin
CAS NO SEQ117921	1344-28-1	ACGIH	TWA(inhalable	
			particulates):10 mg/m3	
CAS NO SEQ117922	1344-28-1	ACGIH	TWA(respirable particles):3	
			mg/m3	

ACGIH : American Conference of Governmental Industrial Hygienists

AIHA : American Industrial Hygiene Association

Australia OELs : Australia. Adopted National Exposure Standards for Atmospheric Contaminants in the Occupational Environment

CMRG : Chemical Manufacturer's Recommended Guidelines

TWA: Time-Weighted-Average

STEL: Short Term Exposure Limit

CEIL: Ceiling

Sen: Sensitiser

Sk: Absorption through the skin may be a significant source of exposure.

8.2. Exposure controls

8.2.1. Engineering controls

Use in a well-ventilated area.

8.2.2. Personal protective equipment (PPE)

Eye/face protection

Select and use eye/face protection to prevent contact based on the results of an exposure assessment. The following eye/face protection(s) are recommended:

Safety glasses with side shields.

Select and use eye protection in accordance with AS/NZS 1336. Eye protection should comply with the performance specifications of AS/NZS 1337.

Skin/hand protection

See Section 7.1 for additional information on skin protection.

Respiratory protection

None required.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

information on basic physical and chemical property	
Physical state	Solid.
Specific Physical Form:	Paste
Colour	White
Odour	Slight Acrylate
Odour threshold	No data available.
рН	Not applicable.
Melting point/Freezing point	No data available.
Boiling point/Initial boiling point/Boiling range	Not applicable.
Flash point	No flash point
Evaporation rate	Not applicable.
Flammability (solid, gas)	Not classified
Flammable Limits(LEL)	No data available.
Flammable Limits(UEL)	No data available.
Vapour pressure	Not applicable.
Vapor Density and/or Relative Vapor Density	Not applicable.
Density	2.1 g/cm3
Relative density	2.1 [<i>Ref Std</i> :WATER=1]
Water solubility	Negligible
Solubility- non-water	No data available.
Partition coefficient: n-octanol/water	Not applicable.
Autoignition temperature	No data available.
Decomposition temperature	No data available.
Viscosity/Kinematic Viscosity	Approximately 300,000 mPa-s
Volatile organic compounds (VOC)	No data available.
Percent volatile	No data available.
VOC less H2O & exempt solvents	No data available.
Molecular weight	No data available.

Nanoparticles

This material contains nanoparticles.

SECTION 10: Stability and reactivity

10.1 Reactivity

This material is considered to be non reactive under normal use conditions

10.2 Chemical stability

Stable.

10.3. Conditions to avoid

None known.

10.4. Possibility of hazardous reactions

Hazardous polymerisation will not occur.

10.5 Incompatible materials

None known.

10.6 Hazardous decomposition products

Substance None known. **Condition**

SECTION 11: Toxicological information

The information below may not be consistent with the material classification in Section 2 if specific ingredient classifications are mandated by a competent authority. In addition, toxicological data on ingredients may not be reflected in the material classification and/or the signs and symptoms of exposure, because an ingredient may be present below the threshold for labelling, an ingredient may not be available for exposure, or the data may not be relevant to the material as a whole.

11.1 Information on Toxicological effects

Signs and Symptoms of Exposure

Based on test data and/or information on the components, this material may produce the following health effects:

Inhalation

This product may have a characteristic odour; however, no adverse health effects are anticipated.

Skin contact

Contact with the skin during product use is not expected to result in significant irritation. Allergic skin reaction (non-photo induced): Signs/symptoms may include redness, swelling, blistering, and itching.

Eye contact

Contact with the eyes during product use is not expected to result in significant irritation.

Ingestion

May be harmful if swallowed.

Gastrointestinal irritation: Signs/symptoms may include abdominal pain, stomach upset, nausea, vomiting and diarrhoea.

Toxicological Data

If a component is disclosed in section 3 but does not appear in a table below, either no data are available for that endpoint or the data are not sufficient for classification.

Acute Toxicity

Name	Route	Species	Value
Overall product	Ingestion		No data available; calculated ATE >2,000 -
			≤5,000 mg/kg
Silane Treated Ceramic	Dermal		LD50 estimated to be > 5,000 mg/kg
Silane Treated Ceramic	Ingestion		LD50 estimated to be 2,000 - 5,000 mg/kg
Bisphenol A Polyethylene Glycol	Dermal	Professional	LD50 estimated to be $>$ 5,000 mg/kg
Diether Dimethacrylate (BISEMA-6)		judgement	
Diurethane Dimethacrylate (UDMA)	Dermal	Professional	LD50 estimated to be $> 5,000 \text{ mg/kg}$
		judgement	
Bisphenol A Polyethylene Glycol	Ingestion	Rat	LD50 > 2,000 mg/kg
Diether Dimethacrylate (BISEMA-6)			
Diurethane Dimethacrylate (UDMA)	Ingestion	Rat	LD50 > 5,000 mg/kg
(1-methylethylidene)bis[4,1-	Dermal	Professional	LD50 estimated to be $> 5,000 \text{ mg/kg}$
phenyleneoxy(2-hydroxy-3,1-		judgement	
propanediyl)] bismethacrylate			
(1-methylethylidene)bis[4,1-	Ingestion	Rat	LD50 > 11,700 mg/kg
phenyleneoxy(2-hydroxy-3,1-			

propanediyl)] bismethacrylate			
Triethylene Glycol Dimethacrylate (TEGDMA)	Dermal	Professional judgement	LD50 estimated to be > 5,000 mg/kg
Triethylene Glycol Dimethacrylate (TEGDMA)	Ingestion	Rat	LD50 10,837 mg/kg
Aluminium oxide	Dermal		LD50 estimated to be > 5,000 mg/kg
Aluminium oxide	Inhalation-Dust/Mist (4 hours)	Rat	LC50 > 2.3 mg/l
Aluminium oxide	Ingestion	Rat	LD50 > 5,000 mg/kg

ATE = acute toxicity estimate

Skin Corrosion/Irritation

Name	Species	Value
Silane Treated Ceramic	similar compounds	No significant irritation
(1-methylethylidene)bis[4,1-phenyleneoxy(2-	Rabbit	No significant irritation
hydroxy-3,1-propanediyl)] bismethacrylate		
Triethylene Glycol Dimethacrylate (TEGDMA)	Guinea pig	Mild irritant
Aluminium oxide	Rabbit	No significant irritation

Serious Eye Damage/Irritation

Name	Species	Value
Silane Treated Ceramic	similar compounds	Mild irritant
(1-methylethylidene)bis[4,1-phenyleneoxy(2-	In vitro data	No significant irritation
hydroxy-3,1-propanediyl)] bismethacrylate		
Triethylene Glycol Dimethacrylate (TEGDMA)	Professional judgement	Moderate irritant
Aluminium oxide	Rabbit	No significant irritation

Skin Sensitisation

Name	Species	Value
Silane Treated Ceramic	similar compounds	Not classified
Bisphenol A Polyethylene Glycol Diether Dimethacrylate (BISEMA-6)	Guinea pig	Not classified
	C in a nin	O
Diurethane Dimethacrylate (UDMA)	Guinea pig	Sensitising
(1-methylethylidene)bis[4,1-phenyleneoxy(2-	Mouse	Not classified
hydroxy-3,1-propanediyl)] bismethacrylate		
Triethylene Glycol Dimethacrylate (TEGDMA)	Human and animal	Sensitising

Respiratory Sensitisation

For the component/components, either no data are currently available or the data are not sufficient for classification.

Germ Cell Mutagenicity

Name	Route	Value
Bisphenol A Polyethylene Glycol Diether	In Vitro	Not mutagenic
Dimethacrylate (BISEMA-6)		
(1-methylethylidene)bis[4,1-phenyleneoxy(2-	In Vitro	Not mutagenic
hydroxy-3,1-propanediyl)] bismethacrylate		
Triethylene Glycol Dimethacrylate (TEGDMA)	In Vitro	Some positive data exist, but the data are not
		sufficient for classification
Aluminium oxide	In Vitro	Not mutagenic

Carcinogenicity

Name	Route	Species	Value
Silane Treated Ceramic	Inhalation	similar compounds	Some positive data exist, but the data

			are not sufficient for classification
Triethylene Glycol Dimethacrylate (TEGDMA)	Dermal	Mouse	Not carcinogenic
Aluminium oxide	Inhalation	Rat	Not carcinogenic

Reproductive Toxicity

Reproductive and/or Developmental Effects

Name	Route	Value	Species	Test result	Exposure Duration
(1- methylethylidene)bis[4,1-phenyleneoxy(2- hydroxy-3,1- propanediyl)] bismethacrylate	Ingestion	Not classified for development	Rat	NOAEL 1,000 mg/kg/day	during gestation
Triethylene Glycol Dimethacrylate (TEGDMA)	Ingestion	Not classified for female reproduction	Mouse	NOAEL 1 mg/kg/day	1 generation
Triethylene Glycol Dimethacrylate (TEGDMA)	Ingestion	Not classified for male reproduction	Mouse	NOAEL 1 mg/kg/day	1 generation
Triethylene Glycol Dimethacrylate (TEGDMA)	Ingestion	Not classified for development	Mouse	NOAEL 1 mg/kg/day	1 generation

Target Organ(s)

Specific Target Organ Toxicity - single exposure

For the component/components, either no data are currently available or the data are not sufficient for classification.

Specific Target Organ Toxicity - repeated exposure

Name	Route	Target Organ(s)	Value	Species	Test result	Exposure Duration
Silane Treated Ceramic	Inhalation	pulmonary fibrosis	Not classified	similar compounds	NOAEL Not available	
(1- methylethylid ene)bis[4,1- phenyleneoxy (2-hydroxy- 3,1- propanediyl)] bismethacryla te	Ingestion	endocrine system hematopoietic system liver heart skin gastrointestinal tract bone, teeth, nails, and/or hair immune system muscles nervous system eyes kidney and/or bladder respiratory system vascular system	Not classified	Rat	NOAEL 1,000 mg/kg/day	90 days
Triethylene Glycol Dimethacrylat e (TEGDMA)	Dermal	kidney and/or bladder blood	Not classified	Mouse	NOAEL 833 mg/kg/day	78 weeks
Aluminium oxide	Inhalation	pneumoconiosis	Some positive data exist, but the data are not sufficient for	Human	NOAEL Not available	occupational exposure

			classification			
Aluminium	Inhalation	pulmonary	Not classified	Human	NOAEL Not	occupational
oxide		fibrosis			available	exposure

Aspiration Hazard

For the component/components, either no data are currently available or the data are not sufficient for classification.

Exposure Levels

Refer Section 8.1 Control Parameters of this Safety Data Sheet.

Interactive Effects

Not determined.

SECTION 12: Ecological information

The information below may not be consistent with the material classification in Section 2 if specific ingredient classifications are mandated by a competent authority. Additional information leading to material classification in Section 2 is available upon request. In addition, environmental fate and effects data on ingredients may not be reflected in this section because an ingredient is present below the threshold for labelling, an ingredient is not expected to be available for exposure, or the data is considered not relevant to the material as a whole.

12.1. Toxicity

Acute aquatic hazard:

Not acutely toxic to aquatic life by GHS criteria.

Chronic aquatic hazard:

Not chronically toxic to aquatic life by GHS criteria.

No product test data available.

Material	CAS Number	Organism	Туре	Exposure	Test endpoint	Test result
Silane Treated Ceramic	444758-98-9		Data not available or insufficient for classification			N/A
(1- methylethylide ne)bis[4,1- phenyleneoxy(2-hydroxy-3,1- propanediyl)] bismethacrylate	1565-94-2	Common Carp	Analogous Compound	96 hours	No tox obs at lmt of water sol	>100 mg/l
(1- methylethylide ne)bis[4,1- phenyleneoxy(2-hydroxy-3,1- propanediyl)] bismethacrylate	1565-94-2	Green Algae	Endpoint not reached	96 hours	EC50	>100 mg/l
(1- methylethylide ne)bis[4,1- phenyleneoxy(2-hydroxy-3,1-	1565-94-2	Green Algae	Experimental	96 hours	EC10	1.1 mg/l

nronan-di-1)1		1	1	1	1	
propanediyl)] bismethacrylate						
Bisphenol A	41637-38-1	Activated	Estimated	3 hours	EC50	>1,000 mg/l
Polyethylene	41037-38-1	sludge	Estimated	5 nours	ECSU	~1,000 mg/1
Glycol Diether		sludge				
Dimethacrylate						
(BISEMA-6)						100 /
Bisphenol A	41637-38-1	Green algae	Estimated	72 hours	No tox obs at	>100 mg/l
Polyethylene					lmt of water sol	
Glycol Diether						
Dimethacrylate						
(BISEMA-6)						
Bisphenol A	41637-38-1	Rainbow trout	Estimated	96 hours	No tox obs at	>100 mg/l
Polyethylene					lmt of water sol	
Glycol Diether						
Dimethacrylate						
(BISEMA-6)						
Bisphenol A	41637-38-1	Green algae	Estimated	72 hours	No tox obs at	>100 mg/l
Polyethylene				_	lmt of water sol	
Glycol Diether						
Dimethacrylate						
(BISEMA-6)						
Diurethane	72869-86-4	Green algae	Endpoint not	72 hours	ErC50	>100 mg/l
Dimethacrylate	72007-00-4	Green argae	reached	72 110013		> 100 mg/1
(UDMA)			reaction			
Diurethane	72869-86-4	Water flea	Experimental	48 hours	EC50	>100 mg/l
Dimethacrylate	72809-80-4	water fiea	Experimental	48 110015	LC.30	>100 mg/1
(UDMA) Diurethane	72869-86-4	Zebra Fish	Experimental	96 hours	LC50	$10.1 m a^{/1}$
	/2809-80-4	Zeora risii	Experimental	90 nours	LC30	10.1 mg/l
Dimethacrylate						
(UDMA)	70060.06.4			70.1	E C10	. 100 //
Diurethane	72869-86-4	Green algae	Endpoint not	72 hours	ErC10	>100 mg/l
Dimethacrylate			reached			
(UDMA)						
Aluminium	1344-28-1	Fish	Experimental	96 hours	LC50	>100 mg/l
oxide						
Aluminium	1344-28-1	Green Algae	Experimental	72 hours	EC50	>100 mg/l
oxide						
Aluminium	1344-28-1	Water flea	Experimental	48 hours	LC50	>100 mg/l
oxide			-			-
Aluminium	1344-28-1	Green Algae	Experimental	72 hours	NOEC	>100 mg/l
oxide		0	F			
Triethylene	109-16-0	Green Algae	Experimental	72 hours	EC50	>100 mg/l
Glycol				2 110 010		1.00 mg/1
Dimethacrylate						
(TEGDMA)						
Triethylene	109-16-0	Zebra Fish	Experimental	96 hours	LC50	16.4 mg/l
Glycol				50 110415		10.1116/1
Dimethacrylate						
(TEGDMA)						
	109-16-0	Groop alaga	Exporimonto ¹	72 hours	NOEC	18.6 mg/l
Triethylene	109-10-0	Green algae	Experimental	72 hours	NUEC	18.6 mg/l
Glycol						
Dimethacrylate						
(TEGDMA)						
Triethylene	109-16-0	Water flea	Experimental	21 days	NOEC	32 mg/l

Glycol			
Dimethacrylate			
(TEGDMA)			

12.2. Persistence and degradability

Material	CAS Number	Test type	Duration	Study Type	Test result	Protocol
Silane Treated Ceramic	444758-98-9	Data not available-			N/A	
Cerainic		insufficient				
(1- methylethylide ne)bis[4,1- phenyleneoxy(2-hydroxy-3,1- propanediyl)] bismethacrylate		Experimental Hydrolysis		Hydrolytic half-life (pH 7)	29 days (t 1/2)	
(1- methylethylide ne)bis[4,1- phenyleneoxy(2-hydroxy-3,1- propanediyl)] bismethacrylate	1565-94-2	Experimental Biodegradation	28 days	BOD	21 % BOD/ThBOD	similar to OECD 301F
Bisphenol A Polyethylene Glycol Diether Dimethacrylate (BISEMA-6)	41637-38-1	Experimental Biodegradation	28 days	BOD	24 % BOD/ThBOD	OECD 301D - Closed bottle test
Diurethane Dimethacrylate (UDMA)	72869-86-4	Experimental Biodegradation	28 days	CO2 evolution	22 %CO2 evolution/THC O2 evolution (does not pass 10-day window)	OECD 301B - Modified sturm or CO2
Aluminium oxide	1344-28-1	Data not available- insufficient			N/A	
Triethylene Glycol Dimethacrylate (TEGDMA)	109-16-0	Experimental Biodegradation	28 days	CO2 evolution	85 % weight	OECD 301B - Modified sturm or CO2

12.3 : Bioaccumulative potential

Material	CAS Number	Test type	Duration	Study Type	Test result	Protocol
Silane Treated	444758-98-9	Data not	N/A	N/A	N/A	N/A
Ceramic		available or				
		insufficient for				
		classification				
(1-	1565-94-2	Experimental		Log Kow	4.63	
methylethylide		Bioconcentrati				
ne)bis[4,1-		on				
phenyleneoxy(
2-hydroxy-3,1-						

1: 1)]		1	1			
propanediyl)]						
bismethacrylate						
Bisphenol A	41637-38-1	Estimated		Bioaccumulatio	6.6	Estimated:
Polyethylene		Bioconcentrati		n factor		Bioconcentration factor
Glycol Diether		on				
Dimethacrylate						
(BISEMA-6)						
Bisphenol A	41637-38-1	Experimental		Log Kow	≥4.66	OECD 117 log Kow
Polyethylene		Bioconcentrati				HPLC method
Glycol Diether		on				
Dimethacrylate						
(BISEMA-6)						
Diurethane	72869-86-4	Experimental		Log Kow	3.39	Non-standard method
Dimethacrylate		Bioconcentrati				
(UDMA)		on				
Aluminium	1344-28-1	Data not	N/A	N/A	N/A	N/A
oxide		available or				
		insufficient for				
		classification				
Triethylene	109-16-0	Experimental		Log Kow	2.3	Non-standard method
Glycol		Bioconcentrati		_		
Dimethacrylate		on				
(TEGDMA)						

12.4. Mobility in soil

Please contact manufacturer for more details

12.5 Other adverse effects

No information available.

SECTION 13: Disposal considerations

13.1. Disposal methods

Dispose of contents/ container in accordance with the local/regional/national/international regulations.

As a disposal alternative, incinerate in a permitted waste incineration facility. Dispose of waste product in a permitted industrial waste facility.

SECTION 14: Transport Information

Australian Dangerous Goods Code (ADG) - Road/Rail Transport

UN No.: Not applicable. Proper shipping name: Not applicable. Class/Division: Not applicable. Sub Risk: Not applicable. Packing Group: Not applicable.

Hazchem Code: Not applicable IERG: Not applicable.

International Air Transport Association (IATA) - Air Transport UN No.: Not applicable. Proper shipping name: Not applicable. Class/Division: Not applicable. **Sub Risk:** Not applicable. **Packing Group:** Not applicable.

International Maritime Dangerous Goods Code (IMDG)- Marine Transport

UN No.: Not applicable.
Proper shipping name: Not applicable.
Class/Division: Not applicable.
Sub Risk: Not applicable.
Packing Group: Not applicable.
Marine Pollutant: Not applicable.

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

Australian Inventory Status:

This product is regulated by the Therapeutics Goods Administration and is exempt from compliance with the Industrial Chemicals (Notification and Assessment) Act 1989 as amended.

SECTION 16: Other information

Revision information:

Complete document review.

DISCLAIMER: The information on this Safety Data Sheet is based on our experience and is correct to the best of our knowledge at the date of publication, but we do not accept any liability for any loss, damage or injury resulting from its use (except as required by law). The information may not be valid for any use not referred to in this Safety Data Sheet or use of the product in combination with other materials. For these reasons, it is important that customers carry out their own test to satisfy themselves as to the suitability of the product for their own intended applications. Greenguard ® is a United States based program. The 'Low VOC' reference related to United States Federal and State

Greenguard ® is a United States based program. The 'Low VOC' reference related to United States Federal and State regulations exemptions for some solvents.

3M Australia SDSs are available at www.3m.com.au