## MiniBlaster

Autoclavable Clinical Sandblaster



## **Installation Options**

## MiniBlaster Autoclavable Clinical Sandblaster



B Air Source: Dental unit 35-60 PSI (standard handpiece tubing)

	Connector	Description	Item No.
Ideally the Sandblaster is connected directly to the compressor at 80-120 PSI, however it can simply be connected to the high speed or low speed tubing just like a hand piece using one of the special adaptors available and can then be operated by the standard foot control of the dental unit.		For 2-3 hole connectors.	505112
		For 4 hole connectors.	505111
		Suitable for Sirona® (Siemens) dental units.	505114
A range of connections are available for the most popular dental units such as Kavo®, Sirona® (Siemens®) & W&H® in addition to standard 2-3 hole and 4 hole connectors.		Suitable for Kavo® dental units.	505113
		Suitable for W&H® dental units.	505117

This connection is supplied with the unit