


Monovek[®]

**Dyed (violet)
monofilament
polydioxanone
(PDO)**

Synthetic
Long-term Absorbable



Suggested procedures

- Orthopaedics
- Plastic
- General / Intestinal
- Paediatric
- Oral and maxillofacial
- Dental
- Ophthalmic

Features


















- Monofilament provides smooth passage and minimal tissue trauma
- Optimum wound support
- Ideal suture when extended wound support is required
- Excellent long term tensile strength
- Superior knot run down and knotting capabilities
- Uniform absorption through hydrolysis provides predictable absorption irrespective of patient health

Tensile Strength	75%	65%	50%
Days	14	28	49
Complete Absorption	6 months		

- Available in:
 - a range of needles including heavy needles
 - sizes 6/0 to USP 2
 - sterile pre-cut lengths of various sizes
- Sutures made with 300 series stainless steel
- Up to five-year shelf life

Commonly used Monovek[®] product codes and specifications.

















Contact us for any other suture specification.

Needle*	Thread							
	USP Metric	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4
3/8 Circle Reverse Cutting								
10mm 	45cm	NAM601	NAM501					
10mm 	45cm	NADM601	NADM501					
12mm 	45cm	NAM602	NAM502					
13mm 	45cm	NAM6040	NAM5040	NAM4040				
	76cm	NAM16040	NAM15040	NAM14040				
35mm 	45cm					NAM2031	NAM031	NAM131
	76cm					NAM12031	NAM1031	NAM1131
40mm 	45cm					NAM2013	NAM013	NAM113
	76cm					NAM12013	NAM1013	NAM1113
45mm 	76cm					NAM12010	NAM1010	NAM1110
3/8 Circle Reverse Premium Cutting								
11mm 	45cm	NAM602D	NAM502D					
13mm 	76cm		NAM15049	NAM14049				
16mm 	45cm	NAM604	NAM504	NAM404	NAM304			
	76cm		NAM1504	NAM1404	NAM1304			
19mm 	45cm		NAM505	NAM405	NAM305			
	76cm		NAM1505	NAM1405	NAM1305			
24mm 	45cm			NAM406	NAM306	NAM206		
	76cm			NAM1406	NAM1306	NAM1206		
26mm 	45cm				NAM3047	NAM2047		
	76cm				NAM13047	NAM12047		
30mm 	45cm					NAM209		
	76cm					NAM1209		
3/8 Circle Conventional Cutting								
13mm 	45cm	NAM6093	NAM5093	NAM4093				
16mm 	76cm	NAM16037	NAM15037	NAM14037				
24mm 	76cm				NAM13023	NAM12023		

* Needle images and sizes are approximations only.

Commonly used Monovek[®] product codes and specifications.

Contact us for any other suture specification.

Needle*	Thread								
	USP Metric	6/0 0.7	5/0 1	4/0 1.5	3/0 2	2/0 3	0 3.5	1 4	2 5
1/2 Circle Reverse Cutting									
25mm 	45cm				NAM3034	NAM2034			
	76cm				NAM13034	NAM12034			
30mm 	45cm				NAM3038	NAM2038			
	76cm				NAM13038	NAM12038			
35mm 	45cm					NAM2017	NAM017		
	76cm				NAM13017	NAM12017	NAM1017	NAM1117	
40mm 	76cm					NAM12012	NAM1012	NAM1112	NAM1212
	90cm					NAM2012S	NAM012S	NAM112S	NAM212S
48mm 	76cm						NAM1045		
	90cm						NAM045S	NAM145S	NAM245S
1/2 Circle Conventional Cutting									
22mm 	45cm				NAM3018				
1/2 Circle Round Bodied Taper									
13mm 	76cm	NADM6055	NADM5055						
16mm 	76cm		NAM5060	NAM4060	NAM3060	NAM2060			
18mm 	76cm		NAM5064	NAM4064	NAM3064	NAM2064			
20mm 	76cm		NAM5065	NAM4065	NAM3065	NAM2065			
22mm 	76cm		NAM50#23	NAM40#23	NAM30#23	NAM20#23			
25mm 	76cm			NAM4068	NAM3068	NAM2068	NAM068		
30mm 	76cm			NAM4069	NAM3069	NAM2069	NAM069	NAM169	NAM269
	90cm			NAM14069	NAM13069	NAM12069	NAM1069	NAM1169	NAM1269
35mm 	76cm				NAM3070	NAM2070	NAM070	NAM170	
	90cm				NAM13070	NAM12070	NAM1070	NAM1170	
40mm 	76cm					NAM2097	NAM097	NAM197	NAM297
	90cm					NAM12097	NAM1097	NAM1197	NAM1297
48mm 	76cm					NAM2098	NAM098	NAM198	NAM298
	90cm					NAM12098	NAM1098	NAM1198	NAM1298

* Needle images and sizes are approximations only.

Proudly Australian owned and operated.

Dynek is a leading, family owned manufacturer of sutures, committed to technical and surgical excellence. Established in 1974, Dynek has a proud history of innovation and quality as the only Australian manufacturer of sutures, with a singular focus of excellent patient outcomes.



dynek.com

AU: 1800 088 172

NZ: 0800 459 649

INT: +61 8 8268 2033

sales@dynek.com



- The information contained in this document is intended as a guide only.
- Always refer to the instructions for use when using Dynek sutures.
- Monovek® is included on the ARTG (AU), Medsafe (NZ), and CE Mark (EU).
- Knot tying requires the standard surgical technique of flat and square ties with additional throws as indicated by surgical circumstance.
- A minimum order quantity may apply to all products.

QA: BE074/2 06/2022