

**Dealer information** 

**Venus Bulk Flow ONE** 



### **Product Description**



# VENUS BULK FLOW ONE ONE SHADING MEETS BULK FILLING.

#### **Venus Bulk Flow ONE simplifies and streamlines the restorative process**

As a further universal shade solution for efficient restorations in the posterior region, it represents the first flowable **one-shade bulk-filling composite** in the dental market that adapts to any tooth shade from A1 to D4.



• Always grab the right shade: With the incredible shadematching properties of Venus Bulk Flow ONE the restoration becomes invisible by blending into the surrounding dentition.



• Efficient handling: Cavities are filled quickly and easily with increments of up to 4 mm without the need for an additional capping layer.



• Long-lasting restorations: Low shrinkage stress, high flexural strength and high wear resistance ensure durable restorations that withstand the high mastication load in the posterior region.

Learn more about our new ONE: watch movie >>>

What makes **Venus Bulk Flow ONE** *unique*?

It is the first flowable bulk fill composite that can be used without additional capping.





Source: Dr. Janine Schweppe, Germany

Light-curing times and maximum layer thickness Venus Bulk Flow ONE:

- Curing depth 4 mm
- Curing time 20 sec

#### Indications for Venus Bulk Flow ONE:

- Class I, II, III and V direct restorations
- First layer lining for class I and II cavities
- Repair of direct and indirect restorations in combination with a suitable adhesive
- Splinting of loosened teeth resulting from trauma or periodontal associated events
- Extended fissure sealing.

## Selling units & pricing: Refills and Value-Packs

Venus Bulk Flow ONE - Refills			RRP AUD
66094723	VENUS BULK FLOW ONE SYR 1 x 2G	1 x 2g syringe 5 x cannulas 1 x IFU	\$55.40
66094724	VENUS BULK FLOW ONE PLT 2x10x0,20G	2 x 10 x 0,2g PLT 1 x IFU	\$118.89

Venus Bulk Flow ONE – Value Packs			RRP AUD
66095335	VENUS BULK FLOW ONE VALUE-KIT SYRINGE	10 x 2g syringes 10 x 5 cannulas 1 x pictorial card, 1 x IFU	\$458.73
66095388	VENUS BULK FLOW ONE VALUE-KIT PLT	80 x 0.20g PLT 1 x pictorial card, 1 x IFU	\$426.98







